

3914 DEL AMO BLVD SUITE 901 TORRANCE CA 90503 www.celprogen.com Phone: 310 542 8822 Fax: 310 542 8028 Email: info@celprogen.com stemcells@celprogen.com

DATA SHEET

Product name: Rat (Adipose Derived) Culture – T225 Plated Cells

Catalog number: 55010-11-T225

Description: (Adipose Derived) Culture T225 Plated Cells. Also available in T25, T75 and

T150 tissue culture flask with plated Cells, shipped at room temperature or in Frozen Vial, shipped in Dry-ice. They were maintained in Celprogen's (Adipose Derived) Culture Complete Growth Medium and sub-cultured every 24 to 48

hours on (Adipose Derived) Culture Extra-cellular Matrix.

Source: Rat Lipo-aspriate

Storage Conditions: Liquid nitrogen vapor phase for frozen Ampule of (Adipose Derived) Culture. For

plated cells in tissue culture flask, upon receipt of the cells wipe the flask with 70% ethanol and transfer to sterile tissue culture hood. In the tissue culture hood remove the media of the cells and wash the cells with 1X PBS sterile solution, for 2-3 minutes, remove the PBS solution and then trypsinize. After trypsinization of the Cells neutralize the trypsin with equal volume of (Adipose Derived) Culture Complete Growth Media with serum and collect the Cell suspension in sterile conical centrifuge tube in the tissue culture hood. Centrifuge the cell suspension at 100g for 7 minutes in centrifuge. Plate cells $5x10^5$ cells per pre-coated flasks with (Adipose Derived) Culture Extra-cellular Matrix for Expansion in (Adipose

Derived) Culture Complete Growth Medium.

Positive Markers: CD34, CD45, CD90, CD117 & Nestin

Morphology & Proliferation:

Mixed population of cells with approximately 95% attached cells and the other 5.0% in suspension, need to change cell culture media every day after 48 hours of

initial cell culture or when the media starts changing color to slight yellow for pink. Fast growing cell culture. Change media with Celprogen's (Adipose

Derived) Culture Complete Growth Medium with the appropriate (Adipose Derived) Culture Extra-cellular Matrix. Temperature 37^oC in 5% CO₂ humidified

incubator.



1871 N GAFFEY STREET SUITE A & B SAN PEDRO CA 90731 www.celprogen.com Phone: 310 547 3975 Fax: 310 547 2975 Email: info@celprogen.com stemcells@celprogen.com

Sub-culturing:

- 1. Thaw the vial with gentle agitation in a 37°C water bath or a dry 37°C shaking incubator. For water bath thawing keep the O-ring out of the water.
- 2. Remove the thawed vial and wipe with 70% ethanol. Then transfer to the tissue culture hood.
- 3. Transfer the vial contents to a sterile centrifuge tube, and gently add pre-warmed (Adipose Derived) Culture Growth Media to the centrifuge tube. Use additional (Adipose Derived) Culture Complete Media to rinse the vial and transfer the liquid to the centrifuge tube repeat this once more to ensure you have all the cells transferred to the 15ml centrifuge tube. Centrifuge the cells at 100g for 7 minutes. Remove the supernatant and re-suspend the cell pellet in 500ul of (Adipose Derived) Culture Complete Growth Medium.
- **4.** Add the 500ul of cells to T225 flask pre-coated with (Adipose Derived) Culture Extra-cellular Matrix with 15ml of (Adipose Derived) Culture Complete Growth Medium.
- **5.** Incubate the cells in the T225 flask in a 37^oC in 5% CO₂ humidified incubator. Perform 100% Media Change every 24 to 48 hours.
- **6.** Medium renewal every other or 2-3 days, sub-culturing ratio: 1:3

Freezing Medium: Available for purchase Cat# M55010-11FM

Trypsin: Available for purchase Cat# T1509-014

IX PBS: Available for Purchase Cat# P1408-013

Storage temperature: Liquid nitrogen vapor phase

Product Orders: Before submitting an order you will be asked to read and accept the terms and

conditions of Celprogen's Material Transfer Agreement (MTA).

Permits/Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory

permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the

permits.

Notices

& Disclaimers: CELPROGEN products are intended for laboratory research purposes only. They

are not intended for use in Humans. Product established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by

CELPROGEN Inc.