



1871 N GAFFEY STREET
SUITE A & B
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

Rat Gastric Stem Cell Data Sheet

Product name: Rat Gastric Stem Cell – Frozen Vial
--

Catalog number: 55106-33

Description: Frozen Ampule (1.2×10^6 Cells) of 1×10^6 Viable cell upon thawing, shipped with dry-ice. Also available in T25 and T75 tissue culture flask with plated cells, shipped at room temperature. They were maintained in Celprogen's Rat Gastric Stem Cell Complete Growth Medium and sub-cultured every 24 to 48 hours on Rat Gastric Stem Cell Extra-cellular Matrix.

Source: Rat Gastric tissue

Storage Conditions: Liquid nitrogen vapor phase for frozen Ampule of Rat Gastric Stem Cell. For plated cells in T75 tissue culture flask, upon receipt of the cells wipe the flask with 70% ethanol and transfer to sterile tissue culture hood. In the tissue culture hood remove the media of the cells and wash the cells with 1X PBS sterile solution, for 2 -3 minutes, remove the PBS solution and then trypsinize. After trypsinization of the Cells neutralize the trypsin with equal volume of Rat Gastric Stem Cell Complete Growth Media with serum and collect the Cell suspension in sterile conical centrifuge tube in the tissue culture hood. Centrifuge the cell suspension at 100g for 7 minutes in centrifuge. Plate cells 5×10^5 cells per pre-coated flasks with Rat Gastric Stem Cell Extra-cellular Matrix for Expansion in Rat Gastric Stem Cell Complete Growth Medium.

Positive markers: CK 19, CD 34, Alkaline Phosphatase, HNF-1B, epithelial membrane antigen (EMA), mucin 5AC, cytokeratin-18 (CK-18), trefoil factor 1 (TFF-1) and Oct-4

Morphology & Proliferation: Mixed population of cells with approximately 95% attached cells and the other 5.0% in suspension, need to change cell culture media every day after 48 hours of initial cell culture or when the media starts changing color to slight yellow for pink. Fast growing cell culture. Change media with Celprogen's Rat Gastric Stem Cell Culture Complete Growth Medium with the appropriate Rat Gastric Stem Cell Extra-cellular Matrix. Temperature 37°C in 5% CO_2 humidified incubator.

Sub-culturing:

1. Thaw the vial with gentle agitation in a 37°C water bath or a dry 37°C shaking incubator. For water bath thawing keep the O-ring out of the water.



1871 N GAFFEY STREET
SUITE A & B
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

2. Remove the thawed vial and wipe with 70% ethanol. Then transfer to the tissue culture hood.
3. Transfer the vial contents to a sterile centrifuge tube, and gently add pre-warmed Rat Gastric Stem Cell Growth Media to the centrifuge tube. Use additional Rat Gastric Stem Cell Culture Complete Media to rinse the vial and transfer the liquid to the centrifuge tube repeat this once more to ensure you have all the cells transferred to the 15ml centrifuge tube. Centrifuge the cells at 100g for 7 minutes. Remove the supernatant and re-suspend the cell pellet in 500ul of Rat Gastric Stem Cell Culture Complete Growth Medium.
4. Add the 500ul of cells to T75 flask pre-coated with Rat Gastric Stem Cell Extra-cellular Matrix with 15ml of Rat Gastric Stem Cell Complete Growth Medium.
5. Incubate the cells in the T75 flask in a 37⁰C in 5% CO₂ humidified incubator. Perform 100% Media Change every 24 to 48 hours.
6. Medium renewal every other or 2-3 days, sub-culturing ratio: 1:3

Freezing Medium: Available for purchase Cat# M55106-33FM

Trypsin: Available for purchase Cat# T1509-014

IX PBS: Available for Purchase Cat# P1408-013

Storage temperature: Liquid nitrogen vapor phase

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits/Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices

& Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans. The product, Rat Gastric Stem Cell, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.*