

3914 DEL AMO BLVD. SUITE 901 TORRANCE, CA 90503 www.celprogen.com Phone: 310 542 8822 Fax: 310 542 8028 Email: info@celprogen.com stemcells@celprogen.com

DATA SHEET

Mouse Retinal Stem Cell Complete Differentiation Media with Serum

Catalog number: M66105-32DS

Description: Mouse Retinal Stem Cell Complete Differentiation Media with fetal bovine Serum and Antibiotics: total volume of 500ml. This media is available in ready to use format. Media is optimized for use with Mouse Retinal Stem Cell . Each lot is qualified to support Mouse Retinal Stem Cell differentiation, growth, and proliferation. Formulated for use in a 5% CO₂ humidified incubator, and use in low O₂ environment: 5% O₂ in a controlled 37^oC humidified incubator.

Storage Conditions: 2-8 ⁰C

Tissue culture: The Mouse Retinal Stem Cell Complete Differentiation Media with Serum was tested and Optimized with Mouse Retinal Stem Cell in tissue culture. The Mouse Retinal Stem Cell was grown on Celprogen's Mouse Retinal Stem Cell Extra-cellular Matrix pre-coated flasks and Mouse Retinal Stem Cell Complete Differentiation Media with Serum. This media was optimized on Extra-cellular Matrix and Media combination Cell Culture system for up to 12 passages within 80-120 population doublings. This media is formulated to be utilized in a 37^oC humidified incubator with 5% CO₂. (Protocol is provided with the purchase of product)

- **Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).
- **Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices &

Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* The product, Mouse Retinal Stem Cell Complete Differentiation Media with Serum, is established and manufactured by CELPROGEN Inc., for Research use only. This product is not for re-sale or may not be transferred to a third party prior to written request and approval by CELPROGEN Inc.