

3914 Del Amo Blvd Suite 901 Torrance, CA 90503 www.celprogen.com

Phone: 310 542 8822 Fax: 310 542 8028 Email: info@celprogen.com stemcells@celprogen.com

DATA SHEET

Product name: Rat Neonatal Oligodendrocyte -Plated T75

Catalog number: 12090-03-T75

Description:

Rat Neonatal Oligodendrocyte T75 Plated Cells. Also available in T25, T150, and T225 tissue culture flasks with plated Cells, shipped at room temperature or Frozen vial shipped in Dry-ice. They were maintained in Celprogen's Rat Neonatal Oligodendrocyte Complete Growth Medium and sub-cultured every 24 to 48 hours on Rat Neonatal Oligodendrocyte Extra-cellular Matrix.

Source: Rat Neonatal

Negative-PCR and mycoplasma agar methods Mycoplasma test: Sterility: Negative for bacteria, yeast, and mold

Storage Conditions:

Liquid nitrogen vapor phase for frozen Ampule of Rat Neonatal Oligodendrocyte. For plated cells in tissue culture flask, upon receipt of the cells wipe the flask with 70% ethanol and transfer to sterile tissue culture hood. In the tissue culture hood remove the media of the cells and wash the cells with 1X PBS sterile solution, for 2-3 minutes, remove the PBS solution and then trypsinize. After trypsinization of the Cells neutralize the trypsin with equal volume of Rat Neonatal Oligodendrocyte Complete Growth Media with serum and collect the Cell suspension in sterile conical centrifuge tube in the tissue culture hood. Centrifuge the cell suspension at 100g for 7 minutes in centrifuge. Plate cells 5×10^5 cells per pre-coated flasks with Rat Neonatal Oligodendrocyte Extra-cellular Matrix for Expansion in Rat Neonatal Oligodendrocyte Complete Growth Medium.

Tests Performed:

Microbial: Negative Incorporation of acetylated Low Density Lipoprotein (LDL): Positive

Positive Markers:

Vimentin, Sox1, O4

Morphology

& Proliferation:

Mixed population of cells with approximately 95% attached cells and the other 5.0% in suspension, need to change cell culture media every day after 48 hours of initial cell culture or when the media starts changing color to slight yellow for pink. Fast growing cell culture. Change media with Celprogen's Rat Neonatal Oligodendrocyte Complete Growth Medium with the appropriate Rat Neonatal Oligodendrocyte



Extra-cellular Matrix. Temperature 37^{0} C in 5% CO₂ humidified incubator.



3914 Del Amo Blvd Suite 901 Torrance, CA 90503 www.celprogen.com Phone: 310 542 8822 Fax: 310 542 8028 Email: info@celprogen.com stemcells@celprogen.com

Sub-culturing:

- 1. That the vial with gentle agitation in a $37^{\circ}C$ water bath or a dry $37^{\circ}C$ shaking incubator. For water bath that the optimized of the water.
- 2. Remove the thawed vial and wipe with 70% ethanol. Then transfer to the tissue culture hood.
- **3.** Transfer the vial contents to a sterile centrifuge tube, and gently add pre-warmed Rat Neonatal Oligodendrocyte Growth Media to the centrifuge tube. Use additional Rat Neonatal Oligodendrocyte Complete Media to rinse the vial and transfer the liquid to the centrifuge tube repeat this once more to ensure you have all the cells transferred to the 15ml centrifuge tube. Centrifuge the cells at

1500 RMP for 5 minutes. Remove the supernatant and re-suspend the cell pellet in 500ul of Rat Neonatal Oligodendrocyte Complete Growth Medium.

- **4.** Add the 500ul of cells to appropriately sized flask pre-coated with Rat Neonatal Oligodendrocyte Extra-cellular Matrix with 15ml of Rat Neonatal Oligodendrocyte Complete Growth Medium.
- 5. Incubate the cells in the T75 flask in a 37^oC in 5% CO₂ humidified incubator. Perform 100% Media Change every 24 to 48 hours.
- 6. Medium renewal every other or 2-3 days, sub-culturing ratio: 1:3

Freezing Medium:	Available for purchase Cat# M12090-03FM
Trypsin:	Available for purchase Cat# T1509-014
IX PBS:	Available for Purchase Cat# P1408-013
Storage temperature:	Liquid nitrogen vapor phase

Product Orders:

Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits/Forms:

In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Biosaftely Level 1

Notices & Disclaimers:

CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.

This Product Rat Neonatal Oligodendrocyte established and manufactured by CELPROGEN Inc., and

is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.