



3914 DEL AMO BLVD  
SUITE 901  
TORRANCE CA 90503  
www.celprogen.com

Phone: 310 542 8822  
Fax: 310 542 8028  
Email: info@celprogen.com  
stemcells@celprogen.com

## DATA SHEET

<b>Product name:</b>	Human Vaginal Epithelial Cell Culture Media Phenol and Serum Free Media
----------------------	---

**Catalog number:** M36078-11PN

**Description:** Human Vaginal Epithelial Cell Culture Media Phenol and Serum Free Media with Antibiotics: total volume of 500ml. This media is available in ready to use format. Media is optimized for use with Human Vaginal Epithelial Cell Culture Media . Each lot is qualified to support Human Vaginal Epithelial Cell Culture Media Growth and Proliferation. Formulated for use in a 5% CO<sub>2</sub> humidified incubator, and use in low O<sub>2</sub> environment: 5% O<sub>2</sub> in a controlled 37<sup>0</sup>C humidified incubator. Used for cell based assays, so phenol reagent would not interfere with endpoint colormatric readings.

**Storage Conditions:** 2- 8 <sup>0</sup>C

**Tissue culture:** The Human Vaginal Epithelial Cell Culture Media Phenol and Serum Free Media was tested and optimized with Human Vaginal Epithelial Cell Culture Media in tissue culture. The Human Vaginal Epithelial Cell Culture Media was grown on Celprogen's Human Vaginal Epithelial Cell Culture Media Extra-cellular Matrix pre-coated flasks. This media was optimized on this Extra-cellular Matrix and Media combination Cell Culture system up to 12 passages within 80-120 population doublings in a 37<sup>0</sup>C humidified incubator with 5% CO<sub>2</sub>. Note: Prior to using serum free media, the cells have to be weaned off their serum conditions. (Protocol is provided with the purchase of product)

**Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* The Product, Human Vaginal Epithelial Cell Culture Media Phenol and Serum Free Media, is established and manufactured by CELPROGEN Inc., for Research use only. This product



3914 DEL AMO BLVD  
SUITE 901  
TORRANCE CA 90503  
www.celprogen.com

Phone: 310 542 8822  
Fax: 310 542 8028  
Email: info@celprogen.com  
stemcells@celprogen.com

is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.