

# Activation of Nuclear Factor (Erythroid-2 Like) Factor 2 by Toxic Bile Acids Provokes Adaptive Defense Responses to Enhance Cell Survival at the Emergence of Oxidative Stress

Kah Poh Tan, Mingdong Yang, and Shinya Ito

*Division of Clinical Pharmacology and Toxicology, Department of Pediatrics (S.I.), Physiology and Experimental Medicine Program, Research Institute (K.P.T., M.Y., S.I.), Hospital for Sick Children; and Department of Pharmacology, Faculty of Medicine, University of Toronto, Ontario, Canada (K.P.T., S.I.)*

Received June 26, 2007; accepted August 27, 2007

## ABSTRACT

Oxidative stress, causing necrotic and apoptotic cell death, is associated with bile acid toxicity. Using liver (HepG2, Hepa1c1c7, and primary human hepatocytes) and intestinal (C2bbe1, a Caco-2 subclone) cells, we demonstrated that toxic bile acids, such as lithocholic acid (LCA) and chenodeoxycholic acid, induced the nuclear factor (erythroid-2 like) factor 2 (Nrf2) target genes, especially the rate-limiting enzyme in glutathione (GSH) biosynthesis [glutamate cysteine ligase modulatory subunit (GCLM) and glutamate cysteine ligase catalytic subunit (GCLC)] and thioredoxin reductase 1. Nrf2 activation and induction of Nrf2 target genes were also evident *in vivo* in the liver of CD-1 mice treated 7 to 8 h or 4 days with LCA. Silencing of Nrf2 via small-interfering RNA suppressed basal and bile acid-induced mRNA levels of the above-mentioned genes. Consistent with this, overexpression of Nrf2 enhanced, but dominant-negative Nrf2 attenuated, Nrf2 target gene induction by bile acids. The activation of

Nrf2-antioxidant responsive element (ARE) transcription machinery by bile acids was confirmed by increased nuclear accumulation of Nrf2, enhanced ARE-reporter activity, and increased Nrf2 binding to ARE. It is noteworthy that Nrf2 silencing increased cell susceptibility to LCA toxicity, as evidenced by reduced cell viability and increased necrosis and apoptosis. Concomitant with GCLC/GCLM induction, cellular GSH was significantly increased in bile acid-treated cells. Cotreatment with *N*-acetyl-L-cysteine, a GSH precursor, ameliorated LCA toxicity, whereas cotreatment with buthionine sulfoximine, a GSH synthesis blocker, exacerbated it. In summary, this study provides molecular evidence linking bile acid toxicity to oxidative stress. Nrf2 is centrally involved in counteracting such oxidative stress by enhancing adaptive antioxidative response, particularly GSH biosynthesis, and hence cell survival.

Exposure to excessive bile acids is toxic to the cells, contributing an etiopathological factor to a number of liver and intestinal diseases such as cholestasis and colorectal cancer (Rao et al., 2001; Debruyne et al., 2002). Among the bile acids, litho-

cholic acid (LCA), a hydrophobic secondary bile acid produced by colonic microflora on chenodeoxycholic acid (CDCA), is the most toxic bile acid, with genotoxic and mutagenesis-enhancing properties (Kawalek et al., 1983; Kozoni et al., 2000). In rodents, it induces intrahepatic cholestasis-like hepatotoxicity (Staudinger et al., 2001), and it promotes chemical-induced colon carcinogenesis (Kozoni et al., 2000). CDCA, the most hydrophobic primary bile acid, is able to cause severe liver injury in species such as rhesus monkey, and it causes mild hepatotoxicity in humans; its chronic administration results in increased colonic production of LCA (Hofmann, 2004).

This work was supported by Canadian Institute for Health Research Grant MT13747 (to S.I.). K.P.T. was supported by the Research Training Competition studentship, Research Institute, Hospital for Sick Children and University of Toronto Fellowship.

Article, publication date, and citation information can be found at <http://molpharm.aspetjournals.org>.  
doi:10.1124/mol.107.039370.

**ABBREVIATIONS:** LCA, lithocholic acid; CDCA, chenodeoxycholic acid; FXR, farnesoid X receptor; VDR, vitamin D receptor (VDR); PXR, pregnane X receptor; ROS, reactive oxygen species; Nrf2, nuclear factor (erythroid-2 like) factor 2; ARE, antioxidant-responsive element; GCL, glutamate cysteine ligase; GCLC, glutamate cysteine ligase catalytic subunit; GCLM, glutamate cysteine ligase modulatory subunit; GSH, glutathione; MEM, minimal essential medium; FBS, fetal bovine serum; DMSO, dimethyl sulfoxide; TBL, total bilirubin; ALT, alanine aminotransferase; AST, aspartate aminotransferase; PCR, polymerase chain reaction; NQO1, NAD(P)H quinone oxidoreductase; TRx1, thioredoxin reductase 1; siRNA, small interfering RNA; siNrf2, small interfering against Nrf2; siCtr, small interfering RNA control; BHQ, *tert*-butylhydroquinone; NTCP, Na(+)-dependent taurocholate cotransporting polypeptide; GCA, glycocholic acid; GCDCA, glycochenodeoxycholic acid; ChIP, chromatin immunoprecipitation; bp, base pair(s); LDH, lactate dehydrogenase; MES, 2-(*N*-morpholino)ethanesulfonic acid; ANOVA, analysis of variance; DCA, deoxycholic acid; GST, glutathione transferase; UDCA, ursodeoxycholic acid; CA, cholic acid; BSO, buthionine sulfoximine; ABC, ATP-binding cassette; NAC, *N*-acetyl-L-cysteine.

The integrity and coordination of efficient hepatic bile flow and intestinal bile extraction are hence critical for species survival. The liver and intestinal cells achieve this through a concerted network involving the nuclear transcription factors, such as farnesoid X receptor (FXR), vitamin D receptor (VDR), retinoid X receptor, liver X receptor, pregnane X receptor (PXR), and/or constitutive androstane receptor. These receptors regulate bile-metabolizing and -conjugation enzymes and bile transporters to prevent excessive accumulation of bile acids (Eloranta and Kullak-Ublick, 2005). Bile acids are regarded as signaling molecules that facilitate synchronization of the above-mentioned regulators in their handling of cellular bile fate. The cross-talk among these receptors is important in maintaining homeostasis of physiological bile extraction, constituting the baseline protection against bile acid toxicity.

Meanwhile, increased cellular production of reactive oxygen species (ROS) and oxidative stress has been implicated in exposure to toxicological concentrations of bile acids. Bile acid-induced oxidative stress results from induction of membrane permeability transition consequent to mitochondrial toxicity and activation of death receptors (CD95), which subsequently lead to apoptosis, via activation of proapoptotic effectors caspases, and necrosis (Palmeira and Rolo, 2004). Whether there exists any regulator to counteract such oxidative stress and the progression of bile acid toxicity is presently unknown. Because of its important role as an oxidative stress sensor and antiapoptosis factor, we hypothesized that the nuclear factor (erythroid 2-like) factor 2 (Nrf2) may play a central role by enhancing adaptive response and cell survival during exposure to excess bile acids.

Nrf2, a basic leucine zipper transcription factor that binds to antioxidant responsive element (ARE), is a chief regulator for many antioxidative, cytoprotective genes (Kensler et al., 2007). Among Nrf2 target genes, the glutamate cysteine ligase (GCL), composed of modulatory (GCLM) and catalytic (GCLC) subunits, is the rate-limiting enzyme for cellular biosynthesis of glutathione (GSH), an important intracellular antioxidant in preserving redox balances. Emerging studies have shown that Nrf2 is a multiorgan protector against various toxic reactive insults; among others are chemical carcinogens (Ramos-Gomez et al., 2001) and acetaminophen (Chan et al., 2001). Hence, robust Nrf2 activation in the cell may be a critical adaptive response to overcome oxidative stress-induced disease processes. However, Nrf2 activation is not merely a cellular response to all circumstances of oxidative stress, because exposure to some oxidative stress inducers such as high-dose UVB ray would in turn result in Nrf2 deactivation (Kannan and Jaiswal, 2006). Presently, it is not known whether toxic bile acids could activate Nrf2.

In this study, we combined *in vitro* and *in vivo* approaches to demonstrate that Nrf2 is activated by cytotoxic bile acids, thereby inducing genes, such as GCL and hence GSH biosynthesis, to protect the cells against bile acid toxicity.

## Materials and Methods

**Cell Culture and Chemicals.** The human hepatoma-derived HepG2 (American Type Culture Collection, Manassas, VA) and mouse hepatoma-derived Hepa1c1c7 (a gift from Dr. Patricia Harper, The Hospital for Sick Children, Toronto, ON, Canada) were maintained in  $\alpha$ -MEM with 10% fetal bovine serum (FBS). C2bbe1,

a subclone of colon carcinoma Caco-2 that displays a more homogeneous brush-border epithelial-like morphology (American Type Culture Collection), was maintained in Dulbecco's modified Eagle's medium supplemented with 10% FBS, 1.5 g/l sodium bicarbonate, and 10 mg/l holo-transferrin. The human primary hepatocytes were purchased from Celprogen (San Pedro, CA), and they were grown in specially formulated serum-free growth media (Celprogen). Experiments of all secondary cell lines were conducted within 10 cell passages. Treatments were given at  $\sim$ 80% confluence for all cell lines except for C2bbe1. Because C2bbe1 cells differentiate to mature colonocytes at confluence, treatments to this cell line were given 2 to 3 days after confluence. All chemicals were purchased from Sigma-Aldrich (St. Louis, MO), unless otherwise indicated. Test bile acids were dissolved in dimethyl sulfoxide (DMSO) [0.2% (v/v)]. Oligonucleotides were synthesized at the Toronto Centre for Applied Genomics (Toronto, ON, Canada) or Integrated DNA Technologies, Inc. (Coralville, IA).

**In Vivo Mouse Experiments.** The animal care and experimental procedures were approved by the Animal Care Committee at the University of Toronto and the Hospital for Sick Children (Toronto, ON, Canada). To examine whether Nrf2 target genes may have been modulated during acute exposure to LCA before the onset of symptomatic liver injury, 9- to 12-week-old CD-1 mice (Charles River Canada, Montreal, QC, Canada) were injected *i.p.* with a single dose of LCA at 125 mg/kg b.wt. dissolved in sterilized DMSO (final amount,  $<$ 1% b.wt.). Mice were killed 7 to 8 h after the treatment. In a separate experiment aiming to investigate changes in similar genes upon extended treatment with LCA, mice were injected the same dose of LCA dissolved in sterilized corn oil (final amount  $\sim$ 2% b.wt.) twice daily for 4 days. This treatment protocol has been used in the past to induce cholestatic liver injury in mice (Staudinger et al., 2001). Mice were killed 16 h after the last dosing. The use of corn oil as a solvent for LCA in the extended treatment protocol was to avoid the possible toxicity with chronic exposure to DMSO. At necropsy, portions of their livers were sampled in RNAlater reagent (Invitrogen, Carlsbad, CA) and neutral-buffered 10% formalin for mRNA and histological analyses, respectively. Nuclear protein extraction of chilled livers was carried out using Nuclear Extraction kit (Panomics, Fremont, CA). Sera of mice were collected for analysis of liver function/injury markers: total bilirubin (TBL), alanine aminotransferase (ALT), aspartate aminotransferase (AST), and  $\gamma$ -glutamyl transferase using established automated methods (Department of Pediatric Laboratory Medicine, The Hospital for Sick Children).

**cDNA Synthesis and Quantitative Reverse-Transcription Polymerase Chain Reaction.** Total RNA was isolated with RNeasy kit (QIAGEN, Valencia, CA) and reverse-transcribed into cDNA using random hexamers and Moloney murine leukemia virus or SuperScript II reverse transcriptase (Invitrogen). Aliquots of cDNA equivalent to 100 ng of RNA were used for real-time PCR performed on Applied Biosystems (Foster City, CA) 7500 Real-Time PCR system or Prism 7700 Sequence Detection system with reaction mode set at 50°C for 2 min, 95°C for 20 s, followed by 40 cycles of 95°C for 15 s and 56 or 60°C for 1 min. The primers for ribosomal 18S,  $\beta$ -actin, tata-box binding protein, glyceraldehyde-3-phosphate dehydrogenase, GCLM, GCLC, and NAD(P)H quinone oxidoreductase 1 (NQO1) were purchased from predesigned and preoptimized Taqman primer-probe sets (Assay-on-Demand Gene Expression probe; Applied Biosystems), whereas custom-made primers for SYBR Green real-time PCR detection were used for the other gene transcripts (primer sequences available upon request). To ensure specificity, primer pairs were designed to span across two neighboring exons and detection of a single peak in dissociation curve analysis. The  $\Delta\Delta C_t$  method (Livak and Schmittgen, 2001) was used to quantify the amplification-fold difference between treatment and vehicle-treated control groups, with the  $C_t$  value of target genes being adjusted to individual housekeeping gene (glyceraldehyde-3-phosphate dehydrogenase,  $\beta$ -actin, tata-box binding protein, and/or 18S),

whichever expression was not affected by treatment protocols. Measurements were done in duplicate or triplicate with variability  $<0.5 C_t$ .

**Immunoblotting.** Whole cell lysate was prepared in radioimmunoprecipitation assay buffer with protease inhibitor cocktail (Roche Applied Science, Indianapolis, IN). Ten to 30  $\mu\text{g}$  of protein was dissolved in 4 to 12% bis-tris gel (NuPage Novex gel system; Invitrogen) and then transferred onto a nitrocellulose membrane (GE Healthcare, Chalfont St. Giles, Buckinghamshire, UK). Primary antibodies (working concentration) used were as follows: rabbit polyclonal anti-GCLC antibody-1 (1:2000) (NeoMarkers, Fremont, CA), rabbit antiserum against GCLM (1:3000) (custom-made; Alpha Diagnostics, San Antonio, TX; see below), rabbit polyclonal anti-Nrf2 c-20 (1:750) (Santa Cruz Biotechnology, Inc., Santa Cruz, CA), rabbit polyclonal anti-thioredoxin reductase 1 (TRx1) (1:3000) (Abcam Inc., Cambridge, MA), mouse monoclonal anti- $\beta$ -actin (1:10,000) (Sigma-Aldrich), and goat polyclonal anti-lamin B c-20 (1:200) (Santa Cruz Biotechnology, Inc.). Based on analyses of hydrophilicity, antigenicity, accessibility, and sequence homology with other related proteins, an antiserum against a peptide (amino acids 29–45) of human GCLM was raised in rabbits. The immunogenicity and specificity were checked by enzyme-linked immunosorbent assay, and its ability to detect an  $\sim 28$ -kDa protein (predicted size of GCLM) with reactivity halted after preabsorption of the antibody in excess immunogen. To ensure equal loading for whole cell lysate and nuclear protein,  $\beta$ -actin and lamin B, respectively, were probed on the same stripped blot membranes after being used for detecting target proteins.

**RNA Interference.** A combo of four gene-specific small-interfering RNAs (siRNA) against human Nrf2 (NM\_006164) was used (Dharmacon SMARTpool siRNA reagent; Dharmacon RNA Technologies, Lafayette, CO). Overnight-seeded HepG2 and C2bBE1 cells at  $\sim 40$  and  $\sim 60\%$  confluence, respectively, were transfected for 48 h with 50 nM siRNA against Nrf2 (siNrf2) or equal molar mismatched siRNA controls (siCtr). These siRNAs were earlier complexed with liposome carrier Dharmafect I (Dharmacon RNA Technologies) at 0.2  $\mu\text{M}$ /nM siRNA concentration in serum-free Opti-MEM (Invitrogen). Under this condition, the transfected cells after 48 h looked normal morphologically, and they did not differ from untransfected cells in cell viability and mRNA levels of inflammatory marker interleukin-6 (data not shown). Treatments with bile acids were then followed for 16 to 18 h. To ensure achieving functional and specific silencing, the mRNA levels of Nrf2, known Nrf2 target genes, and homologous subtypes Nrf1 and Nrf3 were compared between siNrf2 and siCtr groups before and after treatments in all experiments.

**Plasmid Constructs.** The expression vectors for Nrf2 (pEF\_Nrf2), dominant-negative Nrf2 (pEF\_DNrf2), and empty vector (pEF) were kindly provided by Dr. Jawed Alam (Ochsner Clinic Foundation, New Orleans, LA). To make an ARE-reporter construct (pGL3\_ARE), a DNA duplex (CGGGGTACCGCCGCACAAAGCGCTGAGTCACGGGGAGGCAGATCTCC) (core ARE is underlined;  $-3595/-3625$  of hGCLC gene) containing the indispensable ARE motif of GCLC ( $-3604/-3614$ ) (Mulcahy et al., 1997) with KpnI/BglII at 5' and 3' ends, respectively, was constructed by annealing two polyacrylamide gel electrophoresis-purified complementary oligonucleotides. This insert was ligated to similar restriction enzyme sites of pGL3 luciferase reporter plasmid with simian virus 40 promoter (Promega, Madison, WI). A similar reporter construct has been successfully used previously (Mulcahy et al., 1997). The responsiveness and robustness of our ARE reporter to Nrf2 transactivation was confirmed by testing of a panel of Nrf2 activators such as *tert*-butylhydroquinone (BHQ), lipoic acid, and diethyl maleate (data not shown), as well as cotransfection with Nrf2 and dominant-negative Nrf2 expression vectors. To ensure specificity, a mutant ARE reporter construct (pGL3\_mARE) introducing three point mutations on ARE was cloned by PCR-mediated site-directed mutagenesis using Pfu Turbo DNA polymerase (Stratagene, La Jolla, CA) with complementary primers 5'-AGCtaTGAGgCACGGGGAGGCAG-3' (underlined sequence is core ARE; lowercases are mutated nucleotides) on

the template pGL3\_ARE. The PCR condition was 95°C for 5 min followed by 22 cycles of 95°C for 15 s, 55°C for 30 s, 68°C for 10 min, and a final extension of 68°C for 10 min. The template was then digested by DpnI, and mutant clones were transformed in XL-1 blue competent cells (Stratagene). Successful insertion and mutation introduction were confirmed by sequencing. The cDNA clone of human Na(+)-taurocholate cotransporting polypeptide (NTCP) (Origene Technologies, Rockville, MD) was subcloned into the NotI site of pTarget expression vector (Promega), and it was stably transfected into HepG2. Stable clones transfected with NTCP or empty vector (pTarget) were selected using 500  $\mu\text{g}/\text{ml}$  G-418-supplemented growth media.

**Transfection, Reporter Assays, and Overexpression Studies.** HepG2 cells at  $\sim 50\%$  confluence were transiently transfected overnight with 0.1  $\mu\text{g}$  of the firefly luciferase reporter pGL3\_ARE or pGL3\_mARE, 0.02  $\mu\text{g}$  of the *Renilla reniformis* luciferase control reporter pRL-TK with or without cotransfection with 0.2  $\mu\text{g}$  of expression vectors using Lipofectamine 2000 (Invitrogen) as transfection carrier. Treatments with bile acids were then carried out for another 16 to 18 h in all experiments, unless otherwise stated. Conjugated bile acid treatments [glycocholic acid (GCA) and glychenodeoxycholic acid (GCDCA)] were done on NTCP-transfected HepG2. Luciferase activities of the cell extracts were determined with the Dual-Luciferase Reporter Assay system (Promega). Relative luciferase activity (relative light unit) was calculated from firefly luciferase values normalized to those of the control *R. reniformis* luciferase, and activity is expressed as ratios to vehicle-treated empty pGL3 promoter construct, and, if any, cotransfected expression vector. All experiments were done in triplicate, and they were repeated at least twice. For overexpression studies, Hepa1c1c7 cells at 50% confluence in T25 flasks were transfected with 3  $\mu\text{g}$  of Nrf2 or dominant-negative Nrf2 expression vector for 24 h, followed by treatments with bile acids for another 20 to 22 h.

**Quantitative Chromatin Immunoprecipitation.** The assay was performed using the ChIP assay kit (Upstate, Charlottesville, VA) with slight modifications. After 6 h of treatments, chromatin protein-DNA of HepG2 cells was fixed (cross-linked) in neutral-buffered 1% formaldehyde at room temperature for 10 min. Further fixation was stopped by 125 mM glycine buffer. The DNA was sheared by sonication on ice into fragments of  $\sim 500$  bp. An aliquot (one fourth) of sample supernatant was saved as input DNA for later PCR analysis. After preclearing with protein A agarose beads, supernatants were incubated with a ChIP-graded anti-Nrf2 antibody (1:250; Santa Cruz Biotechnology, Inc.) in rotation at 4°C overnight. To control for nonspecific binding of antibody used, an equal amount of the host antibody against an irrelevant protein (rabbit polyclonal anti-CYP1A1) from the same manufacturer (Santa Cruz Biotechnology, Inc.) was included in a separate batch of control supernatants and followed through the remaining protocols. Antibody-chromatin complexes were collected by salmon sperm DNA/protein A beads. DNA was released from cross-linked complexes with proteinase K at 65°C for 4 h followed by 72°C for 10 min. DNA was then extracted and eluted with 120  $\mu\text{l}$  of Tris, pH 8.0, buffer using the DNeasy kit (QIAGEN), and the contaminant RNA was cleaved with RNase A (Invitrogen). For detection of the ARE of GCLM ( $-56/-66$ ) (Erickson et al., 2002) and of GCLC ( $-3604/-3614$ ) (Mulcahy et al., 1997) by real-time PCR, the primer sets and Taqman probe (5'-5-carboxyfluorescein, 3'-5-carboxytetramethylrhodamine) were designed by PrimerQuest software (Integrated DNA Technologies, Inc.), which amplify 5'-region exactly on the core ARE. The primers for detecting the ARE of GCLM were as follows: sense, 5'-CGCGGGATGAGTACGGT-3'; antisense, 5'-GGGAGAGCTGATTCCAAACTGA-3'; and probe, 5'-ACGAAGCACTTCTCGGCTACGAT-3', which amplify a 79-bp product ( $-33/-112$ ). For probing the ARE of GCLC, the primers used were sense, 5'-GGATGAGACTTTGCCCTAAGAAG-3'; antisense, 5'-GCGCAGTTGTGTGATACAG-3'; and probe, 5'-CGCACAAGCGCTGAGTCAC-3', which amplify a 160-bp product ( $-3479/-3609$ ). Quantitation of NRF2 bound to these AREs after the treatments was carried out on 5% of DNA eluates with qPCR anal-

ysis similar to that for the mRNA, except that the  $C_t$  value of amplicon from each sample's input DNA was used as normalization control as described previously (Beresford and Boss, 2001).

**Cytotoxicity, Necrosis, and Apoptosis.** For cytotoxicity analysis, a nontoxic assay, namely, Alamar Blue (BioSource, Nivelles, Belgium), was used. This assay uses a fluorometric indicator to measure the chemical reduction of cell medium, which correlates directly to the metabolic activity of viable cells. The working assay medium [10% (v/v) Alamar Blue in  $\alpha$ -MEM, 2% FBS, and 1% penicillin/streptomycin, 37°C] was first incubated with cells seeded on 24-well culture plate before treatment to obtain baseline/pre-treatment values. The measurement was made at excitation/emission/cut-off  $\lambda = 540/590/570$  nm after 1 h of incubation with the assay medium at 37°C. Immediately after the measurement, the cells were rinsed with prewarmed PBS followed by the treatments. At various time points, treatment media were replaced with fresh assay media to allow for a continuous monitoring of cell viability. The fluorescent unit of each treatment and control was expressed as percentage of change relative to individual baseline/pre-treatment value.

To determine necrosis, cellular release of lactate dehydrogenase (LDH) into treatment media was measured with an LDH detection kit (Roche Applied Science). To control for cell mass and spontaneous release of LDH by viable cells into media, the ratio of LDH activity in the medium to the cells (cell lysate) was determined and then subtracted from those of the vehicle-treated controls. The measurement was made colorimetrically at  $\lambda = 490$  nm. The intra-assay variability of duplicate determinations was 2.2.

To assess apoptosis, the cellular caspases activity was measured using the rhodamine 110-conjugated substrate *N*-benzyloxycarbonyl-Asp-Glu-Val-Asp (Invitrogen). Although traditionally known to detect caspase-3 activity, recent analysis by the manufacturer showed that this substrate is also a target of multiple caspases such as 6, 7, 8, and 10. The caspases activity of cell lysate was quantitated at excitation/emission  $\lambda = 496/520$  nm, and it was normalized to individual protein concentration measured by Bio-Rad protein assay (Bio-Rad Laboratories, Hercules, CA).

**Total GSH Quantitation.** Cellular GSH was quantitated with a GSH assay kit (Cayman Chemical, Ann Arbor, MI) based on an established GSH recycling enzymatic method (Tietze, 1969). After 24-h treatment with bile acids, HepG2 cells were lysed in ice-cold MES buffer after a quick freeze-thaw cycle and then deproteinized by 0.5 g/ml metaphosphoric acid. An aliquot of each sample was saved before deproteinization for determining protein content. The total GSH of deproteinized cell supernatants was measured against an oxidized glutathione standard curve according to the manufacturer's instruction. The measurement unit was expressed as nanomoles per milligram of protein per minute. The intra-assay variability of duplicate determinations for all samples repeated in four experiments was 3.1.

**Statistical Analysis.** Statistical tests were conducted using SigmaStat 3.1 (Systat Software, Inc., Point Richmond, CA) or SPSS 10.1 (SPSS Inc., Chicago, IL). Normality and equal variance tests were first carried out to guide subsequent statistical analyses. Multiple group comparisons were carried out by one-way ANOVA (parametric) or one-way ANOVA on ranks (nonparametric). Once statistical significance was attained ( $p < 0.05$ ), the Dunnett's (parametric) or Dunn's (nonparametric) test comparing between treatment and control groups was initiated. Comparisons between two groups on single variable were accomplished by Student's independent *t* test (parametric) or Mann-Whitney U test (nonparametric). Difference with  $p < 0.05$  was considered statistically significant.

## Results

**Induction of Nrf2 Genes by Bile Acids in Liver and Intestinal Cells.** A dose-response increase in mRNA of GCLM and GCLC after LCA and/or CDCA treatment was

noted for HepG2 and C2bbe1 cells, with a significant >4-fold induction at  $\geq 50$   $\mu$ M LCA and  $\geq 100$   $\mu$ M CDCA (Fig. 1a). Peak inductions of GCL subunit genes occurred at 50 to 75  $\mu$ M LCA and at 100 to 150  $\mu$ M CDCA. Increased bile acid concentrations, i.e., LCA  $\geq 100$   $\mu$ M or CDCA  $\geq 250$   $\mu$ M, resulted in markedly cell death and reduced inductions of GCL genes at 24 h of treatment (data not shown). Significantly higher GCL gene transcripts, although at a lower magnitude (2–4-fold induction), were also noticed for the primary human hepatocytes treated with both bile acids (Fig. 1a). So were similar treatments given to the mouse hepatoma Hepa1c1c7 (data not shown). In C2bbe1 cells, CDCA treatment (100  $\mu$ M) caused a modest increase in GCL genes ( $\sim 2$ -fold), whereas treatment with deoxycholic acid (DCA) ( $\geq 100$   $\mu$ M), a secondary bile acid often associated with toxicity and carcinogenesis in colonic cells, resulted in comparable GCL inductions to those of LCA treatment (data not shown).

Furthermore, mRNA of Nrf2 and a panel of known Nrf2 target genes such as NQO1, TRx1, ferritin light subunit, and heme oxygenase I were also simultaneously induced  $\geq 2$ -fold by bile acids in all test cells (Fig. 1b). Particularly, TRx1, an important seleno-enzyme in cellular thiol and redox maintenance, was increased >4-fold in HepG2 and C2bbe1. Note that glutathione transferase P1, which was induced by bile acids in HepG2 and C2bbe1, was not evident in primary hepatocytes (Fig. 1b). Instead, another GST subtype, glutathione transferase A, was increased by bile acids to  $\sim 2$ -fold in primary hepatocytes (data not shown). This disparity suggests possible cell type specificity in regulation of GSTs by bile acids.

Increased protein levels corresponding to mRNA induction were also noted (Fig. 1c). The apoptosis marker (caspases activity) and cell viability analyses showed that a mild toxicity began to occur in HepG2 cells at 60 to 80  $\mu$ M LCA treatment, followed by a precipitous increase in cell death and caspase activity at >80  $\mu$ M (Fig. 1d). It is noteworthy that the induction of GCL subunits and other antioxidative genes peaked in the range of LCA (60–80  $\mu$ M) during which HepG2 began to experience mild toxicity. These findings suggest that induction of the cytoprotective genes may represent an adaptive cell defense mechanism against LCA toxicity.

**In Vivo Activation of Nrf2 Target Genes.** Acute administration (7–8 h) of LCA to mice at a dose known to induce cholestatic liver injury (Staudinger et al., 2001) resulted in Nrf2 accumulation in the nuclei, a signature event of Nrf2 activation (Fig. 2a). This phenomenon coincided with significant inductions of Nrf2 target genes (Fig. 2b, top) found to be increased in the in vitro studies. It is noteworthy that the increase of TRx1 transcripts rose to  $\sim 50$ -fold at acute exposure to LCA, implying a possible critical role of this enzyme in early toxicity of LCA. At this shorter exposure to LCA, however, analysis of serum liver function and cholestatic markers (ALT, AST,  $\gamma$ -glutamyl transferase, and TBL) as well as liver histology did not indicate liver dysfunction or pathological changes (data not shown).


With prolonged LCA treatment during which symptomatic liver injury (elevated ALT, AST, and TBL, and liver necrosis in histological analysis; data not shown) already occurred, induction of Nrf2 target genes, such as GCL subunit gene transcripts, was found to sustain compared with those treated acutely with LCA (Fig. 2b, bottom). Nqo1 was in-

creased with prolonged treatment, whereas TRx1 induction was subdued. Consistent with the observations from primary human hepatocytes,  $\alpha$  class of mouse Gst (Gsta1/2), rather than glutathione transferase P1, was found to be induced by LCA, with  $\sim 10$ -fold induction at 4 days of treatment (Fig. 2b, bottom). Treatment with vehicle alone did not differ in mRNA of genes under study compared with untreated animals (data not shown).

**Involvement of Nrf2 and Activation of Nrf2-ARE Transcription Machinery.** To examine whether Nrf2 participated in the preceding gene activations, we silenced Nrf2 of HepG2 and C2bbe1 via siRNA. This resulted in significant reductions of  $>60\%$  in Nrf2 mRNA and protein without interfering with other homologous Nrf subtypes (Fig. 3, a and b). Nrf2 silencing significantly decreased the basal levels of known Nrf2 target genes (Fig. 3, c–e), an observation similar to that seen in *in vivo* Nrf2 knockout mice (Lee et al., 2005). In addition, the induction of GCLM, GCLC, and other Nrf2 target genes by bile acids in HepG2 (Fig. 3, c and e) and C2bbe1 (Fig. 3d) has been mitigated. Comparable reduction

in inducible expressions of Gclm occurred with transfection of dominant-negative Nrf2 in Hepa1c1c7 cells, consistent with the enhanced gene induction with Nrf2 overexpression (Fig. 3f). Similar observations were noted for other Nrf2 target genes such as Gclc and Nqo1 (data not shown).

To verify that there was an activation of Nrf2-ARE transcription machinery with exposure to toxic bile acids, we extracted the nuclear proteins of bile-acid-treated HepG2 over different time points across 24 h. Translocation of cytosolic Nrf2 to nucleus represents the prerequisite event of receptor activation. Nrf2 started to be enriched in cell nuclei within 1 to 3 h of bile acid treatments, and it was sustained through 24 h, with CDCA-treated cells showing reduced Nrf2 translocation events with longer time of exposure (24 h) (Fig. 4a). Furthermore, various bile acids were found to increase the activity of an ARE-reporter assay in a dose-dependent manner, suggesting that these bile acids were capable of inducing Nrf2 transactivation (Fig. 4b). The magnitude of luciferase activity of the highest concentration of test bile acids was compatible with those of treatments with antioxi-




**Fig. 1.** Induction of Nrf2 target genes by bile acids. **a**, mRNA levels of GCLM (white bar) and GCLC (black bar) in HepG2 and C2bbe1 cells (left) after 24-h treatment with increasing doses of LCA or CDCA: [LCA] = 6.25, 12.5, 25, 50, and 75  $\mu$ M and [CDCA] = 25, 50, 75, 100, and 150  $\mu$ M. Right, human primary hepatocytes (passages 4–6) were treated with 75  $\mu$ M LCA or 100  $\mu$ M CDCA for 24 h. *y*-axis, -fold change versus vehicle-treated controls. \*,  $p < 0.05$ , significantly different from vehicle-treated control by one-way ANOVA followed by post hoc test for HepG2 and C2bbe1 or \*\*,  $p < 0.01$  by *t* test for primary hepatocyte. Mean and S.E.M. ( $n = 3$ –6). **b**, mRNA levels of other known Nrf2 target genes after 24-h treatment with bile acids in HepG2, C2bbe1, and primary hepatocytes: [LCA] = 75  $\mu$ M and [CDCA] = 100  $\mu$ M. Mean and S.E.M. ( $n = 3$ –6). **c**, representative immunoblots of protein lysate (10  $\mu$ g for HepG2; 30  $\mu$ g for C2bbe1) probed for Nrf2 target genes after 24-h treatment with bile acids. [LCA] = 70  $\mu$ M and [CDCA] = 100  $\mu$ M. **d**, cell viability (by Alamar Blue) and apoptotic marker (caspases activity) measured across increasing concentrations of LCA treatment in HepG2. Viability was measured at 24 h, whereas caspases activity was measured at 6 h of LCA treatment. Note that the induction of GCL genes peaked at 60 to 80  $\mu$ M LCA (shown by closed inverted triangles) during which mild cellular toxicity began to occur. Representative results from four determinations are shown.

dants *tert*-BHQ (100  $\mu$ M) and  $\alpha$ -lipoic acid (200  $\mu$ M) (data not shown), denoting that bile acids are equally potent Nrf2 activators. Note that there was an  $\sim$ 8-fold increase in reporter activity with vehicle DMSO treatment compared with that of the empty vector harboring only the simian virus 40 promoter. This suggests the existence of a strong constitutive Nrf2-ARE transactivational activity in HepG2 cells, an observation in line with the persistent oxidative stress observed in many cancerous cell lines (Brown and Bicknell, 2001). HepG2 cells are deficient in conjugated bile acid transporters such as NTCP, which leads to its resistance to conjugated bile acid-induced oxidative stress (Kullak-Ublick et al., 1996). Transfection of NTCP expression vector hence restores, to some degree, its sensitiveness. In this study, we found that GCDCA, a known cholestatic conjugated bile acid, significantly induced the ARE reporter. This suggests that activation of Nrf2-ARE may be crucial to counteract the toxicity of GCDCA as reported previously (Dent et al., 2005). The potency of bile acids in activating this reporter based upon molarities was LCA > CDCA  $\approx$  DCA > GCDCA  $\geq$  ursode-

oxycholic acid (UDCA) > cholic acid (CA) > GCA. This order is in consensus with the toxicity profile of these bile acids, particularly in terms of their ability to produce oxidative stress (Krähenbühl et al., 1994). Overexpression of Nrf2 further enhanced the reporter activity by bile acids, whereas coexpression of dominant-negative Nrf2 attenuated the activity, and a mutant ARE construct was completely uninducible by bile acids (Fig. 4c). Using the quantitative ChIP assay, we found an increased Nrf2 occupancy to the AREs of both GCL subunits in the native cell context upon 6-h treatment with bile acids (Fig. 4d). Taken together, our data suggest that activation of the Nrf2-ARE machinery underlies induction of Nrf2-target genes by toxicological concentrations of bile acids.

**Protective Role of Nrf2 in Bile Acid Toxicity.** To directly investigate the role of Nrf2 in protection against toxic bile acids, we first silenced Nrf2 of HepG2 via RNA interference upon which the cells were subjected to toxic LCA challenges. Nrf2 knockdown increased cell susceptibility to toxic LCA, with a significantly decreased cell viability starting from 4 h of treatment (Fig. 5a). Without LCA challenge, Nrf2-knockdown cells did not differ in cell viability from those treated with siCtr (data not shown). Significant protective effects of Nrf2 against LCA toxicity was also observed in C2bbe1 cells, and in HepG2 with  $\geq$ 300  $\mu$ M CDCA (data not shown). To investigate which route of cell death was particularly involved in protection by Nrf2, established markers of necrosis and apoptosis were examined. LCA at 90  $\mu$ M was used because at this dose, both apoptosis and necrosis were found to simultaneously occur. Necrotic event, as determined by LDH released into the culture media, remained constantly higher in Nrf2 knockdown cells than those of siCtr starting from 2 h of LCA treatment (Fig. 5b). Nrf2 silencing alone did not affect the cellular release of LDH (data not shown). Likewise, in the assessment of apoptosis, Nrf2-knockdown cells exhibited much higher and prolonged elevation of caspases activity than did siRNA control-treated cells upon LCA treatment (Fig. 5c). Silencing of Nrf2 alone did not result in increased basal caspases activity.

**Role of GSH in Resisting LCA Toxicity.** The increase in GCLM and GCLC, the rate-limiting enzyme in GSH biosynthesis, observed in earlier experiments after LCA (75  $\mu$ M) or CDCA (100  $\mu$ M) treatment was accompanied by a significant increase by >4-fold in cellular GSH levels at 24 h (Fig. 6a). This increase was comparable with treatment with 200  $\mu$ M  $\alpha$ -lipoic acid, a GSH inducer. To determine whether the induced cellular GSH is a protective mechanism against toxic bile acid, we cotreated HepG2 with a toxic dose of LCA and a GSH biosynthesis blocker, buthionine sulfoximine (BSO), which inhibits the activity of GCL subunits and blocks GSH biosynthesis. BSO treatment together with toxic LCA decreased cell resistance toward LCA exposure with more apparent effects in late treatment ( $\sim$ 24 h) (Fig. 6b). Furthermore, depletion of cellular GSH by pretreatment with BSO before LCA challenge markedly lifted cell resistance with a drastic drop in cell viability within 4 h of treatment. Conversely, toxic LCA challenge in the presence of an antioxidant and GSH precursor *N*-acetyl-L-cysteine (NAC) was found to alleviate the toxicity. This suggests that the basal as well as inducible GSH are important determinants of cellular resistance to LCA. Consistent with these findings, NAC cotreatment significantly reduced the oxidative stress-responsive ARE-reporter activity by LCA, indicating an antioxi-




**Fig. 2.** Activation of Nrf2 in mice treated with LCA. a, representative immunoblots of nuclear fractions (30  $\mu$ g) probed against Nrf2 in the liver of mice after 7- to 8-h treatment with cholestatic LCA (125 mg/kg b.wt.). Lamin B was used as loading control for nuclear protein, whereas  $\beta$ -actin was probed to show unintentional contamination of cytosolic proteins in nuclear fraction preparation. Note that the increased nuclear Nrf2 cannot be explained by inclusion of contaminant cytosolic proteins. b, mRNA levels of Nrf2 target genes upon LCA treatment for 7 to 8 h (top) or for 4 days (bottom) in the liver of mice. Because there were differences in the basal gene expression of Nrf2 target genes between sexes, comparisons of all target genes between treated and untreated groups were adjusted for sex. Induction of antioxidative genes by LCA, however, occurred to both sexes. Significantly different from vehicle-treated controls by *t* test. \*,  $p < 0.05$ ; \*\*,  $p < 0.01$ . Mean and S.E.M. ( $n = 5$ –9 for 7- to 8-h treatment;  $n = 3$  for 4-day treatment).

tive effect (Fig. 6C). BSO cotreatment, in contrast, further increased the reporter activity, consistent with a heightened oxidative stress (Fig. 6c).

## Discussion

The discovery of bile acids as key signaling molecules in the enterohepatic circulation system reveals a critical role of hepatic and intestinal xenobiotic nuclear receptors in the metabolism and detoxification of bile acids (Chawla et al.,

2000). Particularly, the cytotoxic hydrophobic bile acids CDCA and LCA have been shown to be ligands and potent inducers of these receptors. LCA, at physiological and non-toxicological concentrations (5–30  $\mu\text{M}$ ), can activate FXR (Makishima et al., 1999, 2002) and VDR (Makishima et al., 2002), indicating their crucial role in physiological handling of this bile acid. The major detoxification routes of LCA, i.e., sulfation by sulfotransferase 2A and 7 $\alpha$ -hydroxylation by CYP3As, are coordinated by VDR (Makishima et al., 2002; Echchgadda et al., 2004). FXR, which induces the hepatic bile




**Fig. 3.** Involvement of Nrf2 in induction of glutamate cysteine ligase subunits by bile acids. a, mRNA levels of all Nrf subtypes after 72-h treatment with siRNA. y-axis, -fold change versus vehicle-treated cells transfected with siCtrl. \*,  $p < 0.001$ , significant difference between siRNA groups by  $t$  test. Mean and S.E.M. ( $n = 3-5$ ). b, representative immunoblots (20  $\mu\text{g}$  of protein lysate) of HepG2 after 48-h treatment with siRNA against Nrf2. c, basal (treated with vehicle DMSO) and inducible (treated with 70  $\mu\text{M}$  LCA or 100  $\mu\text{M}$  CDCA) gene transcripts of GCLM and GCLC in HepG2 (mean and S.E.M.;  $n = 3-5$ ). d, basal and inducible levels of GCLM and GCLC gene transcripts in C2bbe1 treated with vehicle or 70  $\mu\text{M}$  LCA. \*,  $p < 0.01$ , significant difference between siCtrl and siNrf2 with or without treatment by  $t$  test. Mean and S.E.M. ( $n = 3-5$ ). e, basal and inducible expressions of other Nrf2 target genes in HepG2. Mean  $\pm$  S.E.M. f, mRNA levels of GCLM induced by 70  $\mu\text{M}$  LCA or 100  $\mu\text{M}$  LCA in Hepa1c1c7 transfected with empty vector (pEF), Nrf2 (pEF\_Nrf2), or dominant-negative Nrf2 (pEF\_DNrf2) expression vector. y-axis, -fold change versus vehicle-treated cells transfected with empty vector pEF. \*,  $p < 0.05$ , significantly different from LCA-treated cells transfected with pEF control by  $t$  test. Mean and S.E.M. ( $n = 4$ ).

salt export pump BSEP (Ananthanarayanan et al., 2001) and down-regulates the bile-synthesizing enzyme CYP7A1 (Makishima et al., 1999), works to prevent intracellular accumulation of bile acids.

Interestingly, at higher and toxicological concentrations of LCA ( $\geq 50 \mu\text{M}$ ) and CDCA ( $\geq 100 \mu\text{M}$ ), which potentially cause cell injury, PXR (Staudinger et al., 2001; Makishima et al., 2002) and Nrf2, as shown in this study, are found to be activated. The activation of PXR and Nrf2 induces the major hydroxylation enzymes CYP3As and antioxidative genes (Eloranta and Kullak-Ublick, 2005; Kensler et al., 2007), which may represent an important adaptive mechanism of cellular defense against toxic bile acids. Furthermore, we observed that induction of multiple bile salt/conjugate efflux transporters such as ATP-binding cassette (ABC) transporters ABCC2, ABCC3, and ABCG2 by bile acids is dependent on Nrf2 (K. P. Tan, G. Woodland, M. Yang, K. Kosuge, M. Yamamoto, and S. Ito, unpublished data). Hence, the collec-

tive induction of cytoprotective genes by Nrf2 and PXR seems to set off a second line of protection against possible progression of bile acid toxicity toward irreversible cell death.

In this study, we showed for the first time that many bile acids, more potently LCA, CDCA, and DCA, are capable of activating redox-sensitive Nrf2. We also provided in vivo evidence that LCA is able to activate Nrf2, inducing similar target genes observed in in vitro studies. Because the induction of Nrf2 target genes by LCA in vivo was found to precede and sustain through biochemically and histologically overt liver injury, the collective induction of these antioxidative genes may be an integral part of cell defense against bile acid toxicity and hepatic injury. Previous studies have reported an increased intracellular production of detrimental hydroperoxides in isolated rat hepatocytes with hydrophobic bile acid exposure (Sokol et al., 1995), an observation in consensus with the increased oxidative stress byproducts in the liver of patients with cholestasis (Vendemiale et al.,



**Fig. 4.** Bile acids activate Nrf2 transcription machinery. **a**, representative immunoblots of nuclear fraction (10  $\mu\text{g}$ ) extracted from HepG2 treated with bile acids (LCA, 70  $\mu\text{M}$ ; CDCA, 100  $\mu\text{M}$ ) at indicated time points over 24 h. Lamin B was used as equal loading control for nuclear protein, whereas  $\beta$ -actin was probed to show possible contaminant cytosolic proteins in nuclear fraction preparation. **b**, ARE-reporter (luciferase) activity of HepG2 treated with increasing doses of various bile acids for 16 to 18 h. Abbreviations (doses): LCA (50, 70, and 90  $\mu\text{M}$ ), CDCA (50, 100, and 150  $\mu\text{M}$ ), CA (150, 200, and 400  $\mu\text{M}$ ), DCA (50, 100, and 150  $\mu\text{M}$ ), UDCA (50, 100, and 200  $\mu\text{M}$ ), GCA (200, 400, and 800  $\mu\text{M}$ ), and GCDCA (100, 200, and 400  $\mu\text{M}$ ). #, GCA and GCDCA were tested in NTCP-transfected HepG2 for 6 h. **y**-axis, -fold change in the ratio of luciferase activity (relative luciferase unit) (see *Materials and Methods* for detail) from those transfected with basic pGL3 promoter construct and treated with vehicle DMSO. \*,  $p < 0.05$ , significantly different from DMSO-treated pGL3\_ARE by one-way ANOVA followed by post hoc test. Mean and S.E.M. ( $n = 2-4$ ). **c**, ARE-reporter (pGL3\_ARE) activity with coexpression of Nrf2 or dominant-negative Nrf2, and mutant ARE-reporter (pGL3\_mARE) activity in HepG2 treated with bile acids (70  $\mu\text{M}$  LCA; 100  $\mu\text{M}$  CDCA). **y**-axis, -fold change in the ratio of luciferase activity (relative luciferase unit) from those transfected with basic pGL3 promoter construct and/or respective expression vector, and treated with DMSO. \*,  $p < 0.05$ , significant difference between vehicle control and bile acids by one-way ANOVA followed by post hoc test. Mean and S.E.M. ( $n = 3-6$ ). **d**, ChIP analysis examining Nrf2 occupancy to AREs of both GCLM and GCLC genes upon 6-h treatment with bile acids (70  $\mu\text{M}$  LCA; 100  $\mu\text{M}$  CDCA) in HepG2. BHQ (200  $\mu\text{M}$ ), known to transcriptionally activate GCLM and GCLC, was included as positive control. Negligible detection from samples incubated with host IgG (anti-CYP1A1) ruled out contribution of nonspecific binding from antibody. \*,  $p < 0.05$ , significantly different from controls ( $t = 0$ ) and vehicle treatment by  $t$  test. Mean and S.E.M. of triplicate determinations of representative experiments.

2002). Because Nrf2 activation is indicative of cellular anti-oxidative response, our study provides molecular evidence linking mechanism of bile acid toxicity to oxidative stress.


We further showed that induction of hepatic GCL subunits via Nrf2, which provokes GSH biosynthesis, can increase hepatocyte resistance and survival during excessive bile acid exposure. The essential role of GSH in hepatic protection against injury and oxidative xenobiotic insults has been well exemplified (Huang et al., 2001; Glosli et al., 2002). In agreement, *in vivo* knockout of Nrf2 enhances sensitivity of death receptor-induced hepatic apoptosis as a result of decreased

GSH levels (Morito et al., 2003). GSH is also known to protect against mitochondrial injury, a major mechanism of bile acid toxicity (Palmeira and Rolo, 2004). A fraction of cytosolic GSH that becomes mitochondrial GSH is crucial in the defense of oxidant-induced mitochondrial-mediated cell death (Fernandez-Checa and Kaplowitz, 2005). In addition, Nrf2 activation has been shown to protect mitochondria by preventing inhibition of mitochondrial complex II upon exposure to oxidative neurotoxins (Calkins et al., 2005). In intestinal mucosa, cellular GSH has an essential role in maintaining epithelial integrity, transport activity, and metabolism of and susceptibility to luminal toxins (Aw, 2005). Overall, our study, coupled with supportive evidence from recent literature, suggests that protection conferred by hepatic and intestinal Nrf2 against bile acid-induced oxidative stress is, at least partly, achieved by increasing GSH levels.


The simultaneous induction of other Nrf2 target genes may work in concert with GCL subunits to combat bile acid-induced oxidative stress and facilitate adaptive responses. Of particular importance is TRx1, an enzyme engaged in NADPH-dependent catalysis of various redox proteins (Rundlöf and Arner, 2004). It has been shown to act as a key adaptation-promoting mediator for prior exposure to 4-hydroxynonenal, a reactive lipid peroxidation-derived molecule, in inducing cellular tolerance to future oxidative stress attack (Chen et al., 2005). Indeed, intermediate cellular stress has recently been proposed to provide an adaptation advantage by invoking enhanced cellular survival/tolerance mechanisms (Schoemaker et al., 2003; Chen et al., 2005). Activation of nuclear factor  $\kappa$ B as well as Nrf2 has been shown to play an important role in this adaptation process. The drastic induction of TRx1 observed in mice upon acute exposure to toxic LCA in this study may indicate a critical role of this enzyme in adaptation process against LCA toxicity. To address whether and how this process is taking place, future studies are needed.

The precise mechanisms by which toxic bile acids activate Nrf2 remain a subject of future studies. Enormous production of ROS from mitochondrial stress has long been accounted for the main source of oxidative stress induced by bile acids (Palmeira and Rolo, 2004). Insurgence of these ROS potentially targets the cysteine oxidant-sensors of Keap1, an actin-anchored cytosolic sequester that facilitates Nrf2 degradation by ubiquitin-proteasome pathway, which leads to liberation and activation of Nrf2 (Kensler et al., 2007). In addition, subsets of both conjugated and unconjugated bile acids have been shown to activate multiple kinase signaling pathways such as protein kinase C, extracellular signal-regulated kinase 1/2, mitogen-activated protein kinase, p38 mitogen-activated protein kinase, c-Jun NH<sub>2</sub>-terminal kinase, and/or phosphatidylinositol 3-kinase/AKT (Debruyne et al., 2002; Dent et al., 2005). These signaling pathways have been shown as well to influence the stability of Nrf2-Keap1 complex and to post-transcriptionally regulate the Nrf2 target genes (Kensler et al., 2007).

In summary, we characterized a molecular cell defense event associated with bile acid-provoked oxidative stress. Exposure to cytotoxic bile acids in the liver and intestinal cells was shown here to cause Nrf2 activation, thereby up-regulating a battery of cytoprotective genes, particularly GCL subunits, to enhance cell survival at the emergence of oxidative stress.



**Fig. 5.** Nrf2 is a cellular protector against LCA-induced toxicity. a, cell viability of HepG2 after knockdown of Nrf2 via siRNA and treatment with toxic LCA (100 and 120  $\mu$ M). Viability was measured at indicated time points over 24 h. Fluorescent values of each LCA-treated siRNA group was first subtracted by the average values of individual siRNA group treated with vehicle DMSO, and they are expressed as percentage of change to baseline/pre-treatment values. \*,  $p < 0.01$ , significant difference between siCtr and siNrf2 groups with *t* test. Representative results; mean  $\pm$  S.E.M. of four determinations. b, analysis of LDH release ratio, a marker of necrotic event or cell injury, upon LCA challenge (90  $\mu$ M) in HepG2 after Nrf2 silencing. Values of LCA treatment were subtracted by the average values of vehicle treatment for each siRNA group. Representative results were presented; mean of duplicate determinations. c, analysis of caspases activity upon triggered by LCA toxicity (90  $\mu$ M) in HepG2 knockdown of Nrf2. Representative results of duplicate determinations were shown.



**Fig. 6.** Increased GSH production is a cellular protective mechanism against bile acid toxicity. **a**, cellular GSH levels of HepG2 after 24-h treatment with 75  $\mu$ M LCA, 100  $\mu$ M CDCA, or 200  $\mu$ M lipoic acid (LA; positive control). BSO (60  $\mu$ M), a GSH synthesis blocker, was used as negative control. \*,  $p < 0.05$ , significantly different from DMSO control by one-way ANOVA followed by post hoc test. Mean and S.E.M. ( $n = 4$ ). **b**, cell viability of HepG2 treated with toxic LCA (100  $\mu$ M) together with 0.4 mM NAC, or 60  $\mu$ M BSO, or with overnight pretreatment with 60  $\mu$ M BSO (pre-BSO). Fluorescent units were first subtracted by those of respective treatment control (i.e., DMSO, NAC, or BSO alone), and they are expressed as percentage of change to individual baseline/pretreatment values. \*,  $p < 0.05$ , significantly different from LCA treatment by one-way ANOVA followed by post hoc test. Mean  $\pm$  S.E.M. of three independent experiments with four determinations. **c**, ARE-reporter activity in HepG2 treated with LCA with or without 0.4 mM NAC or 60  $\mu$ M BSO. *tert*-BHP (100  $\mu$ M), a peroxide radical generator, was used to show that increased cellular oxidative stress led to increased ARE reporter activity. *y*-axis, -fold change in the ratio of luciferase activity (relative luciferase unit) from those transfected with basic pGL3 promoter construct and treated with vehicle DMSO. \*,  $p < 0.05$ , significant difference. Chemical treatments were for 8 h. Mean and S.E.M. ( $n = 3-4$ ).

## Acknowledgments

We thank Dr. Jawed Alam for providing expression vector plasmids, Dr. Patricia Harper for providing Hepa1c1c7 cells, and Christopher Tierney (Dharmacon RNA Technologies) for technical assistance in RNA interference studies.

## References

- Ananthanarayanan M, Balasubramanian N, Makishima M, Mangelsdorf DJ, and Suchy FJ (2001) Human bile salt export pump promoter is transactivated by the farnesoid X receptor/bile acid receptor. *J Biol Chem* **276**:28857–28865.
- Aw TY (2005) Intestinal glutathione: determinant of mucosal peroxide transport, metabolism, and oxidative susceptibility. *Toxicol Appl Pharmacol* **204**:320–328.
- Beresford GW and Boss JM (2001) CIITA coordinates multiple histone acetylation modifications at the HLA-DRA promoter. *Nat Immunol* **2**:652–657.
- Brown NS and Bicknell R (2001) Hypoxia and oxidative stress in breast cancer. Oxidative stress: its effects on the growth, metastatic potential and response to therapy of breast cancer. *Br Cancer Res* **3**:323–327.
- Calkins MJ, Jakel RJ, Johnson DA, Chan K, Kan YW, and Johnson JA (2005) Protection from mitochondrial complex II inhibition in vitro and in vivo by Nrf2-mediated transcription. *Proc Natl Acad Sci U S A* **102**:244–249.
- Chan K, Han XD, and Kan YW (2001) An important function of Nrf2 in combating oxidative stress: detoxification of acetaminophen. *Proc Natl Acad Sci U S A* **98**:4611–4616.
- Chawla A, Saez E, and Evans RM (2000) "Don't know much bile-ology". *Cell* **103**:1–4.
- Chen ZH, Saito Y, Yoshida Y, Sekine A, Noguchi N, and Niki E (2005) 4-Hydroxynon-ene induces adaptive response and enhances PC12 cell tolerance primarily through induction of thioredoxin reductase 1 via activation of Nrf2. *J Biol Chem* **280**:41921–41927.
- Debruyne PR, Bruyneel EA, Karaguni IM, Li X, Flatau G, Muller O, Zimmer A, Gespach C, and Mareel MM (2002) Bile acids stimulate invasion and haptotaxis in human colorectal cancer cells through activation of multiple oncogenic signaling pathways. *Oncogene* **21**:6740–6750.
- Dent P, Fang Y, Gupta S, Studer E, Mitchell C, Spiegel S, and Hylemon PB (2005) Conjugated bile acids promote ERK1/2 and AKT activation via a pertussis toxin-sensitive mechanism in murine and human hepatocytes. *Hepatology* **42**:1291–1299.
- Echchgadda I, Song CS, Roy AK, and Chatterjee B (2004) Dehydroepiandrosterone sulfotransferase is a target for transcriptional induction by the vitamin D receptor. *Mol Pharmacol* **65**:720–729.
- Eloranta JJ and Kullak-Ublick GA (2005) Coordinate transcriptional regulation of bile acid homeostasis and drug metabolism. *Arch Biochem Biophys* **433**:397–412.
- Erickson AM, Nevarea Z, Gipp JJ, and Mulcahy RT (2002) Identification of a variant antioxidant response element in the promoter of the human glutamate-cysteine ligase modifier subunit gene. Revision of the ARE consensus sequence. *J Biol Chem* **277**:30730–30737.
- Fernandez-Checa JC and Kaplowitz N (2005) Hepatic mitochondrial glutathione: transport and role in disease and toxicity. *Toxicol Appl Pharmacol* **204**:263–273.
- Glosli H, Tronstad KJ, Wergedal H, Muller F, Svardal A, Aukrust P, Berge RK, and Prydz H (2002) Human TNF- $\alpha$  in transgenic mice induces differential changes in redox status and glutathione-regulating enzymes. *FASEB J* **16**:1450–1452.
- Hofmann AF (2004) Detoxification of lithocholic acid, a toxic bile acid: relevance to drug hepatotoxicity. *Drug Metab Rev* **36**:703–722.
- Huang ZZ, Chen C, Zeng Z, Yang H, Oh J, Chen L, and Lu SC (2001) Mechanism and significance of increased glutathione level in human hepatocellular carcinoma and liver regeneration. *FASEB J* **15**:19–21.
- Kannan S and Jaiswal AK (2006) Low and High dose UVB regulation of transcription factor NF-E2-related factor 2. *Cancer Res* **66**:8421–8429.
- Kawalek JC, Hallmark RK, and Andrews AW (1983) Effect of lithocholic acid on the mutagenicity of some substituted aromatic amines. *J Natl Cancer Inst* **71**:293–298.
- Kensler TW, Wakabayashi N, and Biswal S (2007) Cell survival responses to environmental stress via the Keap1-Nrf2-ARE pathway. *Annu Rev Pharmacol Toxicol* **47**:89–116.
- Kozoni V, Tsioulas G, Shiff S, and Rigas B (2000) The effect of lithocholic acid on proliferation and apoptosis during the early stages of colon carcinogenesis: differential effect on apoptosis in the presence of a colon carcinogen. *Carcinogenesis* **21**:999–1005.
- Krähenbühl S, Fischer S, Talos C, and Reichen J (1994) Ursodeoxycholate protects oxidative mitochondrial metabolism from bile acid toxicity: dose-response study in isolated rat liver mitochondria. *Hepatology* **20**:1595–1601.
- Kullak-Ublick GA, Beuers U, and Paumgartner G (1996) Molecular and functional characterization of bile acid transport in human hepatoblastoma HepG2 cells. *Hepatology* **23**:1053–1060.
- Lee JM, Li J, Johnson DA, Stein TD, Kraft AD, Calkins MJ, Jakel RJ, and Johnson JA (2005) Nrf2, a multi-organ protector? *FASEB J* **19**:1061–1066.
- Livak KJ and Schmittgen TD (2001) Analysis of relative gene expression data using real-time quantitative PCR and the  $2^{-\Delta\Delta C(T)}$  method. *Methods* **25**:402–408.
- Makishima M, Okamoto AY, Repa JJ, Tu H, Learned RM, Luk A, Hull MV, Lustig KD, Mangelsdorf DJ, and Shan B (1999) Identification of a nuclear receptor for bile acids. *Science* **284**:1362–1365.
- Makishima M, Lu TT, Xie W, Whitfield GK, Domoto H, Evans RM, Haussler MR, and Mangelsdorf DJ (2002) Vitamin D receptor as an intestinal bile acid sensor. *Science* **296**:1313–1316.
- Morito N, Yoh K, Itoh K, Hirayama A, Koyama A, Yamamoto M, and Takahashi S (2003) Nrf2 regulates the sensitivity of death receptor signals by affecting intracellular glutathione levels. *Oncogene* **22**:9275–9281.
- Mulcahy RT, Wartman MA, Bailey HH, and Gipp JJ (1997) Constitutive and betanaphthoflavone-induced expression of the human gamma-glutamylcysteine syn-

- thetase heavy subunit gene is regulated by a distal antioxidant response element/TRE sequence. *J Biol Chem* **272**:7445–7454.
- Palmeira CM and Rolo AP (2004) Mitochondrially-mediated toxicity of bile acids. *Toxicology* **203**:1–15.
- Ramos-Gomez M, Kwak M-K, Dolan PM, Itoh K, Yamamoto M, Talalay P, and Kensler TW (2001) Sensitivity to carcinogenesis is increased and chemoprotective efficacy of enzyme inducers is lost in nrf2 transcription factor-deficient mice. *Proc Natl Acad Sci U S A* **98**:3410–3415.
- Rao CV, Hirose Y, Indranie C, and Reddy BS (2001) Modulation of experimental colon tumorigenesis by types and amounts of dietary fatty acids. *Cancer Res* **61**:1927–1933.
- Rundlöf AK and Arner ES (2004) Regulation of the mammalian selenoprotein thioredoxin reductase 1 in relation to cellular phenotype, growth, and signaling events. *Antioxid Redox Signal* **6**:41–52.
- Schoemaker MH, Gommans WM, Conde de al Rosa L, Homan M, Klok P, Teautwein C, van Goor H, Poelstra K, Haisma HJ, Jansen PL, et al. (2003) Resistance of rat hepatocytes against bile acid-induced apoptosis in cholestatic liver injury is due to nuclear factor-kappa B activation. *J Hepatol* **39**:153–161.

- Sokol RJ, Winkhofer-Roob BM, Devereaux MW, and McKim JM Jr (1995) Generation of hydroperoxides in isolated rat hepatocytes and hepatic mitochondria exposed to hydrophobic bile acids. *Gastroenterology* **109**:1249–1256.
- Staudinger JL, Goodwin B, Jones SA, Hawkins-Brown D, MacKenzie KI, LaTour A, Liu Y, Klaassen CD, Brown KK, Reinhard J, et al. (2001) The nuclear receptor PXR is a lithocholic acid sensor that protects against liver toxicity. *Proc Natl Acad Sci U S A* **98**:3369–3374.
- Tietze F (1969) Enzymic method for quantitative determination of nanogram amounts of total and oxidized glutathione: applications to mammalian blood and other tissues. *Anal Biochem* **27**:502–522.
- Vendemiale G, Gratiagiano I, Lupo L, Memeo V, and Altomare E (2002) Hepatic oxidative alterations in patients with extra-hepatic cholestasis. Effects of surgical drainage. *J Hepatol* **37**:601–605.

---

**Address correspondence to:** Dr. Shinya Ito, Division of Clinical Pharmacology and Toxicology, Hospital for Sick Children, 555 University Ave., Toronto, Ontario M5G 1X8, Canada. E-mail: shinya.ito@sickkids.ca

---