



3914 DEL AMO BLVD.
SUITE 901
TORRANCE CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Product name: Human Embryonic Stem Cell Total RNA -50µg

Catalog number: 36101RNA

Description: Human Embryonic Stem Cell Total RNA (50µg) extracted from Human Embryonic derived Stem Cells.

Storage Conditions: -80 °C

Tests performed:

1. Integrity of the RNA was examined by visual inspection for the presence of intact bands of 18s and 28s ribosomal RNA when electrophoreses on a denaturing agarose gel. Quality and purity of total RNA tested by spectrophotometer. A260/280 is between 1.8 and 2.0 (detected in 10 mM Tris-Cl, pH 7.5.). The ratio of 28S/18S is ≥ 1 .
2. The RNA was treated by DNase I, and is tested as DNA free RNA.

Expiration: 12 months when stored at -80°C

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits/Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* All products are established by CELPROGEN Inc., and is for Research Use Only. This product is not for resale or may not be transferred to Third party prior to written request and approval by CELPROGEN Inc.