



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Pancreatic Stem Cell Extracellular Differentiation Matrix (ECM)

Catalog number: D36004-05

Description: Human Pancreatic Stem Cell Extracellular Differentiation Matrix for growth, differentiation, and proliferation of Human Pancreatic Stem Cell in tissue culture with their appropriate media. The Human Pancreatic Stem Cell Differentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
D36004-05-T25	Human Pancreatic Stem Cell Differentiation Extracellular Matrix - T25 Flask (10/Pk)
D36004-05-T75	Human Pancreatic Stem Cell Differentiation Extracellular Matrix - T75 Flask (5/Pk)
D36004-05-T150	Human Pancreatic Stem Cell Differentiation Extracellular Matrix - T150 Flasks (5/pk)
D36004-05-T225	Human Pancreatic Stem Cell Differentiation Extracellular Matrix - T225 Flasks (5/pk)
D36004-05-4Well	Human Pancreatic Stem Cell Differentiation ECM - 4 Well Plates (5/Pk)
D36004-05-4MCS	Human Pancreatic Stem Cell Differentiation ECM - 4 Well Microscope Chamber Slides
D36004-05-MS4	Human Pancreatic Stem Cell Differentiation ECM 4 Well Microscope Slide
D36004-05-6Well	Human Pancreatic Stem Cell Differentiation ECM - 6 Well Plates (5/Pk)
D36004-05-8MCS	Human Pancreatic Stem Cell Differentiation ECM - 8 Well Microscope Chamber Slides
D36004-05-MS8	Human Pancreatic Stem Cell Differentiation ECM 8 well Microscope Slide
D36004-05-MS10	Human Pancreatic Stem Cell Differentiation ECM 10 Well Microscope Slides
D36004-05-12Well	Human Pancreatic Stem Cell Differentiation ECM - 12 Well Plates (5/Pk)
D36004-05-MS18	Human Pancreatic Stem Cell Differentiation ECM - 18 Well Microscope Slides
D36004-05-24Well	Human Pancreatic Stem Cell Differentiation ECM - 24 Well Plates (5/Pk)
D36004-05-48Well	Human Pancreatic Stem Cell Differentiation ECM - 48 Well Plates (5/Pk)
D36004-05-96Well	Human Pancreatic Stem Cell Differentiation ECM - 96 Well Plates (5/Pk)
D36004-05-96WW	Human Pancreatic Stem Cell Differentiation ECM - 96 White Well Plates (5/pk)
D36004-05-96BW	Human Pancreatic Stem Cell Differentiation ECM - 96 Black Well Plates (5/pk)
D36004-05-384WW	Human Pancreatic Stem Cell Differentiation ECM - 384 White Well Plates (5/pk)
D36004-05-384BW	Human Pancreatic Stem Cell Differentiation ECM - 384 Black Well Plates (5/pk)
D36004-05-CS12	Human Pancreatic Stem Cell Differentiation ECM - 12mm coverslips (10/pk)
D36004-05-CS15	Human Pancreatic Stem Cell Differentiation ECM - 15mm coverslips (10/pk)
D36004-05-CS18	Human Pancreatic Stem Cell Differentiation ECM - 18mm coverslips (10/pk)
D36004-05-CS22	Human Pancreatic Stem Cell Differentiation ECM - 22mm coverslips (10/pk)
D36004-05-CS25	Human Pancreatic Stem Cell Differentiation ECM - 25mm coverslips (10/pk)



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

D36004-05-PD6	Human Pancreatic Stem Cell Differentiation ECM - 6cm Petri Dish (5/pk)
D36004-05-PD10	Human Pancreatic Stem Cell Differentiation ECM - 10cm Petri Dish (5/pk)
D36004-05-6W-GB	Human Pancreatic Stem Cell Differentiation ECM 6 Well Glass Bottom Plates (5/pk)
D36004-05-96W-GB	Human Pancreatic Stem Cell Differentiation ECM 96 Well Glass Bottom Plates (5/pk)
D36004-05-3D-6Well	Human Pancreatic Stem 3D Cell Culture System Differentiation ECM 6 Well (1/pk)
D36004-05-3D-12Well	Human Pancreatic Stem 3D Cell Culture System Differentiation ECM 12 Well (1/pk)
D36004-05-3D-24Well	Human Pancreatic Stem 3D Cell Culture System Differentiation ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months
 Room temperature for short term storage 12 months

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Pancreatic Stem Cell Differentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.