



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

Human Neuronal Stem Cell (Hippocampus) Extracellular Differentiation Matrix (ECM)
---

**Catalog number:** D36119-47

**Description:** Human Neuronal Stem Cell (Hippocampus) Extracellular Differentiation Matrix for growth, differentiation, and proliferation of Human Neuronal Stem Cell (Hippocampus) in tissue culture with their appropriate media. The Human Neuronal Stem Cell (Hippocampus) Differentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
D36119-47-T25	Human Neuronal Stem Cell (Hippocampus) Differentiation Extracellular Matrix - T25 Flask (10/Pk)
D36119-47-T75	Human Neuronal Stem Cell (Hippocampus) Differentiation Extracellular Matrix - T75 Flask (5/Pk)
D36119-47-T150	Human Neuronal Stem Cell (Hippocampus) Differentiation Extracellular Matrix - T150 Flasks (5/pk)
D36119-47-T225	Human Neuronal Stem Cell (Hippocampus) Differentiation Extracellular Matrix - T225 Flasks (5/pk)
D36119-47-4Well	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 4 Well Plates (5/Pk)
D36119-47-4MCS	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 4 Well Microscope Chamber Slides
D36119-47-MS4	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM 4 Well Microscope Slide
D36119-47-6Well	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 6 Well Plates (5/Pk)
D36119-47-8MCS	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 8 Well Microscope Chamber Slides
D36119-47-MS8	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM 8 well Microscope Slide
D36119-47-MS10	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM 10 Well Microscope Slides
D36119-47-12Well	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 12 Well Plates (5/Pk)
D36119-47-MS18	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 18 Well Microscope Slides
D36119-47-24Well	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 24 Well Plates (5/Pk)
D36119-47-48Well	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 48 Well Plates (5/Pk)
D36119-47-96Well	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 96 Well Plates (5/Pk)



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

D36119-47-96WW	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 96 White Well Plates (5/pk)
D36119-47-96BW	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 96 Black Well Plates (5/pk)
D36119-47-384WW	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 384 White Well Plates (5/pk)
D36119-47-384BW	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 384 Black Well Plates (5/pk)
D36119-47-CS12	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 12mm coverslips (10/pk)
D36119-47-CS15	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 15mm coverslips (10/pk)
D36119-47-CS18	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 18mm coverslips (10/pk)
D36119-47-CS22	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 22mm coverslips (10/pk)
D36119-47-CS25	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 25mm coverslips (10/pk)
D36119-47-PD6	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 6cm Petri Dish (5/pk)
D36119-47-PD10	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM - 10cm Petri Dish (5/pk)
D36119-47-6W-GB	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM 6 Well Glass Bottom Plates (5/pk)
D36119-47-96W-GB	Human Neuronal Stem Cell (Hippocampus) Differentiation ECM 96 Well Glass Bottom Plates (5/pk)
D36119-47-3D-6Well	Human Neuronal SC 3D Cell Culture System Differentiation ECM 6 Well (1/pk)
D36119-47-3D-12Well	Human Neuronal SC 3D Cell Culture System Differentiation ECM 12 Well (1/pk)
D36119-47-3D-24Well	Human Neuronal SC 3D Cell Culture System Differentiation ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months  
 Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Neuronal Stem Cell (Hippocampus) Differentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use



1871 N GAFFEY STREET  
SUITE A & B  
SAN PEDRO CA 90731  
[www.celprogen.com](http://www.celprogen.com)

Phone: 310 547 3975  
Fax: 310 547 2975  
Email: [info@celprogen.com](mailto:info@celprogen.com)  
[stemcells@celprogen.com](mailto:stemcells@celprogen.com)

Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.