



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

### Rat Bone Marrow Derived Stem Cell Extracellular Differentiation Matrix (ECM)

**Catalog number:** D55093-20

**Description:** Rat Bone Marrow Derived Stem Cell Extracellular Differentiation Matrix for growth, differentiation, and proliferation of Rat Bone Marrow Derived Stem Cell in tissue culture with their appropriate media. The Rat Bone Marrow Derived Stem Cell Differentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
D55093-20-T25	Rat Bone Marrow Derived Stem Cell Differentiation Extracellular Matrix - T25 Flask (10/Pk)
D55093-20-T75	Rat Bone Marrow Derived Stem Cell Differentiation Extracellular Matrix - T75 Flask (5/Pk)
D55093-20-T150	Rat Bone Marrow Derived Stem Cell Differentiation Extracellular Matrix - T150 Flasks (5/pk)
D55093-20-T225	Rat Bone Marrow Derived Stem Cell Differentiation Extracellular Matrix - T225 Flasks (5/pk)
D55093-20-4Well	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 4 Well Plates (5/Pk)
D55093-20-4MCS	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 4 Well Microscope Chamber Slides
D55093-20-MS4	Rat Bone Marrow Derived Stem Cell Differentiation ECM 4 Well Microscope Slide
D55093-20-6Well	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 6 Well Plates (5/Pk)
D55093-20-8MCS	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 8 Well Microscope Chamber Slides
D55093-20-MS8	Rat Bone Marrow Derived Stem Cell Differentiation ECM 8 well Microscope Slide
D55093-20-MS10	Rat Bone Marrow Derived Stem Cell Differentiation ECM 10 Well Microscope Slides
D55093-20-12Well	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 12 Well Plates (5/Pk)
D55093-20-MS18	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 18 Well Microscope Slides
D55093-20-24Well	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 24 Well Plates (5/Pk)
D55093-20-48Well	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 48 Well Plates (5/Pk)
D55093-20-96Well	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 96 Well Plates (5/Pk)
D55093-20-96WW	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 96 White Well Plates (5/pk)
D55093-20-96BW	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 96 Black Well Plates (5/pk)



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

D55093-20-384WW	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 384 White Well Plates (5/pk)
D55093-20-384BW	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 384 Black Well Plates (5/pk)
D55093-20-CS12	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 12mm coverslips (10/pk)
D55093-20-CS15	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 15mm coverslips (10/pk)
D55093-20-CS18	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 18mm coverslips (10/pk)
D55093-20-CS22	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 22mm coverslips (10/pk)
D55093-20-CS25	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 25mm coverslips (10/pk)
D55093-20-PD6	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 6cm Petri Dish (5/pk)
D55093-20-PD10	Rat Bone Marrow Derived Stem Cell Differentiation ECM - 10cm Petri Dish (5/pk)
D55093-20-6W-GB	Rat Bone Marrow Derived Stem Cell Differentiation ECM 6 Well Glass Bottom Plates (5/pk)
D55093-20-96W-GB	Rat Bone Marrow Derived Stem Cell Differentiation ECM 96 Well Glass Bottom Plates (5/pk)
D55093-20-3D-6Well	Rat Bone Marrow Derived Stem 3D Cell Culture System Differentiation ECM 6 Well (1/pk)
D55093-20-3D-12Well	Rat Bone Marrow Derived Stem 3D Cell Culture System Differentiation ECM 12 Well (1/pk)
D55093-20-3D-24Well	Rat Bone Marrow Derived Stem 3D Cell Culture System Differentiation ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Rat Bone Marrow Derived Stem Cell Differentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.