



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Mouse Keratinocyte Stem Cell Extracellular Differentiation Matrix (ECM)

Catalog number: D66008-09

Description: Mouse Keratinocyte Stem Cell Extracellular Differentiation Matrix for growth, differentiation, and proliferation of Mouse Keratinocyte Stem Cell in tissue culture with their appropriate media. The Mouse Keratinocyte Stem Cell Differentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
D66008-09-T25	Mouse Keratinocyte Stem Cell Differentiation Extracellular Matrix - T25 Flask (10/Pk)
D66008-09-T75	Mouse Keratinocyte Stem Cell Differentiation Extracellular Matrix - T75 Flask (5/Pk)
D66008-09-T150	Mouse Keratinocyte Stem Cell Differentiation Extracellular Matrix - T150 Flasks (5/pk)
D66008-09-T225	Mouse Keratinocyte Stem Cell Differentiation Extracellular Matrix - T225 Flasks (5/pk)
D66008-09-4Well	Mouse Keratinocyte Stem Cell Differentiation ECM - 4 Well Plates (5/Pk)
D66008-09-4MCS	Mouse Keratinocyte Stem Cell Differentiation ECM - 4 Well Microscope Chamber Slides
D66008-09-MS4	Mouse Keratinocyte Stem Cell Differentiation ECM 4 Well Microscope Slide
D66008-09-6Well	Mouse Keratinocyte Stem Cell Differentiation ECM - 6 Well Plates (5/Pk)
D66008-09-8MCS	Mouse Keratinocyte Stem Cell Differentiation ECM - 8 Well Microscope Chamber Slides
D66008-09-MS8	Mouse Keratinocyte Stem Cell Differentiation ECM 8 well Microscope Slide
D66008-09-MS10	Mouse Keratinocyte Stem Cell Differentiation ECM 10 Well Microscope Slides
D66008-09-12Well	Mouse Keratinocyte Stem Cell Differentiation ECM - 12 Well Plates (5/Pk)
D66008-09-MS18	Mouse Keratinocyte Stem Cell Differentiation ECM - 18 Well Microscope Slides
D66008-09-24Well	Mouse Keratinocyte Stem Cell Differentiation ECM - 24 Well Plates (5/Pk)
D66008-09-48Well	Mouse Keratinocyte Stem Cell Differentiation ECM - 48 Well Plates (5/Pk)
D66008-09-96Well	Mouse Keratinocyte Stem Cell Differentiation ECM - 96 Well Plates (5/Pk)
D66008-09-96WW	Mouse Keratinocyte Stem Cell Differentiation ECM - 96 White Well Plates (5/pk)
D66008-09-96BW	Mouse Keratinocyte Stem Cell Differentiation ECM - 96 Black Well Plates (5/pk)
D66008-09-384WW	Mouse Keratinocyte Stem Cell Differentiation ECM - 384 White Well Plates (5/pk)
D66008-09-384BW	Mouse Keratinocyte Stem Cell Differentiation ECM - 384 Black Well Plates (5/pk)
D66008-09-CS12	Mouse Keratinocyte Stem Cell Differentiation ECM - 12mm coverslips (10/pk)
D66008-09-CS15	Mouse Keratinocyte Stem Cell Differentiation ECM - 15mm coverslips (10/pk)
D66008-09-CS18	Mouse Keratinocyte Stem Cell Differentiation ECM - 18mm coverslips (10/pk)
D66008-09-CS22	Mouse Keratinocyte Stem Cell Differentiation ECM - 22mm coverslips (10/pk)
D66008-09-CS25	Mouse Keratinocyte Stem Cell Differentiation ECM - 25mm coverslips (10/pk)



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

D66008-09-PD6	Mouse Keratinocyte Stem Cell Differentiation ECM - 6cm Petri Dish (5/pk)
D66008-09-PD10	Mouse Keratinocyte Stem Cell Differentiation ECM - 10cm Petri Dish (5/pk)
D66008-09-6W-GB	Mouse Keratinocyte Stem Cell Differentiation ECM 6 Well Glass Bottom Plates (5/pk)
D66008-09-96W-GB	Mouse Keratinocyte Stem Cell Differentiation ECM 96 Well Glass Bottom Plates (5/pk)
D66008-09-3D-6Well	Mouse Keratinocyte Stem Cell3D Cell Culture System Differentiation ECM 6 Well (1/pk)
D66008-09-3D-12Well	Mouse Keratinocyte Stem Cell3D Cell Culture System Differentiation ECM 12 Well (1/pk)
D66008-09-3D-24Well	Mouse Keratinocyte Stem Cell3D Cell Culture System Differentiation ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers:

CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans. This product, Mouse Keratinocyte Stem Cell Differentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.