



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

|   |
|---|
| Mouse Epithelial Stem Cell Extracellular Differentiation Matrix (ECM) |
|---|

**Catalog number:** D66104-31

**Description:** Mouse Epithelial Stem Cell Extracellular Differentiation Matrix for growth, differentiation, and proliferation of Mouse Epithelial Stem Cell in tissue culture with their appropriate media. The Mouse Epithelial Stem Cell Differentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

| Catalog          | Description  |
|------------------|--|
| D66104-31-T25    | Mouse Epithelial Stem Cell Differentiation Extracellular Matrix - T25 Flask (10/Pk)  |
| D66104-31-T75    | Mouse Epithelial Stem Cell Differentiation Extracellular Matrix - T75 Flask (5/Pk)   |
| D66104-31-T150   | Mouse Epithelial Stem Cell Differentiation Extracellular Matrix - T150 Flasks (5/pk) |
| D66104-31-T225   | Mouse Epithelial Stem Cell Differentiation Extracellular Matrix - T225 Flasks (5/pk) |
| D66104-31-4Well  | Mouse Epithelial Stem Cell Differentiation ECM - 4 Well Plates (5/Pk)                |
| D66104-31-4MCS   | Mouse Epithelial Stem Cell Differentiation ECM - 4 Well Microscope Chamber Slides    |
| D66104-31-MS4    | Mouse Epithelial Stem Cell Differentiation ECM 4 Well Microscope Slide               |
| D66104-31-6Well  | Mouse Epithelial Stem Cell Differentiation ECM - 6 Well Plates (5/Pk)                |
| D66104-31-8MCS   | Mouse Epithelial Stem Cell Differentiation ECM - 8 Well Microscope Chamber Slides    |
| D66104-31-MS8    | Mouse Epithelial Stem Cell Differentiation ECM 8 well Microscope Slide               |
| D66104-31-MS10   | Mouse Epithelial Stem Cell Differentiation ECM 10 Well Microscope Slides             |
| D66104-31-12Well | Mouse Epithelial Stem Cell Differentiation ECM - 12 Well Plates (5/Pk)               |
| D66104-31-MS18   | Mouse Epithelial Stem Cell Differentiation ECM - 18 Well Microscope Slides           |
| D66104-31-24Well | Mouse Epithelial Stem Cell Differentiation ECM - 24 Well Plates (5/Pk)               |
| D66104-31-48Well | Mouse Epithelial Stem Cell Differentiation ECM - 48 Well Plates (5/Pk)               |
| D66104-31-96Well | Mouse Epithelial Stem Cell Differentiation ECM - 96 Well Plates (5/Pk)               |
| D66104-31-96WW   | Mouse Epithelial Stem Cell Differentiation ECM - 96 White Well Plates (5/pk)         |
| D66104-31-96BW   | Mouse Epithelial Stem Cell Differentiation ECM - 96 Black Well Plates (5/pk)         |
| D66104-31-384WW  | Mouse Epithelial Stem Cell Differentiation ECM - 384 White Well Plates (5/pk)        |
| D66104-31-384BW  | Mouse Epithelial Stem Cell Differentiation ECM - 384 Black Well Plates (5/pk)        |
| D66104-31-CS12   | Mouse Epithelial Stem Cell Differentiation ECM - 12mm coverslips (10/pk)             |
| D66104-31-CS15   | Mouse Epithelial Stem Cell Differentiation ECM - 15mm coverslips (10/pk)             |
| D66104-31-CS18   | Mouse Epithelial Stem Cell Differentiation ECM - 18mm coverslips (10/pk)             |
| D66104-31-CS22   | Mouse Epithelial Stem Cell Differentiation ECM - 22mm coverslips (10/pk)             |
| D66104-31-CS25   | Mouse Epithelial Stem Cell Differentiation ECM - 25mm coverslips (10/pk)             |



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

|                     |   |
|---------------------|---|
| D66104-31-PD6       | Mouse Epithelial Stem Cell Differentiation ECM - 6cm Petri Dish (5/pk)            |
| D66104-31-PD10      | Mouse Epithelial Stem Cell Differentiation ECM - 10cm Petri Dish (5/pk)           |
| D66104-31-6W-GB     | Mouse Epithelial Stem Cell Differentiation ECM 6 Well Glass Bottom Plates (5/pk)  |
| D66104-31-96W-GB    | Mouse Epithelial Stem Cell Differentiation ECM 96 Well Glass Bottom Plates (5/pk) |
| D66104-31-3D-6Well  | Mouse Epithelial Stem 3D Cell Culture System Differentiation ECM 6 Well (1/pk)    |
| D66104-31-3D-12Well | Mouse Epithelial Stem 3D Cell Culture System Differentiation ECM 12 Well (1/pk)   |
| D66104-31-3D-24Well | Mouse Epithelial Stem 3D Cell Culture System Differentiation ECM 24 Well (1/pk)   |

**Storage Conditions:** 2- 8 °C long term storage 18 months  
 Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Mouse Epithelial Stem Cell Differentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.