



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

Human Hepatocyte iPC Differentiated Cell (Adenovirus) ECM
-----------------------------------------------------------

**Catalog number:** D933003-02

**Description:** Human Hepatocyte iPC Stem Cell Extracellular Differentiation Matrix for growth, differentiation, and proliferation of Human Hepatocyte iPC Stem cells in tissue culture with their appropriate media. The Human Hepatocyte iPC Stem Cell Differentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
D933003-02-T25	Human Hepatocyte iPC Differentiation Extracellular Matrix - T25 Flask (10/Pk)
D933003-02-T75	Human Hepatocyte iPC Differentiation Extracellular Matrix - T75 Flask (5/Pk)
D933003-02-T150	Human Hepatocyte iPC Differentiation Extracellular Matrix - T150 Flasks (5/pk)
D933003-02-T225	Human Hepatocyte iPC Differentiation Extracellular Matrix - T225 Flasks (5/pk)
D933003-02-4Well	Human Hepatocyte iPC Differentiation ECM - 4 Well Plates (5/Pk)
D933003-02-4MCS	Human Hepatocyte iPC Differentiation ECM - 4 Well Microscope Chamber Slides
D933003-02-MS4	Human Hepatocyte iPC Differentiation ECM 4 Well Microscope Slide
D933003-02-6Well	Human Hepatocyte iPC Differentiation ECM - 6 Well Plates (5/Pk)
D933003-02-8MCS	Human Hepatocyte iPC Differentiation ECM - 8 Well Microscope Chamber Slides
D933003-02-MS8	Human Hepatocyte iPC Differentiation ECM 8 well Microscope Slide
D933003-02-MS10	Human Hepatocyte iPC Differentiation ECM 10 Well Microscope Slides
D933003-02-12Well	Human Hepatocyte iPC Differentiation ECM - 12 Well Plates (5/Pk)
D933003-02-MS18	Human Hepatocyte iPC Differentiation ECM - 18 Well Microscope Slides
D933003-02-24Well	Human Hepatocyte iPC Differentiation ECM - 24 Well Plates (5/Pk)
D933003-02-48Well	Human Hepatocyte iPC Differentiation ECM - 48 Well Plates (5/Pk)
D933003-02-96Well	Human Hepatocyte iPC Differentiation ECM - 96 Well Plates (5/Pk)
D933003-02-96WW	Human Hepatocyte iPC Differentiation ECM - 96 White Well Plates (5/pk)
D933003-02-96BW	Human Hepatocyte iPC Differentiation ECM - 96 Black Well Plates (5/pk)
D933003-02-384WW	Human Hepatocyte iPC Differentiation ECM - 384 White Well Plates (5/pk)
D933003-02-384BW	Human Hepatocyte iPC Differentiation ECM - 384 Black Well Plates (5/pk)
D933003-02-CS12	Human Hepatocyte iPC Differentiation ECM - 12mm coverslips (10/pk)
D933003-02-CS15	Human Hepatocyte iPC Differentiation ECM - 15mm coverslips (10/pk)
D933003-02-CS18	Human Hepatocyte iPC Differentiation ECM - 18mm coverslips (10/pk)



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

D933003-02-CS22	Human Hepatocyte iPC Differentiation ECM - 22mm coverslips (10/pk)
D933003-02-CS25	Human Hepatocyte iPC Differentiation ECM - 25mm coverslips (10/pk)
D933003-02-PD6	Human Hepatocyte iPC Differentiation ECM - 6cm Petri Dish (5/pk)
D933003-02-PD10	Human Hepatocyte iPC Differentiation ECM - 10cm Petri Dish (5/pk)
D933003-02-6W-GB	Human Hepatocyte iPC Differentiation ECM 6 Well Glass Bottom Plates (5/pk)
D933003-02-96W-GB	Human Hepatocyte iPC Differentiation ECM 96 Well Glass Bottom Plates (5/pk)
D933003-02-3D-6Well	Human Hepatocyte iPC3D Cell Culture System Differentiation ECM 6 Well (1/pk)
D933003-02-3D-12Well	Human Hepatocyte iPC3D Cell Culture System Differentiation ECM 12 Well (1/pk)
D933003-02-3D-24Well	Human Hepatocyte iPC3D Cell Culture System Differentiation ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months  
 Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen’s Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Hepatocyte iPC Stem Cell Differentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.