



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

| |
|--|
| Human Fibroblast iPC Derived Stem Cell Differentiated Cell (Non-Viral) ECM |
|--|

Catalog number: D936041-01

Description: Human Fibroblast iPC Derived Stem Cell Extracellular Differentiation Matrix for growth, differentiation, and proliferation of Human Fibroblast iPC Derived Stem Cell in tissue culture with their appropriate media. The Human Fibroblast iPC Derived Stem Cell Differentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

| Catalog | Description |
|-------------------|--|
| D936041-01-T25 | Human Fibroblast iPC Differentiation Extracellular Matrix - T25 Flask (10/Pk) |
| D936041-01-T75 | Human Fibroblast iPC Differentiation Extracellular Matrix - T75 Flask (5/Pk) |
| D936041-01-T150 | Human Fibroblast iPC Differentiation Extracellular Matrix - T150 Flasks (5/pk) |
| D936041-01-T225 | Human Fibroblast iPC Differentiation Extracellular Matrix - T225 Flasks (5/pk) |
| D936041-01-4Well | Human Fibroblast iPC Differentiation ECM - 4 Well Plates (5/Pk) |
| D936041-01-4MCS | Human Fibroblast iPC Differentiation ECM - 4 Well Microscope Chamber Slides |
| D936041-01-MS4 | Human Fibroblast iPC Differentiation ECM 4 Well Microscope Slide |
| D936041-01-6Well | Human Fibroblast iPC Differentiation ECM - 6 Well Plates (5/Pk) |
| D936041-01-8MCS | Human Fibroblast iPC Differentiation ECM - 8 Well Microscope Chamber Slides |
| D936041-01-MS8 | Human Fibroblast iPC Differentiation ECM 8 well Microscope Slide |
| D936041-01-MS10 | Human Fibroblast iPC Differentiation ECM 10 Well Microscope Slides |
| D936041-01-12Well | Human Fibroblast iPC Differentiation ECM - 12 Well Plates (5/Pk) |
| D936041-01-MS18 | Human Fibroblast iPC Differentiation ECM - 18 Well Microscope Slides |
| D936041-01-24Well | Human Fibroblast iPC Differentiation ECM - 24 Well Plates (5/Pk) |
| D936041-01-48Well | Human Fibroblast iPC Differentiation ECM - 48 Well Plates (5/Pk) |
| D936041-01-96Well | Human Fibroblast iPC Differentiation ECM - 96 Well Plates (5/Pk) |
| D936041-01-96WW | Human Fibroblast iPC Differentiation ECM - 96 White Well Plates (5/pk) |
| D936041-01-96BW | Human Fibroblast iPC Differentiation ECM - 96 Black Well Plates (5/pk) |
| D936041-01-384WW | Human Fibroblast iPC Differentiation ECM - 384 White Well Plates (5/pk) |
| D936041-01-384BW | Human Fibroblast iPC Differentiation ECM - 384 Black Well Plates (5/pk) |
| D936041-01-CS12 | Human Fibroblast iPC Differentiation ECM - 12mm coverslips (10/pk) |
| D936041-01-CS15 | Human Fibroblast iPC Differentiation ECM - 15mm coverslips (10/pk) |
| D936041-01-CS18 | Human Fibroblast iPC Differentiation ECM - 18mm coverslips (10/pk) |



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

| | |
|----------------------|--|
| D936041-01-CS22 | Human Fibroblast iPC Differentiation ECM - 22mm coverslips (10/pk) |
| D936041-01-CS25 | Human Fibroblast iPC Differentiation ECM - 25mm coverslips (10/pk) |
| D936041-01-PD6 | Human Fibroblast iPC Differentiation ECM - 6cm Petri Dish (5/pk) |
| D936041-01-PD10 | Human Fibroblast iPC Differentiation ECM - 10cm Petri Dish (5/pk) |
| D936041-01-6W-GB | Human Fibroblast iPC Differentiation ECM 6 Well Glass Bottom Plates (5/pk) |
| D936041-01-96W-GB | Human Fibroblast iPC Differentiation ECM 96 Well Glass Bottom Plates (5/pk) |
| D936041-01-3D-6Well | Human Fibroblast iPC 3D Cell Culture System Differentiation ECM 6 Well (1/pk) |
| D936041-01-3D-12Well | Human Fibroblast iPC 3D Cell Culture System Differentiation ECM 12 Well (1/pk) |
| D936041-01-3D-24Well | Human Fibroblast iPC 3D Cell Culture System Differentiation ECM 24 Well (1/pk) |

Storage Conditions: 2- 8 °C long term storage 18 months
 Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen’s Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Fibroblast iPC Derived Stem Cell Differentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.