



1871 N GAFFEY STREET
SUITE A & B
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Human Cardiomyocyte iPC Derived Stem Cell Differentiated Cell (Lentivirus) ECM

Catalog number: D936044-03

Description: Human Cardiomyocyte iPC Derived Stem Cell Extracellular Differentiation Matrix for growth, differentiation, and proliferation of Human Cardiomyocyte iPC Derived Stem Cells in tissue culture with their appropriate media. The Human Cardiomyocyte iPC Derived Stem Cell Differentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
D936044-03-T25	Human Cardiomyocyte iPC Differentiation Extracellular Matrix - T25 Flask (10/Pk)
D936044-03-T75	Human Cardiomyocyte iPC Differentiation Extracellular Matrix - T75 Flask (5/Pk)
D936044-03-T150	Human Cardiomyocyte iPC Differentiation Extracellular Matrix - T150 Flasks (5/pk)
D936044-03-T225	Human Cardiomyocyte iPC Differentiation Extracellular Matrix - T225 Flasks (5/pk)
D936044-03-4Well	Human Cardiomyocyte iPC Differentiation ECM - 4 Well Plates (5/Pk)
D936044-03-4MCS	Human Cardiomyocyte iPC Differentiation ECM - 4 Well Microscope Chamber Slides
D936044-03-MS4	Human Cardiomyocyte iPC Differentiation ECM 4 Well Microscope Slide
D936044-03-6Well	Human Cardiomyocyte iPC Differentiation ECM - 6 Well Plates (5/Pk)
D936044-03-8MCS	Human Cardiomyocyte iPC Differentiation ECM - 8 Well Microscope Chamber Slides
D936044-03-MS8	Human Cardiomyocyte iPC Differentiation ECM 8 well Microscope Slide
D936044-03-MS10	Human Cardiomyocyte iPC Differentiation ECM 10 Well Microscope Slides
D936044-03-12Well	Human Cardiomyocyte iPC Differentiation ECM - 12 Well Plates (5/Pk)
D936044-03-MS18	Human Cardiomyocyte iPC Differentiation ECM - 18 Well Microscope Slides
D936044-03-24Well	Human Cardiomyocyte iPC Differentiation ECM - 24 Well Plates (5/Pk)
D936044-03-48Well	Human Cardiomyocyte iPC Differentiation ECM - 48 Well Plates (5/Pk)
D936044-03-96Well	Human Cardiomyocyte iPC Differentiation ECM - 96 Well Plates (5/Pk)
D936044-03-96WW	Human Cardiomyocyte iPC Differentiation ECM - 96 White Well Plates (5/pk)
D936044-03-96BW	Human Cardiomyocyte iPC Differentiation ECM - 96 Black Well Plates (5/pk)
D936044-03-384WW	Human Cardiomyocyte iPC Differentiation ECM - 384 White Well Plates (5/pk)
D936044-03-384BW	Human Cardiomyocyte iPC Differentiation ECM - 384 Black Well Plates (5/pk)
D936044-03-CS12	Human Cardiomyocyte iPC Differentiation ECM - 12mm coverslips (10/pk)
D936044-03-CS15	Human Cardiomyocyte iPC Differentiation ECM - 15mm coverslips (10/pk)
D936044-03-CS18	Human Cardiomyocyte iPC Differentiation ECM - 18mm coverslips (10/pk)



1871 N GAFFEY STREET
SUITE A & B
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

D936044-03-CS22	Human Cardiomyocyte iPC Differentiation ECM - 22mm coverslips (10/pk)
D936044-03-CS25	Human Cardiomyocyte iPC Differentiation ECM - 25mm coverslips (10/pk)
D936044-03-PD6	Human Cardiomyocyte iPC Differentiation ECM - 6cm Petri Dish (5/pk)
D936044-03-PD10	Human Cardiomyocyte iPC Differentiation ECM - 10cm Petri Dish (5/pk)
D936044-03-6W-GB	Human Cardiomyocyte iPC Differentiation ECM 6 Well Glass Bottom Plates (5/pk)
D936044-03-96W-GB	Human Cardiomyocyte iPC Differentiation ECM 96 Well Glass Bottom Plates (5/pk)
D936044-03-3D-6Well	Human Cardiomyocyte iPC 3D Cell Culture System Differentiation ECM 6 Well (1/pk)
D936044-03-3D-12Well	Human Cardiomyocyte iPC 3D Cell Culture System Differentiation ECM 12 Well (1/pk)
D936044-03-3D-24Well	Human Cardiomyocyte iPC 3D Cell Culture System Differentiation ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices &
Disclaimers:**

CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans. This product, Human Cardiomyocyte iPC Derived Stem Cell Differentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.