



1871 N GAFFEY STREET
 SUITE A & EMBRYONIC HYPOTHALAMIC NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Email: info@celprogen.com
 stemcells@celprogen.com

Fax: 310 547 2975

DATA SHEET

Rat Embryonic Hypothalamic Neuron Cell Culture Extracellular Matrix (ECM)

Catalog number: E12040-01

Description: Rat Embryonic Hypothalamic Neuron Cell Culture Extracellular Matrix for growth, expansion, and proliferation of Rat Embryonic Hypothalamic Neuron Cell Culture in tissue culture with their appropriate media. The Rat Embryonic Hypothalamic Neuron Cell Culture ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
E12040-01-T25	Rat Embryonic Hypothalamic Neuron Cell Culture Extracellular Matrix - T25 Flask (10/Pk)
E12040-01-T75	Rat Embryonic Hypothalamic Neuron Cell Culture Extracellular Matrix - T75 Flask (5/Pk)
E12040-01-T150	Rat Embryonic Hypothalamic Neuron Cell Culture Extracellular Matrix - T150 Flasks (5/pk)
E12040-01-T225	Rat Embryonic Hypothalamic Neuron Cell Culture Extracellular Matrix - T225 Flasks (5/pk)
E12040-01-4Well	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 4 Well Plates (5/Pk)
E12040-01-4MCS	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 4 Well Microscope Chamber Slides
E12040-01-MS4	Rat Embryonic Hypothalamic Neuron Cell Culture ECM 4 Well Microscope Slide
E12040-01-6Well	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 6 Well Plates (5/Pk)
E12040-01-8MCS	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 8 Well Microscope Chamber Slides
E12040-01-MS8	Rat Embryonic Hypothalamic Neuron Cell Culture ECM 8 well Microscope Slide
E12040-01-MS10	Rat Embryonic Hypothalamic Neuron Cell Culture ECM 10 Well Microscope Slides
E12040-01-12Well	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 12 Well Plates (5/Pk)
E12040-01-MS18	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 18 Well Microscope Slides
E12040-01-24Well	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 24 Well Plates (5/Pk)
E12040-01-48Well	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 48 Well Plates (5/Pk)
E12040-01-96Well	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 96 Well Plates (5/Pk)
E12040-01-96WW	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 96 White Well Plates (5/pk)
E12040-01-96BW	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 96 Black Well Plates (5/pk)
E12040-01-384WW	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 384 White Well Plates (5/pk)



1871 N GAFFEY STREET
 SUITE A & EMBRYONIC HYPOTHALAMIC NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975

Fax: 310 547 2975

Email: info@celprogen.com
 stemcells@celprogen.com

E12040-01-384BW	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 384 Black Well Plates (5/pk)
E12040-01-CS12	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 12mm coverslips (10/pk)
E12040-01-CS15	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 15mm coverslips (10/pk)
E12040-01-CS18	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 18mm coverslips (10/pk)
E12040-01-CS22	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 22mm coverslips (10/pk)
E12040-01-CS25	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 25mm coverslips (10/pk)
E12040-01-PD6	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 6cm Petri Dish (5/pk)
E12040-01-PD10	Rat Embryonic Hypothalamic Neuron Cell Culture ECM - 10cm Petri Dish (5/pk)
E12040-01-6W-GB	Rat Embryonic Hypothalamic Neuron Cell Culture ECM 6 Well Glass Bottom Plates (5/pk)
E12040-01-96W-GB	Rat Embryonic Hypothalamic Neuron Cell Culture ECM 96 Well Glass Bottom Plates (5/pk)
E12040-01-3D-6Well	Rat Embryonic Hypothalamic Neuron Cell Culture 3D Culture System ECM 6 Well (1/pk)
E12040-01-3D-12Well	Rat Embryonic Hypothalamic Neuron Cell Culture 3D Culture System ECM 12 Well (1/pk)
E12040-01-3D-24Well	Rat Embryonic Hypothalamic Neuron Cell Culture 3D Culture System ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Rat Embryonic Hypothalamic Neuron Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.