



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Rat Aortic Artery Smooth Muscle Cell Culture Extracellular Matrix (ECM)

Catalog number: E12112-04

Description: Rat Aortic Artery Smooth Muscle Cell Culture Extracellular Matrix for growth, expansion, and proliferation of Rat Aortic Artery Smooth Muscle Cell Culture in tissue culture with their appropriate media. The Rat Aortic Artery Smooth Muscle Cell Culture ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
E12112-04-T25	Rat Aortic Artery Smooth Muscle Cell Culture Extracellular Matrix - T25 Flask (10/Pk)
E12112-04-T75	Rat Aortic Artery Smooth Muscle Cell Culture Extracellular Matrix - T75 Flask (5/Pk)
E12112-04-T150	Rat Aortic Artery Smooth Muscle Cell Culture Extracellular Matrix - T150 Flasks (5/pk)
E12112-04-T225	Rat Aortic Artery Smooth Muscle Cell Culture Extracellular Matrix - T225 Flasks (5/pk)
E12112-04-4Well	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 4 Well Plates (5/Pk)
E12112-04-4MCS	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 4 Well Microscope Chamber Slides
E12112-04-MS4	Rat Aortic Artery Smooth Muscle Cell Culture ECM 4 Well Microscope Slide
E12112-04-6Well	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 6 Well Plates (5/Pk)
E12112-04-8MCS	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 8 Well Microscope Chamber Slides
E12112-04-MS8	Rat Aortic Artery Smooth Muscle Cell Culture ECM 8 well Microscope Slide
E12112-04-MS10	Rat Aortic Artery Smooth Muscle Cell Culture ECM 10 Well Microscope Slides
E12112-04-12Well	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 12 Well Plates (5/Pk)
E12112-04-MS18	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 18 Well Microscope Slides
E12112-04-24Well	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 24 Well Plates (5/Pk)
E12112-04-48Well	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 48 Well Plates (5/Pk)
E12112-04-96Well	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 96 Well Plates (5/Pk)
E12112-04-96WW	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 96 White Well Plates (5/pk)
E12112-04-96BW	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 96 Black Well Plates (5/pk)
E12112-04-384WW	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 384 White Well Plates (5/pk)
E12112-04-384BW	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 384 Black Well Plates (5/pk)
E12112-04-CS12	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 12mm coverslips (10/pk)
E12112-04-CS15	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 15mm coverslips (10/pk)
E12112-04-CS18	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 18mm coverslips (10/pk)



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

E12112-04-CS22	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 22mm coverslips (10/pk)
E12112-04-CS25	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 25mm coverslips (10/pk)
E12112-04-PD6	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 6cm Petri Dish (5/pk)
E12112-04-PD10	Rat Aortic Artery Smooth Muscle Cell Culture ECM - 10cm Petri Dish (5/pk)
E12112-04-6W-GB	Rat Aortic Artery Smooth Muscle Cell Culture ECM 6 Well Glass Bottom Plates (5/pk)
E12112-04-96W-GB	Rat Aortic Artery Smooth Muscle Cell Culture ECM 96 Well Glass Bottom Plates (5/pk)
E12112-04-3D-6Well	Rat Aortic Artery Smooth Muscle Cell Culture 3D Culture System ECM 6 Well (1/pk)
E12112-04-3D-12Well	Rat Aortic Artery Smooth Muscle Cell Culture 3D Culture System ECM 12 Well (1/pk)
E12112-04-3D-24Well	Rat Aortic Artery Smooth Muscle Cell Culture 3D Culture System ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Rat Aortic Artery Smooth Muscle Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.