



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

### Human Pancreatic Islets of Langerhans Extracellular Matrix (ECM)

**Catalog number:** E35002-04

**Description:** Human Pancreatic Islets of Langerhans Extracellular Matrix for growth, expansion, and proliferation of Human Pancreatic Islets of Langerhans in tissue culture with their appropriate media. The Human Pancreatic Islets of Langerhans ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
E35002-04 -T25	Human Pancreatic Islets of Langerhans Extracellular Matrix - T25 Flask (10/Pk)
E35002-04 -T75	Human Pancreatic Islets of Langerhans Extracellular Matrix - T75 Flask (5/Pk)
E35002-04 -T150	Human Pancreatic Islets of Langerhans Extracellular Matrix - T150 Flasks (5/pk)
E35002-04 -T225	Human Pancreatic Islets of Langerhans Extracellular Matrix - T225 Flasks (5/pk)
E35002-04 -4Well	Human Pancreatic Islets of Langerhans ECM - 4 Well Plates (5/Pk)
E35002-04 -4MCS	Human Pancreatic Islets of Langerhans ECM - 4 Well Microscope Chamber Slides
E35002-04 -MS4	Human Pancreatic Islets of Langerhans ECM 4 Well Microscope Slide
E35002-04 -6Well	Human Pancreatic Islets of Langerhans ECM - 6 Well Plates (5/Pk)
E35002-04 -8MCS	Human Pancreatic Islets of Langerhans ECM - 8 Well Microscope Chamber Slides
E35002-04 -MS8	Human Pancreatic Islets of Langerhans ECM 8 well Microscope Slide
E35002-04 -MS10	Human Pancreatic Islets of Langerhans ECM 10 Well Microscope Slides
E35002-04 -12Well	Human Pancreatic Islets of Langerhans ECM - 12 Well Plates (5/Pk)
E35002-04 -MS18	Human Pancreatic Islets of Langerhans ECM - 18 Well Microscope Slides
E35002-04 -24Well	Human Pancreatic Islets of Langerhans ECM - 24 Well Plates (5/Pk)
E35002-04 -48Well	Human Pancreatic Islets of Langerhans ECM - 48 Well Plates (5/Pk)
E35002-04 -96Well	Human Pancreatic Islets of Langerhans ECM - 96 Well Plates (5/Pk)
E35002-04 -96WW	Human Pancreatic Islets of Langerhans ECM - 96 White Well Plates (5/pk)
E35002-04 -96BW	Human Pancreatic Islets of Langerhans ECM - 96 Black Well Plates (5/pk)
E35002-04 -384WW	Human Pancreatic Islets of Langerhans ECM - 384 White Well Plates (5/pk)
E35002-04 -384BW	Human Pancreatic Islets of Langerhans ECM - 384 Black Well Plates (5/pk)
E35002-04 -CS12	Human Pancreatic Islets of Langerhans ECM - 12mm coverslips (10/pk)
E35002-04 -CS15	Human Pancreatic Islets of Langerhans ECM - 15mm coverslips (10/pk)
E35002-04 -CS18	Human Pancreatic Islets of Langerhans ECM - 18mm coverslips (10/pk)
E35002-04 -CS22	Human Pancreatic Islets of Langerhans ECM - 22mm coverslips (10/pk)
E35002-04 -CS25	Human Pancreatic Islets of Langerhans ECM - 25mm coverslips (10/pk)



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

E35002-04 -PD6	Human Pancreatic Islets of Langerhans ECM - 6cm Petri Dish (5/pk)
E35002-04 -PD10	Human Pancreatic Islets of Langerhans ECM - 10cm Petri Dish (5/pk)
E35002-04 -6W-GB	Human Pancreatic Islets of Langerhans ECM 6 Well Glass Bottom Plates (5/pk)
E35002-04 -96W-GB	Human Pancreatic Islets of Langerhans ECM 96 Well Glass Bottom Plates (5/pk)
E35002-04 -3D-6Well	Human Pancreatic Islets of Langerhans 3D Culture System ECM 6 Well (1/pk)
E35002-04 -3D-12Well	Human Pancreatic Islets of Langerhans 3D Culture System ECM 12 Well (1/pk)
E35002-04 -3D-24Well	Human Pancreatic Islets of Langerhans 3D Culture System ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices &**

**Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Pancreatic Islets of Langerhans Primary Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.