



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Pancreatic Normal Primary Cell Culture Extracellular Matrix (ECM)
--

Catalog number: E36004-06

Description: Human Pancreatic Normal Primary Cell Culture Extracellular Matrix for growth, expansion, and proliferation of Human Pancreatic Normal Primary Cell in tissue culture with their appropriate media. The Human Pancreatic Normal Primary Cell ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop down menu. Each format is tissue culture tested with their appropriate media.

Catalog	Description
E36004-06-T25	Human Pancreatic Normal Primary Cell Extracellular Matrix - T25 Flask (10/Pk)
E36004-06-T75	Human Pancreatic Normal Primary Cell Extracellular Matrix - T75 Flask (5/Pk)
E36004-06-T150	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix - T150 Flasks (5/pk)
E36004-06-T225	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix - T225 Flasks (5/pk)
E36004-06-4Well	Human Pancreatic Normal Primary Cell Culture ECM - 4 Well Plates (5/Pk)
E36004-06-4MCS	Human Pancreatic Normal Primary Cell Extracellular Matrix - 4 Well Microscope Chamber Slides
E36004-06-MS4	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix 4 Well Microscope Slide
E36004-06-6Well	Human Pancreatic Normal Primary Cell Extracellular Matrix - 6 Well Plates (5/Pk)
E36004-06-8MCS	Human Pancreatic Normal Primary Cell Extracellular Matrix - 8 Well Microscope Chamber Slides
E36004-06-MS8	Human Pancreatic Normal Primary Cell Extracellular Matrix 8 well Microscope Slide
E36004-06-MS10	Human Pancreatic Normal Primary Cell Culture - 10 Well Microscope Slides
E36004-06-12Well	Human Pancreatic Normal Primary Cell Extracellular Matrix - 12 Well Plates (5/Pk)
E36004-06-MS18	Human Pancreatic Normal Primary Cell Culture - 18 Well Microscope Slides
E36004-06-24Well	Human Pancreatic Normal Primary Cell Extracellular Matrix - 24 Well Plates (5/Pk)
E36004-06-48Well	Human Pancreatic Normal Primary Cell Extracellular Matrix - 48 Well Plates (5/Pk)
E36004-06-96Well	Human Pancreatic Normal Primary Cell Extracellular Matrix - 96 Well Plates (5/Pk)
E36004-06-96WW	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix - 96 White Well Plates (5/pk)
E36004-06-96BW	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix - 96 Black Well Plates (5/pk)
E36004-06-384WW	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix - 384 White Well Plates (5/pk)



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

E36004-06-384BW	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix - 384 Black Well Plates (5/pk)
E36004-06-CS12	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix- 12mm coverslips (10/pk)
E36004-06-CS15	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix- 15mm coverslips (10/pk)
E36004-06-CS18	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix- 18mm coverslips (10/pk)
E36004-06-CS22	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix- 22mm coverslips (10/pk)
E36004-06-CS25	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix- 25mm coverslips (10/pk)
E36004-06-PD6	Human Pancreatic Normal Primary Cell Culture Extracellular Expansion Matrix - 6cm Petri Dish (5/pk)
E36004-06-PD10	Human Pancreatic Normal Primary Cell Culture Extracellular Expansion Matrix - 10cm Petri Dish (5/pk)
E36004-06-6W-GB	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix 6 Well Glass Bottom Plates (5/pk)
E36004-06-96W-GB	Human Pancreatic Normal Primary Cell Culture Extracellular Matrix 96 Well Glass Bottom Plates (5/pk)
E36004-06-3D-6Well	Human Pancreatic Normal Primary 3D Cell Culture System 6 Well Format (1/pk)
E36004-06-3D-12Well	Human Pancreatic Normal Primary 3D Cell Culture System 12 Well Format (1/pk)
E36004-06-3D-24Well	Human Pancreatic Normal Primary 3D Cell Culture System 24 Well Format (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months
 Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Pancreatic Normal Primary Cell Expansion Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.