



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Corneal Epithelial Extracellular Matrix (ECM)

Catalog number: E36045-16

Description: Human Corneal Epithelial Extracellular Matrix for growth, expansion, and proliferation of Human Corneal Epithelial in tissue culture with their appropriate media. The Human Corneal Epithelial ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
E36045-16-T25	Human Corneal Epithelial Extracellular Matrix - T25 Flask (10/Pk)
E36045-16-T75	Human Corneal Epithelial Extracellular Matrix - T75 Flask (5/Pk)
E36045-16-T150	Human Corneal Epithelial Extracellular Matrix - T150 Flasks (5/pk)
E36045-16-T225	Human Corneal Epithelial Extracellular Matrix - T225 Flasks (5/pk)
E36045-16-4Well	Human Corneal Epithelial ECM - 4 Well Plates (5/Pk)
E36045-16-4MCS	Human Corneal Epithelial ECM - 4 Well Microscope Chamber Slides
E36045-16-MS4	Human Corneal Epithelial ECM 4 Well Microscope Slide
E36045-16-6Well	Human Corneal Epithelial ECM - 6 Well Plates (5/Pk)
E36045-16-8MCS	Human Corneal Epithelial ECM - 8 Well Microscope Chamber Slides
E36045-16-MS8	Human Corneal Epithelial ECM 8 well Microscope Slide
E36045-16-MS10	Human Corneal Epithelial ECM 10 Well Microscope Slides
E36045-16-12Well	Human Corneal Epithelial ECM - 12 Well Plates (5/Pk)
E36045-16-MS18	Human Corneal Epithelial ECM - 18 Well Microscope Slides
E36045-16-24Well	Human Corneal Epithelial ECM - 24 Well Plates (5/Pk)
E36045-16-48Well	Human Corneal Epithelial ECM - 48 Well Plates (5/Pk)
E36045-16-96Well	Human Corneal Epithelial ECM - 96 Well Plates (5/Pk)
E36045-16-96WW	Human Corneal Epithelial ECM - 96 White Well Plates (5/pk)
E36045-16-96BW	Human Corneal Epithelial ECM - 96 Black Well Plates (5/pk)
E36045-16-384WW	Human Corneal Epithelial ECM - 384 White Well Plates (5/pk)
E36045-16-384BW	Human Corneal Epithelial ECM - 384 Black Well Plates (5/pk)
E36045-16-CS12	Human Corneal Epithelial ECM - 12mm coverslips (10/pk)
E36045-16-CS15	Human Corneal Epithelial ECM - 15mm coverslips (10/pk)
E36045-16-CS18	Human Corneal Epithelial ECM - 18mm coverslips (10/pk)
E36045-16-CS22	Human Corneal Epithelial ECM - 22mm coverslips (10/pk)
E36045-16-CS25	Human Corneal Epithelial ECM - 25mm coverslips (10/pk)



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

E36045-16-PD6	Human Corneal Epithelial ECM - 6cm Petri Dish (5/pk)
E36045-16-PD10	Human Corneal Epithelial ECM - 10cm Petri Dish (5/pk)
E36045-16-6W-GB	Human Corneal Epithelial ECM 6 Well Glass Bottom Plates (5/pk)
E36045-16-96W-GB	Human Corneal Epithelial ECM 96 Well Glass Bottom Plates (5/pk)
E36045-16-3D-6Well	Human Corneal Epithelial 3D Culture System ECM 6 Well (1/pk)
E36045-16-3D-12Well	Human Corneal Epithelial 3D Culture System ECM 12 Well (1/pk)
E36045-16-3D-24Well	Human Corneal Epithelial 3D Culture System ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices &

Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Corneal Epithelial Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.