



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

### Human Mesenchymal (Adipose Derived) Stem Cell Extracellular Matrix (ECM)

**Catalog number:** E36096-23

**Description:** Human Mesenchymal (Adipose Derived) Stem Cell Extracellular Matrix for growth, expansion, and proliferation of Human Mesenchymal (Adipose Derived) Stem Cell in tissue culture with their appropriate media. The Human Mesenchymal (Adipose Derived) Stem Cell ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
E36096-23-T25	Human Mesenchymal (Adipose Derived) Stem Cell Extracellular Matrix - T25 Flask (10/Pk)
E36096-23-T75	Human Mesenchymal (Adipose Derived) Stem Cell Extracellular Matrix - T75 Flask (5/Pk)
E36096-23-T150	Human Mesenchymal (Adipose Derived) Stem Cell Extracellular Matrix - T150 Flasks (5/pk)
E36096-23-T225	Human Mesenchymal (Adipose Derived) Stem Cell Extracellular Matrix - T225 Flasks (5/pk)
E36096-23-4Well	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 4 Well Plates (5/Pk)
E36096-23-4MCS	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 4 Well Microscope Chamber Slides
E36096-23-MS4	Human Mesenchymal (Adipose Derived) Stem Cell ECM 4 Well Microscope Slide
E36096-23-6Well	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 6 Well Plates (5/Pk)
E36096-23-8MCS	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 8 Well Microscope Chamber Slides
E36096-23-MS8	Human Mesenchymal (Adipose Derived) Stem Cell ECM 8 well Microscope Slide
E36096-23-MS10	Human Mesenchymal (Adipose Derived) Stem Cell ECM 10 Well Microscope Slides
E36096-23-12Well	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 12 Well Plates (5/Pk)
E36096-23-MS18	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 18 Well Microscope Slides
E36096-23-24Well	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 24 Well Plates (5/Pk)
E36096-23-48Well	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 48 Well Plates (5/Pk)
E36096-23-96Well	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 96 Well Plates (5/Pk)
E36096-23-96WW	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 96 White Well Plates (5/pk)
E36096-23-96BW	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 96 Black Well Plates (5/pk)



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

E36096-23-384WW	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 384 White Well Plates (5/pk)
E36096-23-384BW	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 384 Black Well Plates (5/pk)
E36096-23-CS12	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 12mm coverslips (10/pk)
E36096-23-CS15	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 15mm coverslips (10/pk)
E36096-23-CS18	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 18mm coverslips (10/pk)
E36096-23-CS22	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 22mm coverslips (10/pk)
E36096-23-CS25	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 25mm coverslips (10/pk)
E36096-23-PD6	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 6cm Petri Dish (5/pk)
E36096-23-PD10	Human Mesenchymal (Adipose Derived) Stem Cell ECM - 10cm Petri Dish (5/pk)
E36096-23-6W-GB	Human Mesenchymal (Adipose Derived) Stem Cell ECM 6 Well Glass Bottom Plates (5/pk)
E36096-23-96W-GB	Human Mesenchymal (Adipose Derived) Stem Cell ECM 96 Well Glass Bottom Plates (5/pk)
E36096-23-3D-6Well	Human Mesenchymal (Adipose Derived) Stem Cell 3D Culture System ECM 6 Well (1/pk)
E36096-23-3D-12Well	Human Mesenchymal (Adipose Derived) Stem Cell 3D Culture System ECM 12 Well (1/pk)
E36096-23-3D-24Well	Human Mesenchymal (Adipose Derived) Stem Cell 3D Culture System ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Mesenchymal (Adipose Derived) Stem Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.