



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Embryonic Morula Derived Stem Cell Extracellular Matrix (ECM)

Catalog number: E36123-50

Description: Human Embryonic Morula Derived Stem Cell Extracellular Matrix for growth, expansion, and proliferation of Human Embryonic Morula Derived Stem Cell in tissue culture with their appropriate media. The Human Embryonic Morula Derived Stem Cell ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
E36123-50-T25	Human Embryonic Morula Derived Stem Cell Extracellular Matrix - T25 Flask (10/Pk)
E36123-50-T75	Human Embryonic Morula Derived Stem Cell Extracellular Matrix - T75 Flask (5/Pk)
E36123-50-T150	Human Embryonic Morula Derived Stem Cell Extracellular Matrix - T150 Flasks (5/pk)
E36123-50-T225	Human Embryonic Morula Derived Stem Cell Extracellular Matrix - T225 Flasks (5/pk)
E36123-50-4Well	Human Embryonic Morula Derived Stem Cell ECM - 4 Well Plates (5/Pk)
E36123-50-4MCS	Human Embryonic Morula Derived Stem Cell ECM - 4 Well Microscope Chamber Slides
E36123-50-MS4	Human Embryonic Morula Derived Stem Cell ECM 4 Well Microscope Slide
E36123-50-6Well	Human Embryonic Morula Derived Stem Cell ECM - 6 Well Plates (5/Pk)
E36123-50-8MCS	Human Embryonic Morula Derived Stem Cell ECM - 8 Well Microscope Chamber Slides
E36123-50-MS8	Human Embryonic Morula Derived Stem Cell ECM 8 well Microscope Slide
E36123-50-MS10	Human Embryonic Morula Derived Stem Cell ECM 10 Well Microscope Slides
E36123-50-12Well	Human Embryonic Morula Derived Stem Cell ECM - 12 Well Plates (5/Pk)
E36123-50-MS18	Human Embryonic Morula Derived Stem Cell ECM - 18 Well Microscope Slides
E36123-50-24Well	Human Embryonic Morula Derived Stem Cell ECM - 24 Well Plates (5/Pk)
E36123-50-48Well	Human Embryonic Morula Derived Stem Cell ECM - 48 Well Plates (5/Pk)
E36123-50-96Well	Human Embryonic Morula Derived Stem Cell ECM - 96 Well Plates (5/Pk)
E36123-50-96WW	Human Embryonic Morula Derived Stem Cell ECM - 96 White Well Plates (5/pk)
E36123-50-96BW	Human Embryonic Morula Derived Stem Cell ECM - 96 Black Well Plates (5/pk)
E36123-50-384WW	Human Embryonic Morula Derived Stem Cell ECM - 384 White Well Plates (5/pk)
E36123-50-384BW	Human Embryonic Morula Derived Stem Cell ECM - 384 Black Well Plates (5/pk)
E36123-50-CS12	Human Embryonic Morula Derived Stem Cell ECM - 12mm coverslips (10/pk)
E36123-50-CS15	Human Embryonic Morula Derived Stem Cell ECM - 15mm coverslips (10/pk)
E36123-50-CS18	Human Embryonic Morula Derived Stem Cell ECM - 18mm coverslips (10/pk)



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

E36123-50-CS22	Human Embryonic Morula Derived Stem Cell ECM - 22mm coverslips (10/pk)
E36123-50-CS25	Human Embryonic Morula Derived Stem Cell ECM - 25mm coverslips (10/pk)
E36123-50-PD6	Human Embryonic Morula Derived Stem Cell ECM - 6cm Petri Dish (5/pk)
E36123-50-PD10	Human Embryonic Morula Derived Stem Cell ECM - 10cm Petri Dish (5/pk)
E36123-50-6W-GB	Human Embryonic Morula Derived Stem Cell ECM 6 Well Glass Bottom Plates (5/pk)
E36123-50-96W-GB	Human Embryonic Morula Derived Stem Cell ECM 96 Well Glass Bottom Plates (5/pk)
E36123-50-3D-6Well	Human Embryonic Morula Derived Stem Cell 3D Culture System ECM 6 Well (1/pk)
E36123-50-3D-12Well	Human Embryonic Morula Derived Stem Cell 3D Culture System ECM 12 Well (1/pk)
E36123-50-3D-24Well	Human Embryonic Morula Derived Stem Cell 3D Culture System ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months
 Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen’s Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Embryonic Morula Derived Stem Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.