



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Retinal Müller Cell Culture Extracellular Matrix (ECM)

Catalog number: E37090-01

Description: Human Retinal Müller Cell Culture Extracellular Matrix for growth, expansion, and proliferation of Human Retinal Müller Cell in tissue culture with their appropriate media. The Human Retinal Müller Cell ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop down menu. Each format is tissue culture tested with their appropriate media.

| Catalog | Description |
|------------------|---|
| E37090-01-T25 | Human Retinal Müller Cell Extracellular Matrix - T25 Flask (10/Pk) |
| E37090-01-T75 | Human Retinal Müller Cell Extracellular Matrix - T75 Flask (5/Pk) |
| E37090-01-T150 | Human Retinal Müller Cell Culture Extracellular Matrix - T150 Flasks (5/pk) |
| E37090-01-T225 | Human Retinal Müller Cell Culture Extracellular Matrix - T225 Flasks (5/pk) |
| E37090-01-4Well | Human Retinal Müller Cell Culture ECM - 4 Well Plates (5/Pk) |
| E37090-01-4MCS | Human Retinal Müller Cell Extracellular Matrix - 4 Well Microscope Chamber Slides |
| E37090-01-MS4 | Human Retinal Müller Cell Culture Extracellular Matrix 4 Well Microscope Slide |
| E37090-01-6Well | Human Retinal Müller Cell Extracellular Matrix - 6 Well Plates (5/Pk) |
| E37090-01-8MCS | Human Retinal Müller Cell Extracellular Matrix - 8 Well Microscope Chamber Slides |
| E37090-01-MS8 | Human Retinal Müller Cell Extracellular Matrix 8 well Microscope Slide |
| E37090-01-MS10 | Human Retinal Müller Cell Culture - 10 Well Microscope Slides |
| E37090-01-12Well | Human Retinal Müller Cell Extracellular Matrix - 12 Well Plates (5/Pk) |
| E37090-01-MS18 | Human Retinal Müller Cell Culture - 18 Well Microscope Slides |
| E37090-01-24Well | Human Retinal Müller Cell Extracellular Matrix - 24 Well Plates (5/Pk) |
| E37090-01-48Well | Human Retinal Müller Cell Extracellular Matrix - 48 Well Plates (5/Pk) |
| E37090-01-96Well | Human Retinal Müller Cell Extracellular Matrix - 96 Well Plates (5/Pk) |
| E37090-01-96WW | Human Retinal Müller Cell Culture Extracellular Matrix - 96 White Well Plates (5/pk) |
| E37090-01-96BW | Human Retinal Müller Cell Culture Extracellular Matrix - 96 Black Well Plates (5/pk) |
| E37090-01-384WW | Human Retinal Müller Cell Culture Extracellular Matrix - 384 White Well Plates (5/pk) |
| E37090-01-384BW | Human Retinal Müller Cell Culture Extracellular Matrix - 384 Black Well Plates (5/pk) |
| E37090-01-CS12 | Human Retinal Müller Cell Culture Extracellular Matrix- 12mm coverslips (10/pk) |
| E37090-01-CS15 | Human Retinal Müller Cell Culture Extracellular Matrix- 15mm coverslips (10/pk) |
| E37090-01-CS18 | Human Retinal Müller Cell Culture Extracellular Matrix- 18mm coverslips (10/pk) |
| E37090-01-CS22 | Human Retinal Müller Cell Culture Extracellular Matrix- 22mm coverslips (10/pk) |
| E37090-01-CS25 | Human Retinal Müller Cell Culture Extracellular Matrix- 25mm coverslips (10/pk) |



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

| | |
|---------------------|---|
| E37090-01-PD6 | Human Retinal Müller Cell Culture Extracellular Expansion Matrix - 6cm Petri Dish (5/pk) |
| E37090-01-PD10 | Human Retinal Müller Cell Culture Extracellular Expansion Matrix - 10cm Petri Dish (5/pk) |
| E37090-01-6W-GB | Human Retinal Müller Cell Culture Extracellular Matrix 6 Well Glass Bottom Plates (5/pk) |
| E37090-01-96W-GB | Human Retinal Müller Cell Culture Extracellular Matrix 96 Well Glass Bottom Plates (5/pk) |
| E37090-01-3D-6Well | Human Retinal Müller3D Cell Culture System 6 Well Format (1/pk) |
| E37090-01-3D-12Well | Human Retinal Müller3D Cell Culture System 12 Well Format (1/pk) |
| E37090-01-3D-24Well | Human Retinal Müller3D Cell Culture System 24 Well Format (1/pk) |

Storage Conditions: 2- 8 °C long term storage 18 months
 Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen’s Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Mouses.* This product, Human Retinal Müller Cell Expansion Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.