



1871 N GAFFEY STREET
SUITE A & CEREBELLAR GRANULE NEURON
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Mouse Endothelial Progenitor Stem Cell Extracellular Matrix (ECM)

Catalog number: E66110-37

Description: Mouse Endothelial Progenitor Stem Cell Extracellular Matrix for growth, expansion, and proliferation of Mouse Endothelial Progenitor Stem Cell in tissue culture with their appropriate media. The Mouse Endothelial Progenitor Stem Cell ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
E66110-37-T25	Mouse Endothelial Progenitor Stem Cell Extracellular Matrix - T25 Flask (10/Pk)
E66110-37-T75	Mouse Endothelial Progenitor Stem Cell Extracellular Matrix - T75 Flask (5/Pk)
E66110-37-T150	Mouse Endothelial Progenitor Stem Cell Extracellular Matrix - T150 Flasks (5/pk)
E66110-37-T225	Mouse Endothelial Progenitor Stem Cell Extracellular Matrix - T225 Flasks (5/pk)
E66110-37-4Well	Mouse Endothelial Progenitor Stem Cell ECM - 4 Well Plates (5/Pk)
E66110-37-4MCS	Mouse Endothelial Progenitor Stem Cell ECM - 4 Well Microscope Chamber Slides
E66110-37-MS4	Mouse Endothelial Progenitor Stem Cell ECM 4 Well Microscope Slide
E66110-37-6Well	Mouse Endothelial Progenitor Stem Cell ECM - 6 Well Plates (5/Pk)
E66110-37-8MCS	Mouse Endothelial Progenitor Stem Cell ECM - 8 Well Microscope Chamber Slides
E66110-37-MS8	Mouse Endothelial Progenitor Stem Cell ECM 8 well Microscope Slide
E66110-37-MS10	Mouse Endothelial Progenitor Stem Cell ECM 10 Well Microscope Slides
E66110-37-12Well	Mouse Endothelial Progenitor Stem Cell ECM - 12 Well Plates (5/Pk)
E66110-37-MS18	Mouse Endothelial Progenitor Stem Cell ECM - 18 Well Microscope Slides
E66110-37-24Well	Mouse Endothelial Progenitor Stem Cell ECM - 24 Well Plates (5/Pk)
E66110-37-48Well	Mouse Endothelial Progenitor Stem Cell ECM - 48 Well Plates (5/Pk)
E66110-37-96Well	Mouse Endothelial Progenitor Stem Cell ECM - 96 Well Plates (5/Pk)
E66110-37-96WW	Mouse Endothelial Progenitor Stem Cell ECM - 96 White Well Plates (5/pk)
E66110-37-96BW	Mouse Endothelial Progenitor Stem Cell ECM - 96 Black Well Plates (5/pk)
E66110-37-384WW	Mouse Endothelial Progenitor Stem Cell ECM - 384 White Well Plates (5/pk)
E66110-37-384BW	Mouse Endothelial Progenitor Stem Cell ECM - 384 Black Well Plates (5/pk)
E66110-37-CS12	Mouse Endothelial Progenitor Stem Cell ECM - 12mm coverslips (10/pk)
E66110-37-CS15	Mouse Endothelial Progenitor Stem Cell ECM - 15mm coverslips (10/pk)
E66110-37-CS18	Mouse Endothelial Progenitor Stem Cell ECM - 18mm coverslips (10/pk)
E66110-37-CS22	Mouse Endothelial Progenitor Stem Cell ECM - 22mm coverslips (10/pk)
E66110-37-CS25	Mouse Endothelial Progenitor Stem Cell ECM - 25mm coverslips (10/pk)



1871 N GAFFEY STREET
SUITE A & CEREBELLAR GRANULE NEURON
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

E66110-37-PD6	Mouse Endothelial Progenitor Stem Cell ECM - 6cm Petri Dish (5/pk)
E66110-37-PD10	Mouse Endothelial Progenitor Stem Cell ECM - 10cm Petri Dish (5/pk)
E66110-37-6W-GB	Mouse Endothelial Progenitor Stem Cell ECM 6 Well Glass Bottom Plates (5/pk)
E66110-37-96W-GB	Mouse Endothelial Progenitor Stem Cell ECM 96 Well Glass Bottom Plates (5/pk)
E66110-37-3D-6Well	Mouse Endothelial Progenitor Stem Cell 3D Culture System ECM 6 Well (1/pk)
E66110-37-3D-12Well	Mouse Endothelial Progenitor Stem Cell 3D Culture System ECM 12 Well (1/pk)
E66110-37-3D-24Well	Mouse Endothelial Progenitor Stem Cell 3D Culture System ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices &

Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Mouse Endothelial Progenitor Stem Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.