



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

**Human Fibroblast iPC Derived Embryoid Body (Non-Viral) Extracellular Matrix (ECM)**

**Catalog number:** E836041-01

**Description:** Human Fibroblast iPC Derived Embryoid Body (Non-Viral) Extracellular Matrix for growth, expansion, and proliferation of Human Fibroblast iPC Derived Embryoid Body (Non-Viral) in tissue culture with their appropriate media. The Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

<b>Catalog</b>	<b>Description</b>
E836041-01-T25	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) Extracellular Matrix - T25 Flask (10/Pk)
E836041-01-T75	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) Extracellular Matrix - T75 Flask (5/Pk)
E836041-01-T150	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) Extracellular Matrix - T150 Flasks (5/pk)
E836041-01-T225	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) Extracellular Matrix - T225 Flasks (5/pk)
E836041-01-4Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 4 Well Plates (5/Pk)
E836041-01-4MCS	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 4 Well Microscope Chamber Slides
E836041-01-MS4	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM 4 Well Microscope Slide
E836041-01-6Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 6 Well Plates (5/Pk)
E836041-01-8MCS	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 8 Well Microscope Chamber Slides
E836041-01-MS8	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM 8 well Microscope Slide
E836041-01-MS10	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM 10 Well Microscope Slides
E836041-01-12Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 12 Well Plates (5/Pk)
E836041-01-MS18	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 18 Well Microscope Slides
E836041-01-24Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 24 Well Plates (5/Pk)
E836041-01-48Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 48 Well Plates (5/Pk)



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

E836041-01-96Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 96 Well Plates (5/Pk)
E836041-01-96WW	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 96 White Well Plates (5/pk)
E836041-01-96BW	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 96 Black Well Plates (5/pk)
E836041-01-384WW	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 384 White Well Plates (5/pk)
E836041-01-384BW	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 384 Black Well Plates (5/pk)
E836041-01-CS12	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 12mm coverslips (10/pk)
E836041-01-CS15	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 15mm coverslips (10/pk)
E836041-01-CS18	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 18mm coverslips (10/pk)
E836041-01-CS22	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 22mm coverslips (10/pk)
E836041-01-CS25	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 25mm coverslips (10/pk)
E836041-01-PD6	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 6cm Petri Dish (5/pk)
E836041-01-PD10	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM - 10cm Petri Dish (5/pk)
E836041-01-6W-GB	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM 6 Well Glass Bottom Plates (5/pk)
E836041-01-96W-GB	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) ECM 96 Well Glass Bottom Plates (5/pk)
E836041-01-3D-6Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) 3D Culture System ECM 6 Well (1/pk)
E836041-01-3D-12Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) 3D Culture System ECM 12 Well (1/pk)
E836041-01-3D-24Well	Human Fibroblast iPC Derived Embryoid Body (Non-Viral) 3D Culture System ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices &**



1871 N GAFFEY STREET  
SUITE A & CEREBELLAR GRANULE NEURON  
SAN PEDRO CA 90731  
www.celprogen.com

Phone: 310 547 3975  
Fax: 310 547 2975  
Email: info@celprogen.com  
stemcells@celprogen.com

**Disclaimers:**

***CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.*** This product, Human Fibroblast iPC Derived Embryoid Body (Non-Viral) Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.