



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) Extracellular Matrix (ECM)

Catalog number: E836044-01

Description: Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) Extracellular Matrix for growth, expansion, and proliferation of Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) in tissue culture with their appropriate media. The Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
E836044-01-T25	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) Extracellular Matrix - T25 Flask (10/Pk)
E836044-01-T75	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) Extracellular Matrix - T75 Flask (5/Pk)
E836044-01-T150	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) Extracellular Matrix - T150 Flasks (5/pk)
E836044-01-T225	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) Extracellular Matrix - T225 Flasks (5/pk)
E836044-01-4Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 4 Well Plates (5/Pk)
E836044-01-4MCS	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 4 Well Microscope Chamber Slides
E836044-01-MS4	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM 4 Well Microscope Slide
E836044-01-6Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 6 Well Plates (5/Pk)
E836044-01-8MCS	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 8 Well Microscope Chamber Slides
E836044-01-MS8	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM 8 well Microscope Slide
E836044-01-MS10	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM 10 Well Microscope Slides
E836044-01-12Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 12 Well Plates (5/Pk)



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

E836044-01-MS18	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 18 Well Microscope Slides
E836044-01-24Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 24 Well Plates (5/Pk)
E836044-01-48Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 48 Well Plates (5/Pk)
E836044-01-96Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 96 Well Plates (5/Pk)
E836044-01-96WW	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 96 White Well Plates (5/pk)
E836044-01-96BW	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 96 Black Well Plates (5/pk)
E836044-01-384WW	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 384 White Well Plates (5/pk)
E836044-01-384BW	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 384 Black Well Plates (5/pk)
E836044-01-CS12	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 12mm coverslips (10/pk)
E836044-01-CS15	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 15mm coverslips (10/pk)
E836044-01-CS18	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 18mm coverslips (10/pk)
E836044-01-CS22	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 22mm coverslips (10/pk)
E836044-01-CS25	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 25mm coverslips (10/pk)
E836044-01-PD6	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 6cm Petri Dish (5/pk)
E836044-01-PD10	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM - 10cm Petri Dish (5/pk)
E836044-01-6W-GB	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM 6 Well Glass Bottom Plates (5/pk)
E836044-01-96W-GB	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) ECM 96 Well Glass Bottom Plates (5/pk)
E836044-01-3D-6Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) 3D Culture System ECM 6 Well (1/pk)
E836044-01-3D-12Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) 3D Culture System ECM 12 Well (1/pk)
E836044-01-3D-24Well	Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) 3D Culture System ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).



1871 N GAFFEY STREET
SUITE A & CEREBELLAR GRANULE NEURON
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Cardiomyocyte iPC derived Embryoid Body (Non-Viral) Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.