



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

|  |
|--|
| Human Cardiomyocytes Undifferentiated Extracellular Matrix (ECM) |
|--|

**Catalog number:** U36003-04

**Description:** Human Cardiomyocytes Undifferentiated Extracellular Matrix for growth, expansion, and proliferation of Human Cardiomyocytes Undifferentiated in tissue culture with their appropriate media. The Human Cardiomyocytes Undifferentiated ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

| Catalog          | Description   |
|------------------|---|
| U36003-04-T25    | Human Cardiomyocytes Undifferentiated Extracellular Matrix - T25 Flask (10/Pk)  |
| U36003-04-T75    | Human Cardiomyocytes Undifferentiated Extracellular Matrix - T75 Flask (5/Pk)   |
| U36003-04-T150   | Human Cardiomyocytes Undifferentiated Extracellular Matrix - T150 Flasks (5/pk) |
| U36003-04-T225   | Human Cardiomyocytes Undifferentiated Extracellular Matrix - T225 Flasks (5/pk) |
| U36003-04-4Well  | Human Cardiomyocytes Undifferentiated ECM - 4 Well Plates (5/Pk)                |
| U36003-04-4MCS   | Human Cardiomyocytes Undifferentiated ECM - 4 Well Microscope Chamber Slides    |
| U36003-04-MS4    | Human Cardiomyocytes Undifferentiated ECM 4 Well Microscope Slide               |
| U36003-04-6Well  | Human Cardiomyocytes Undifferentiated ECM - 6 Well Plates (5/Pk)                |
| U36003-04-8MCS   | Human Cardiomyocytes Undifferentiated ECM - 8 Well Microscope Chamber Slides    |
| U36003-04-MS8    | Human Cardiomyocytes Undifferentiated ECM 8 well Microscope Slide               |
| U36003-04-MS10   | Human Cardiomyocytes Undifferentiated ECM 10 Well Microscope Slides             |
| U36003-04-12Well | Human Cardiomyocytes Undifferentiated ECM - 12 Well Plates (5/Pk)               |
| U36003-04-MS18   | Human Cardiomyocytes Undifferentiated ECM - 18 Well Microscope Slides           |
| U36003-04-24Well | Human Cardiomyocytes Undifferentiated ECM - 24 Well Plates (5/Pk)               |
| U36003-04-48Well | Human Cardiomyocytes Undifferentiated ECM - 48 Well Plates (5/Pk)               |
| U36003-04-96Well | Human Cardiomyocytes Undifferentiated ECM - 96 Well Plates (5/Pk)               |
| U36003-04-96WW   | Human Cardiomyocytes Undifferentiated ECM - 96 White Well Plates (5/pk)         |
| U36003-04-96BW   | Human Cardiomyocytes Undifferentiated ECM - 96 Black Well Plates (5/pk)         |
| U36003-04-384WW  | Human Cardiomyocytes Undifferentiated ECM - 384 White Well Plates (5/pk)        |
| U36003-04-384BW  | Human Cardiomyocytes Undifferentiated ECM - 384 Black Well Plates (5/pk)        |
| U36003-04-CS12   | Human Cardiomyocytes Undifferentiated ECM - 12mm coverslips (10/pk)             |
| U36003-04-CS15   | Human Cardiomyocytes Undifferentiated ECM - 15mm coverslips (10/pk)             |
| U36003-04-CS18   | Human Cardiomyocytes Undifferentiated ECM - 18mm coverslips (10/pk)             |
| U36003-04-CS22   | Human Cardiomyocytes Undifferentiated ECM - 22mm coverslips (10/pk)             |
| U36003-04-CS25   | Human Cardiomyocytes Undifferentiated ECM - 25mm coverslips (10/pk)             |



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

|                     |  |
|---------------------|--|
| U36003-04-PD6       | Human Cardiomyocytes Undifferentiated ECM - 6cm Petri Dish (5/pk)            |
| U36003-04-PD10      | Human Cardiomyocytes Undifferentiated ECM - 10cm Petri Dish (5/pk)           |
| U36003-04-6W-GB     | Human Cardiomyocytes Undifferentiated ECM 6 Well Glass Bottom Plates (5/pk)  |
| U36003-04-96W-GB    | Human Cardiomyocytes Undifferentiated ECM 96 Well Glass Bottom Plates (5/pk) |
| U36003-04-3D-6Well  | Human Cardiomyocytes Undifferentiated 3D Culture System ECM 6 Well (1/pk)    |
| U36003-04-3D-12Well | Human Cardiomyocytes Undifferentiated 3D Culture System ECM 12 Well (1/pk)   |
| U36003-04-3D-24Well | Human Cardiomyocytes Undifferentiated 3D Culture System ECM 24 Well (1/pk)   |

**Storage Conditions:** 2- 8 °C long term storage 18 months  
 Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Cardiomyocytes Undifferentiated Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.