



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

### Human Adipose Derived Undifferentiated Extracellular Matrix (ECM)

**Catalog number:** U36010-11

**Description:** Human Adipose Derived Undifferentiated Extracellular Matrix for growth, expansion, and proliferation of Human Adipose Derived Undifferentiated in tissue culture with their appropriate media. The Human Adipose Derived Undifferentiated ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
U36010-11-T25	Human Adipose Derived Undifferentiated Extracellular Matrix - T25 Flask (10/Pk)
U36010-11-T75	Human Adipose Derived Undifferentiated Extracellular Matrix - T75 Flask (5/Pk)
U36010-11-T150	Human Adipose Derived Undifferentiated Extracellular Matrix - T150 Flasks (5/pk)
U36010-11-T225	Human Adipose Derived Undifferentiated Extracellular Matrix - T225 Flasks (5/pk)
U36010-11-4Well	Human Adipose Derived Undifferentiated ECM - 4 Well Plates (5/Pk)
U36010-11-4MCS	Human Adipose Derived Undifferentiated ECM - 4 Well Microscope Chamber Slides
U36010-11-MS4	Human Adipose Derived Undifferentiated ECM 4 Well Microscope Slide
U36010-11-6Well	Human Adipose Derived Undifferentiated ECM - 6 Well Plates (5/Pk)
U36010-11-8MCS	Human Adipose Derived Undifferentiated ECM - 8 Well Microscope Chamber Slides
U36010-11-MS8	Human Adipose Derived Undifferentiated ECM 8 well Microscope Slide
U36010-11-MS10	Human Adipose Derived Undifferentiated ECM 10 Well Microscope Slides
U36010-11-12Well	Human Adipose Derived Undifferentiated ECM - 12 Well Plates (5/Pk)
U36010-11-MS18	Human Adipose Derived Undifferentiated ECM - 18 Well Microscope Slides
U36010-11-24Well	Human Adipose Derived Undifferentiated ECM - 24 Well Plates (5/Pk)
U36010-11-48Well	Human Adipose Derived Undifferentiated ECM - 48 Well Plates (5/Pk)
U36010-11-96Well	Human Adipose Derived Undifferentiated ECM - 96 Well Plates (5/Pk)
U36010-11-96WW	Human Adipose Derived Undifferentiated ECM - 96 White Well Plates (5/pk)
U36010-11-96BW	Human Adipose Derived Undifferentiated ECM - 96 Black Well Plates (5/pk)
U36010-11-384WW	Human Adipose Derived Undifferentiated ECM - 384 White Well Plates (5/pk)
U36010-11-384BW	Human Adipose Derived Undifferentiated ECM - 384 Black Well Plates (5/pk)
U36010-11-CS12	Human Adipose Derived Undifferentiated ECM - 12mm coverslips (10/pk)
U36010-11-CS15	Human Adipose Derived Undifferentiated ECM - 15mm coverslips (10/pk)
U36010-11-CS18	Human Adipose Derived Undifferentiated ECM - 18mm coverslips (10/pk)
U36010-11-CS22	Human Adipose Derived Undifferentiated ECM - 22mm coverslips (10/pk)



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

U36010-11-CS25	Human Adipose Derived Undifferentiated ECM - 25mm coverslips (10/pk)
U36010-11-PD6	Human Adipose Derived Undifferentiated ECM - 6cm Petri Dish (5/pk)
U36010-11-PD10	Human Adipose Derived Undifferentiated ECM - 10cm Petri Dish (5/pk)
U36010-11-6W-GB	Human Adipose Derived Undifferentiated ECM 6 Well Glass Bottom Plates (5/pk)
U36010-11-96W-GB	Human Adipose Derived Undifferentiated ECM 96 Well Glass Bottom Plates (5/pk)
U36010-11-3D-6Well	Human Adipose Derived Undifferentiated 3D Culture System ECM 6 Well (1/pk)
U36010-11-3D-12Well	Human Adipose Derived Undifferentiated 3D Culture System ECM 12 Well (1/pk)
U36010-11-3D-24Well	Human Adipose Derived Undifferentiated 3D Culture System ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months  
 Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen’s Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Adipose Derived Undifferentiated Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.