



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Neuronal Progenitor Stem Cell Undifferentiation Extracellular Matrix (ECM)

Catalog number: U36057-02

Description: Human Neuronal Progenitor Stem Cell Undifferentiation Extracellular Matrix for growth, expansion, and proliferation of Human Neuronal Progenitor Stem Cell Undifferentiation in tissue culture with their appropriate media. The Human Neuronal Progenitor Stem Cell Undifferentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
U36057-02-T25	Human Neuronal Progenitor Stem Cell Undifferentiation Extracellular Matrix - T25 Flask (10/Pk)
U36057-02-T75	Human Neuronal Progenitor Stem Cell Undifferentiation Extracellular Matrix - T75 Flask (5/Pk)
U36057-02-T150	Human Neuronal Progenitor Stem Cell Undifferentiation Extracellular Matrix - T150 Flasks (5/pk)
U36057-02-T225	Human Neuronal Progenitor Stem Cell Undifferentiation Extracellular Matrix - T225 Flasks (5/pk)
U36057-02-4Well	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 4 Well Plates (5/Pk)
U36057-02-4MCS	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 4 Well Microscope Chamber Slides
U36057-02-MS4	Human Neuronal Progenitor Stem Cell Undifferentiation ECM 4 Well Microscope Slide
U36057-02-6Well	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 6 Well Plates (5/Pk)
U36057-02-8MCS	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 8 Well Microscope Chamber Slides
U36057-02-MS8	Human Neuronal Progenitor Stem Cell Undifferentiation ECM 8 well Microscope Slide
U36057-02-MS10	Human Neuronal Progenitor Stem Cell Undifferentiation ECM 10 Well Microscope Slides
U36057-02-12Well	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 12 Well Plates (5/Pk)
U36057-02-MS18	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 18 Well Microscope Slides
U36057-02-24Well	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 24 Well Plates (5/Pk)
U36057-02-48Well	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 48 Well Plates (5/Pk)
U36057-02-96Well	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 96 Well Plates (5/Pk)
U36057-02-96WW	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 96 White Well Plates



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

	(5/pk)
U36057-02-96BW	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 96 Black Well Plates (5/pk)
U36057-02-384WW	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 384 White Well Plates (5/pk)
U36057-02-384BW	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 384 Black Well Plates (5/pk)
U36057-02-CS12	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 12mm coverslips (10/pk)
U36057-02-CS15	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 15mm coverslips (10/pk)
U36057-02-CS18	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 18mm coverslips (10/pk)
U36057-02-CS22	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 22mm coverslips (10/pk)
U36057-02-CS25	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 25mm coverslips (10/pk)
U36057-02-PD6	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 6cm Petri Dish (5/pk)
U36057-02-PD10	Human Neuronal Progenitor Stem Cell Undifferentiation ECM - 10cm Petri Dish (5/pk)
U36057-02-6W-GB	Human Neuronal Progenitor Stem Cell Undifferentiation ECM 6 Well Glass Bottom Plates (5/pk)
U36057-02-96W-GB	Human Neuronal Progenitor Stem Cell Undifferentiation ECM 96 Well Glass Bottom Plates (5/pk)
U36057-02-3D-6Well	Human Neuronal Progenitor Stem Cell Undifferentiation 3D Culture System ECM 6 Well (1/pk)
U36057-02-3D-12Well	Human Neuronal Progenitor Stem Cell Undifferentiation 3D Culture System ECM 12 Well (1/pk)
U36057-02-3D-24Well	Human Neuronal Progenitor Stem Cell Undifferentiation 3D Culture System ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months
 Room temperature for short term storage 12 months
Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Neuronal Progenitor Stem Cell Undifferentiation Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.