



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

### Rat Cardiomyocytes Undifferentiated Extracellular Matrix (ECM)

**Catalog number:** U55003-04

**Description:** Rat Cardiomyocytes Undifferentiated Extracellular Matrix for growth, expansion, and proliferation of Rat Cardiomyocytes Undifferentiated in tissue culture with their appropriate media. The Rat Cardiomyocytes Undifferentiated ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
U55003-04-T25	Rat Cardiomyocytes Undifferentiated Extracellular Matrix - T25 Flask (10/Pk)
U55003-04-T75	Rat Cardiomyocytes Undifferentiated Extracellular Matrix - T75 Flask (5/Pk)
U55003-04-T150	Rat Cardiomyocytes Undifferentiated Extracellular Matrix - T150 Flasks (5/pk)
U55003-04-T225	Rat Cardiomyocytes Undifferentiated Extracellular Matrix - T225 Flasks (5/pk)
U55003-04-4Well	Rat Cardiomyocytes Undifferentiated ECM - 4 Well Plates (5/Pk)
U55003-04-4MCS	Rat Cardiomyocytes Undifferentiated ECM - 4 Well Microscope Chamber Slides
U55003-04-MS4	Rat Cardiomyocytes Undifferentiated ECM 4 Well Microscope Slide
U55003-04-6Well	Rat Cardiomyocytes Undifferentiated ECM - 6 Well Plates (5/Pk)
U55003-04-8MCS	Rat Cardiomyocytes Undifferentiated ECM - 8 Well Microscope Chamber Slides
U55003-04-MS8	Rat Cardiomyocytes Undifferentiated ECM 8 well Microscope Slide
U55003-04-MS10	Rat Cardiomyocytes Undifferentiated ECM 10 Well Microscope Slides
U55003-04-12Well	Rat Cardiomyocytes Undifferentiated ECM - 12 Well Plates (5/Pk)
U55003-04-MS18	Rat Cardiomyocytes Undifferentiated ECM - 18 Well Microscope Slides
U55003-04-24Well	Rat Cardiomyocytes Undifferentiated ECM - 24 Well Plates (5/Pk)
U55003-04-48Well	Rat Cardiomyocytes Undifferentiated ECM - 48 Well Plates (5/Pk)
U55003-04-96Well	Rat Cardiomyocytes Undifferentiated ECM - 96 Well Plates (5/Pk)
U55003-04-96WW	Rat Cardiomyocytes Undifferentiated ECM - 96 White Well Plates (5/pk)
U55003-04-96BW	Rat Cardiomyocytes Undifferentiated ECM - 96 Black Well Plates (5/pk)
U55003-04-384WW	Rat Cardiomyocytes Undifferentiated ECM - 384 White Well Plates (5/pk)
U55003-04-384BW	Rat Cardiomyocytes Undifferentiated ECM - 384 Black Well Plates (5/pk)
U55003-04-CS12	Rat Cardiomyocytes Undifferentiated ECM - 12mm coverslips (10/pk)
U55003-04-CS15	Rat Cardiomyocytes Undifferentiated ECM - 15mm coverslips (10/pk)
U55003-04-CS18	Rat Cardiomyocytes Undifferentiated ECM - 18mm coverslips (10/pk)
U55003-04-CS22	Rat Cardiomyocytes Undifferentiated ECM - 22mm coverslips (10/pk)
U55003-04-CS25	Rat Cardiomyocytes Undifferentiated ECM - 25mm coverslips (10/pk)



1871 N GAFFEY STREET  
 SUITE A & CEREBELLAR GRANULE NEURON  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

U55003-04-PD6	Rat Cardiomyocytes Undifferentiated ECM - 6cm Petri Dish (5/pk)
U55003-04-PD10	Rat Cardiomyocytes Undifferentiated ECM - 10cm Petri Dish (5/pk)
U55003-04-6W-GB	Rat Cardiomyocytes Undifferentiated ECM 6 Well Glass Bottom Plates (5/pk)
U55003-04-96W-GB	Rat Cardiomyocytes Undifferentiated ECM 96 Well Glass Bottom Plates (5/pk)
U55003-04-3D-6Well	Rat Cardiomyocytes Undifferentiated 3D Culture System ECM 6 Well (1/pk)
U55003-04-3D-12Well	Rat Cardiomyocytes Undifferentiated 3D Culture System ECM 12 Well (1/pk)
U55003-04-3D-24Well	Rat Cardiomyocytes Undifferentiated 3D Culture System ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices &**

**Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Rat Cardiomyocytes Undifferentiated Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.