



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

## DATA SHEET

<b>Rat Bone Marrow Derived Stem Cell Culture Extracellular Undifferentiation Matrix (ECM)</b>
---

**Catalog number:** U55093-20

**Description:** Rat Bone Marrow Derived Stem Cell Extracellular Undifferentiation Matrix for growth, Undifferentiation, and proliferation of Rat Bone Marrow Derived Stem Cell in tissue culture with their appropriate media. The Rat Bone Marrow Derived Stem Cell Undifferentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at [www.celprogen.com](http://www.celprogen.com) and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

<b>Catalog</b>	<b>Description</b>
U55093-20-T25	Rat Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix - T25 Flask (10/Pk)
U55093-20-T75	Rat Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix - T75 Flask (5/Pk)
U55093-20-T150	Rat Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix - T150 Flasks (5/pk)
U55093-20-T225	Rat Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix - T225 Flasks (5/pk)
U55093-20-4Well	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 4 Well Plates (5/Pk)
U55093-20-4MCS	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 4 Well Microscope Chamber Slides
U55093-20-MS4	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM 4 Well Microscope Slide
U55093-20-6Well	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 6 Well Plates (5/Pk)
U55093-20-8MCS	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 8 Well Microscope Chamber Slides
U55093-20-MS8	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM 8 well Microscope Slide
U55093-20-MS10	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM 10 Well Microscope Slides
U55093-20-12Well	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 12 Well Plates (5/Pk)
U55093-20-MS18	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 18 Well Microscope Slides
U55093-20-24Well	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 24 Well Plates (5/Pk)
U55093-20-48Well	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 48 Well Plates (5/Pk)
U55093-20-96Well	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 96 Well Plates (5/Pk)
U55093-20-96WW	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 96 White Well Plates (5/pk)
U55093-20-96BW	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 96 Black Well Plates (5/pk)



1871 N GAFFEY STREET  
 SUITE A & B  
 SAN PEDRO CA 90731  
 www.celprogen.com

Phone: 310 547 3975  
 Fax: 310 547 2975  
 Email: info@celprogen.com  
 stemcells@celprogen.com

U55093-20-384WW	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 384 White Well Plates (5/pk)
U55093-20-384BW	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 384 Black Well Plates (5/pk)
U55093-20-CS12	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 12mm coverslips (10/pk)
U55093-20-CS15	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 15mm coverslips (10/pk)
U55093-20-CS18	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 18mm coverslips (10/pk)
U55093-20-CS22	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 22mm coverslips (10/pk)
U55093-20-CS25	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 25mm coverslips (10/pk)
U55093-20-PD6	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 6cm Petri Dish (5/pk)
U55093-20-PD10	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM - 10cm Petri Dish (5/pk)
U55093-20-6W-GB	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM 6 Well Glass Bottom Plates (5/pk)
U55093-20-96W-GB	Rat Bone Marrow Derived Stem Cell Undifferentiation ECM 96 Well Glass Bottom Plates (5/pk)
U55093-20-3D-6Well	Rat Bone Marrow Derived Stem 3D Cell Culture System Undifferentiation ECM 6 Well (1/pk)
U55093-20-3D-12Well	Rat Bone Marrow Derived Stem 3D Cell Culture System Undifferentiation ECM 12 Well (1/pk)
U55093-20-3D-24Well	Rat Bone Marrow Derived Stem 3D Cell Culture System Undifferentiation ECM 24 Well (1/pk)

**Storage Conditions:** 2- 8 °C long term storage 18 months  
 Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

**Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Rat Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.