



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Mouse Pancreatic Undifferentiated Extracellular Matrix (ECM)

Catalog number: U66004-05

Description: Mouse Pancreatic Undifferentiated Extracellular Matrix for growth, expansion, and proliferation of Mouse Pancreatic Undifferentiated in tissue culture with their appropriate media. The Mouse Pancreatic Undifferentiated ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

| Catalog | Description |
|------------------|---|
| U66004-05-T25 | Mouse Pancreatic Undifferentiated Extracellular Matrix - T25 Flask (10/Pk) |
| U66004-05-T75 | Mouse Pancreatic Undifferentiated Extracellular Matrix - T75 Flask (5/Pk) |
| U66004-05-T150 | Mouse Pancreatic Undifferentiated Extracellular Matrix - T150 Flasks (5/pk) |
| U66004-05-T225 | Mouse Pancreatic Undifferentiated Extracellular Matrix - T225 Flasks (5/pk) |
| U66004-05-4Well | Mouse Pancreatic Undifferentiated ECM - 4 Well Plates (5/Pk) |
| U66004-05-4MCS | Mouse Pancreatic Undifferentiated ECM - 4 Well Microscope Chamber Slides |
| U66004-05-MS4 | Mouse Pancreatic Undifferentiated ECM 4 Well Microscope Slide |
| U66004-05-6Well | Mouse Pancreatic Undifferentiated ECM - 6 Well Plates (5/Pk) |
| U66004-05-8MCS | Mouse Pancreatic Undifferentiated ECM - 8 Well Microscope Chamber Slides |
| U66004-05-MS8 | Mouse Pancreatic Undifferentiated ECM 8 well Microscope Slide |
| U66004-05-MS10 | Mouse Pancreatic Undifferentiated ECM 10 Well Microscope Slides |
| U66004-05-12Well | Mouse Pancreatic Undifferentiated ECM - 12 Well Plates (5/Pk) |
| U66004-05-MS18 | Mouse Pancreatic Undifferentiated ECM - 18 Well Microscope Slides |
| U66004-05-24Well | Mouse Pancreatic Undifferentiated ECM - 24 Well Plates (5/Pk) |
| U66004-05-48Well | Mouse Pancreatic Undifferentiated ECM - 48 Well Plates (5/Pk) |
| U66004-05-96Well | Mouse Pancreatic Undifferentiated ECM - 96 Well Plates (5/Pk) |
| U66004-05-96WW | Mouse Pancreatic Undifferentiated ECM - 96 White Well Plates (5/pk) |
| U66004-05-96BW | Mouse Pancreatic Undifferentiated ECM - 96 Black Well Plates (5/pk) |
| U66004-05-384WW | Mouse Pancreatic Undifferentiated ECM - 384 White Well Plates (5/pk) |
| U66004-05-384BW | Mouse Pancreatic Undifferentiated ECM - 384 Black Well Plates (5/pk) |
| U66004-05-CS12 | Mouse Pancreatic Undifferentiated ECM - 12mm coverslips (10/pk) |
| U66004-05-CS15 | Mouse Pancreatic Undifferentiated ECM - 15mm coverslips (10/pk) |
| U66004-05-CS18 | Mouse Pancreatic Undifferentiated ECM - 18mm coverslips (10/pk) |
| U66004-05-CS22 | Mouse Pancreatic Undifferentiated ECM - 22mm coverslips (10/pk) |
| U66004-05-CS25 | Mouse Pancreatic Undifferentiated ECM - 25mm coverslips (10/pk) |



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

| | |
|---------------------|--|
| U66004-05-PD6 | Mouse Pancreatic Undifferentiated ECM - 6cm Petri Dish (5/pk) |
| U66004-05-PD10 | Mouse Pancreatic Undifferentiated ECM - 10cm Petri Dish (5/pk) |
| U66004-05-6W-GB | Mouse Pancreatic Undifferentiated ECM 6 Well Glass Bottom Plates (5/pk) |
| U66004-05-96W-GB | Mouse Pancreatic Undifferentiated ECM 96 Well Glass Bottom Plates (5/pk) |
| U66004-05-3D-6Well | Mouse Pancreatic Undifferentiated 3D Culture System ECM 6 Well (1/pk) |
| U66004-05-3D-12Well | Mouse Pancreatic Undifferentiated 3D Culture System ECM 12 Well (1/pk) |
| U66004-05-3D-24Well | Mouse Pancreatic Undifferentiated 3D Culture System ECM 24 Well (1/pk) |

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices &

Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Mouse Pancreatic Undifferentiated Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.