



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Mouse Bone Marrow Derived Stem Cell Culture Extracellular Undifferentiation Matrix (ECM)

Catalog number: U66093-20

Description: Mouse Bone Marrow Derived Stem Cell Extracellular Undifferentiation Matrix for growth, Undifferentiation, and proliferation of Mouse Bone Marrow Derived Stem Cell in tissue culture with their appropriate media. The Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
U66093-20-T25	Mouse Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix - T25 Flask (10/Pk)
U66093-20-T75	Mouse Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix - T75 Flask (5/Pk)
U66093-20-T150	Mouse Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix - T150 Flasks (5/pk)
U66093-20-T225	Mouse Bone Marrow Derived Stem Cell Undifferentiation Extracellular Matrix - T225 Flasks (5/pk)
U66093-20-4Well	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 4 Well Plates (5/Pk)
U66093-20-4MCS	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 4 Well Microscope Chamber Slides
U66093-20-MS4	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM 4 Well Microscope Slide
U66093-20-6Well	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 6 Well Plates (5/Pk)
U66093-20-8MCS	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 8 Well Microscope Chamber Slides
U66093-20-MS8	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM 8 well Microscope Slide
U66093-20-MS10	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM 10 Well Microscope Slides
U66093-20-12Well	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 12 Well Plates (5/Pk)
U66093-20-MS18	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 18 Well Microscope Slides
U66093-20-24Well	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 24 Well Plates (5/Pk)
U66093-20-48Well	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 48 Well Plates (5/Pk)
U66093-20-96Well	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 96 Well Plates (5/Pk)



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

U66093-20-96WW	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 96 White Well Plates (5/pk)
U66093-20-96BW	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 96 Black Well Plates (5/pk)
U66093-20-384WW	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 384 White Well Plates (5/pk)
U66093-20-384BW	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 384 Black Well Plates (5/pk)
U66093-20-CS12	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 12mm coverslips (10/pk)
U66093-20-CS15	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 15mm coverslips (10/pk)
U66093-20-CS18	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 18mm coverslips (10/pk)
U66093-20-CS22	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 22mm coverslips (10/pk)
U66093-20-CS25	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 25mm coverslips (10/pk)
U66093-20-PD6	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 6cm Petri Dish (5/pk)
U66093-20-PD10	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM - 10cm Petri Dish (5/pk)
U66093-20-6W-GB	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM 6 Well Glass Bottom Plates (5/pk)
U66093-20-96W-GB	Mouse Bone Marrow Derived Stem Cell Undifferentiation ECM 96 Well Glass Bottom Plates (5/pk)
U66093-20-3D-6Well	Mouse Bone Marrow Derived Stem 3D Cell Culture System Undifferentiation ECM 6 Well (1/pk)
U66093-20-3D-12Well	Mouse Bone Marrow Derived Stem 3D Cell Culture System Undifferentiation ECM 12 Well (1/pk)
U66093-20-3D-24Well	Mouse Bone Marrow Derived Stem 3D Cell Culture System Undifferentiation ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Mouse Bone



1871 N GAFFEY STREET
SUITE A & B
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

Marrow Derived Stem Cell Undifferentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.