



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 Stem Cells@celprogen.com

DATA SHEET

Mouse Pancreatic Stem Cell Culture Extracellular Undifferentiation Matrix (ECM)

Catalog number: U66097-24

Description: Mouse Pancreatic Stem Cell Extracellular Undifferentiation Matrix for growth, Undifferentiation, and proliferation of Mouse Pancreatic Stem Cell in tissue culture with their appropriate media. The Mouse Pancreatic Stem Cell Undifferentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
U66097-24-T25	Mouse Pancreatic Stem Cell Undifferentiation Extracellular Matrix - T25 Flask (10/Pk)
U66097-24-T75	Mouse Pancreatic Stem Cell Undifferentiation Extracellular Matrix - T75 Flask (5/Pk)
U66097-24-T150	Mouse Pancreatic Stem Cell Undifferentiation Extracellular Matrix - T150 Flasks (5/pk)
U66097-24-T225	Mouse Pancreatic Stem Cell Undifferentiation Extracellular Matrix - T225 Flasks (5/pk)
U66097-24-4Well	Mouse Pancreatic Stem Cell Undifferentiation ECM - 4 Well Plates (5/Pk)
U66097-24-4MCS	Mouse Pancreatic Stem Cell Undifferentiation ECM - 4 Well Microscope Chamber Slides
U66097-24-MS4	Mouse Pancreatic Stem Cell Undifferentiation ECM 4 Well Microscope Slide
U66097-24-6Well	Mouse Pancreatic Stem Cell Undifferentiation ECM - 6 Well Plates (5/Pk)
U66097-24-8MCS	Mouse Pancreatic Stem Cell Undifferentiation ECM - 8 Well Microscope Chamber Slides
U66097-24-MS8	Mouse Pancreatic Stem Cell Undifferentiation ECM 8 well Microscope Slide
U66097-24-MS10	Mouse Pancreatic Stem Cell Undifferentiation ECM 10 Well Microscope Slides
U66097-24-12Well	Mouse Pancreatic Stem Cell Undifferentiation ECM - 12 Well Plates (5/Pk)
U66097-24-MS18	Mouse Pancreatic Stem Cell Undifferentiation ECM - 18 Well Microscope Slides
U66097-24-24Well	Mouse Pancreatic Stem Cell Undifferentiation ECM - 24 Well Plates (5/Pk)
U66097-24-48Well	Mouse Pancreatic Stem Cell Undifferentiation ECM - 48 Well Plates (5/Pk)
U66097-24-96Well	Mouse Pancreatic Stem Cell Undifferentiation ECM - 96 Well Plates (5/Pk)
U66097-24-96WW	Mouse Pancreatic Stem Cell Undifferentiation ECM - 96 White Well Plates (5/pk)
U66097-24-96BW	Mouse Pancreatic Stem Cell Undifferentiation ECM - 96 Black Well Plates (5/pk)
U66097-24-384WW	Mouse Pancreatic Stem Cell Undifferentiation ECM - 384 White Well Plates (5/pk)
U66097-24-384BW	Mouse Pancreatic Stem Cell Undifferentiation ECM - 384 Black Well Plates (5/pk)
U66097-24-CS12	Mouse Pancreatic Stem Cell Undifferentiation ECM - 12mm coverslips (10/pk)
U66097-24-CS15	Mouse Pancreatic Stem Cell Undifferentiation ECM - 15mm coverslips (10/pk)
U66097-24-CS18	Mouse Pancreatic Stem Cell Undifferentiation ECM - 18mm coverslips (10/pk)
U66097-24-CS22	Mouse Pancreatic Stem Cell Undifferentiation ECM - 22mm coverslips (10/pk)
U66097-24-CS25	Mouse Pancreatic Stem Cell Undifferentiation ECM - 25mm coverslips (10/pk)



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 Stem Cells@celprogen.com

U66097-24-PD6	Mouse Pancreatic Stem Cell Undifferentiation ECM - 6cm Petri Dish (5/pk)
U66097-24-PD10	Mouse Pancreatic Stem Cell Undifferentiation ECM - 10cm Petri Dish (5/pk)
U66097-24-6W-GB	Mouse Pancreatic Stem Cell Undifferentiation ECM 6 Well Glass Bottom Plates (5/pk)
U66097-24-96W-GB	Mouse Pancreatic Stem Cell Undifferentiation ECM 96 Well Glass Bottom Plates (5/pk)
U66097-24-3D-6Well	Mouse Pancreatic Stem3D Cell Culture System Undifferentiation ECM 6 Well (1/pk)
U66097-24-3D-12Well	Mouse Pancreatic Stem3D Cell Culture System Undifferentiation ECM 12 Well (1/pk)
U66097-24-3D-24Well	Mouse Pancreatic Stem3D Cell Culture System Undifferentiation ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months
 Room temperature for short term storage 12 months

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Mouse Pancreatic Stem Cell Undifferentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.