



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Endothelial Undifferentiated Extracellular Matrix (ECM)

Catalog number: UD36006-07

Description: Human Endothelial Undifferentiated Extracellular Matrix for growth, expansion, and proliferation of Human Endothelial Undifferentiated in tissue culture with their appropriate media. The Human Endothelial Undifferentiated ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
UD36006-07-T25	Human Endothelial Undifferentiated Extracellular Matrix - T25 Flask (10/Pk)
UD36006-07-T75	Human Endothelial Undifferentiated Extracellular Matrix - T75 Flask (5/Pk)
UD36006-07-T150	Human Endothelial Undifferentiated Extracellular Matrix - T150 Flasks (5/pk)
UD36006-07-T225	Human Endothelial Undifferentiated Extracellular Matrix - T225 Flasks (5/pk)
UD36006-07-4Well	Human Endothelial Undifferentiated ECM - 4 Well Plates (5/Pk)
UD36006-07-4MCS	Human Endothelial Undifferentiated ECM - 4 Well Microscope Chamber Slides
UD36006-07-MS4	Human Endothelial Undifferentiated ECM 4 Well Microscope Slide
UD36006-07-6Well	Human Endothelial Undifferentiated ECM - 6 Well Plates (5/Pk)
UD36006-07-8MCS	Human Endothelial Undifferentiated ECM - 8 Well Microscope Chamber Slides
UD36006-07-MS8	Human Endothelial Undifferentiated ECM 8 well Microscope Slide
UD36006-07-MS10	Human Endothelial Undifferentiated ECM 10 Well Microscope Slides
UD36006-07-12Well	Human Endothelial Undifferentiated ECM - 12 Well Plates (5/Pk)
UD36006-07-MS18	Human Endothelial Undifferentiated ECM - 18 Well Microscope Slides
UD36006-07-24Well	Human Endothelial Undifferentiated ECM - 24 Well Plates (5/Pk)
UD36006-07-48Well	Human Endothelial Undifferentiated ECM - 48 Well Plates (5/Pk)
UD36006-07-96Well	Human Endothelial Undifferentiated ECM - 96 Well Plates (5/Pk)
UD36006-07-96WW	Human Endothelial Undifferentiated ECM - 96 White Well Plates (5/pk)
UD36006-07-96BW	Human Endothelial Undifferentiated ECM - 96 Black Well Plates (5/pk)
UD36006-07-384WW	Human Endothelial Undifferentiated ECM - 384 White Well Plates (5/pk)
UD36006-07-384BW	Human Endothelial Undifferentiated ECM - 384 Black Well Plates (5/pk)
UD36006-07-CS12	Human Endothelial Undifferentiated ECM - 12mm coverslips (10/pk)
UD36006-07-CS15	Human Endothelial Undifferentiated ECM - 15mm coverslips (10/pk)
UD36006-07-CS18	Human Endothelial Undifferentiated ECM - 18mm coverslips (10/pk)
UD36006-07-CS22	Human Endothelial Undifferentiated ECM - 22mm coverslips (10/pk)
UD36006-07-CS25	Human Endothelial Undifferentiated ECM - 25mm coverslips (10/pk)



1871 N GAFFEY STREET
 SUITE A & CEREBELLAR GRANULE NEURON
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

UD36006-07-PD6	Human Endothelial Undifferentiated ECM - 6cm Petri Dish (5/pk)
UD36006-07-PD10	Human Endothelial Undifferentiated ECM - 10cm Petri Dish (5/pk)
UD36006-07-6W-GB	Human Endothelial Undifferentiated ECM 6 Well Glass Bottom Plates (5/pk)
UD36006-07-96W-GB	Human Endothelial Undifferentiated ECM 96 Well Glass Bottom Plates (5/pk)
UD36006-07-3D-6Well	Human Endothelial Undifferentiated 3D Culture System ECM 6 Well (1/pk)
UD36006-07-3D-12Well	Human Endothelial Undifferentiated 3D Culture System ECM 12 Well (1/pk)
UD36006-07-3D-24Well	Human Endothelial Undifferentiated 3D Culture System ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months
 Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen’s Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans.* This product, Human Endothelial Undifferentiated Cell Culture Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.