



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Cardiomyocyte iPC Derived Stem Cell Dedifferentiation Cell (Lentivirus) ECM

Catalog number: dD736044-03

Description: Human Cardiomyocyte iPC Derived Stem Cell Extracellular Dedifferentiation Matrix for growth, Dedifferentiation, and proliferation of Human Cardiomyocyte iPC Derived Stem Cell in tissue culture with their appropriate media. The Human Cardiomyocyte iPC Derived Stem Cell Dedifferentiation ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop-down menu. Each format is tissue culture tested with their appropriate media and cells.

Catalog	Description
dD736044-03-T25	Human Cardiomyocyte iPC Dedifferentiation Extracellular Matrix - T25 Flask (10/Pk)
dD736044-03-T75	Human Cardiomyocyte iPC Dedifferentiation Extracellular Matrix - T75 Flask (5/Pk)
dD736044-03-T150	Human Cardiomyocyte iPC Dedifferentiation Extracellular Matrix - T150 Flasks (5/pk)
dD736044-03-T225	Human Cardiomyocyte iPC Dedifferentiation Extracellular Matrix - T225 Flasks (5/pk)
dD736044-03-4Well	Human Cardiomyocyte iPC Dedifferentiation ECM - 4 Well Plates (5/Pk)
dD736044-03-4MCS	Human Cardiomyocyte iPC Dedifferentiation ECM - 4 Well Microscope Chamber Slides
dD736044-03-MS4	Human Cardiomyocyte iPC Dedifferentiation ECM 4 Well Microscope Slide
dD736044-03-6Well	Human Cardiomyocyte iPC Dedifferentiation ECM - 6 Well Plates (5/Pk)
dD736044-03-8MCS	Human Cardiomyocyte iPC Dedifferentiation ECM - 8 Well Microscope Chamber Slides
dD736044-03-MS8	Human Cardiomyocyte iPC Dedifferentiation ECM 8 well Microscope Slide
dD736044-03-MS10	Human Cardiomyocyte iPC Dedifferentiation ECM 10 Well Microscope Slides
dD736044-03-12Well	Human Cardiomyocyte iPC Dedifferentiation ECM - 12 Well Plates (5/Pk)
dD736044-03-MS18	Human Cardiomyocyte iPC Dedifferentiation ECM - 18 Well Microscope Slides
dD736044-03-24Well	Human Cardiomyocyte iPC Dedifferentiation ECM - 24 Well Plates (5/Pk)
dD736044-03-48Well	Human Cardiomyocyte iPC Dedifferentiation ECM - 48 Well Plates (5/Pk)
dD736044-03-96Well	Human Cardiomyocyte iPC Dedifferentiation ECM - 96 Well Plates (5/Pk)
dD736044-03-96WW	Human Cardiomyocyte iPC Dedifferentiation ECM - 96 White Well Plates (5/pk)
dD736044-03-96BW	Human Cardiomyocyte iPC Dedifferentiation ECM - 96 Black Well Plates (5/pk)
dD736044-03-384WW	Human Cardiomyocyte iPC Dedifferentiation ECM - 384 White Well Plates (5/pk)
dD736044-03-384BW	Human Cardiomyocyte iPC Dedifferentiation ECM - 384 Black Well Plates (5/pk)
dD736044-03-CS12	Human Cardiomyocyte iPC Dedifferentiation ECM - 12mm coverslips (10/pk)
dD736044-03-CS15	Human Cardiomyocyte iPC Dedifferentiation ECM - 15mm coverslips (10/pk)
dD736044-03-CS18	Human Cardiomyocyte iPC Dedifferentiation ECM - 18mm coverslips (10/pk)



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

dD736044-03-CS22	Human Cardiomyocyte iPC Dedifferentiation ECM - 22mm coverslips (10/pk)
dD736044-03-CS25	Human Cardiomyocyte iPC Dedifferentiation ECM - 25mm coverslips (10/pk)
dD736044-03-PD6	Human Cardiomyocyte iPC Dedifferentiation ECM - 6cm Petri Dish (5/pk)
dD736044-03-PD10	Human Cardiomyocyte iPC Dedifferentiation ECM - 10cm Petri Dish (5/pk)
dD736044-03-6W-GB	Human Cardiomyocyte iPC Dedifferentiation ECM 6 Well Glass Bottom Plates (5/pk)
dD736044-03-96W-GB	Human Cardiomyocyte iPC Dedifferentiation ECM 96 Well Glass Bottom Plates (5/pk)
dD736044-03-3D-6Well	Human Cardiomyocyte iPC 3D Cell Culture System Dedifferentiation ECM 6 Well (1/pk)
dD736044-03-3D-12Well	Human Cardiomyocyte iPC 3D Cell Culture System Dedifferentiation ECM 12 Well (1/pk)
dD736044-03-3D-24Well	Human Cardiomyocyte iPC 3D Cell Culture System Dedifferentiation ECM 24 Well (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers:

CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans. This product, Human Cardiomyocyte iPC Derived Stem Cell Dedifferentiation Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.