



3914 Del Amo Blvd
SUITE 901
TORRANCE, CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Human Breast Cancer Stem Cell Freezing Media

Catalog number:	M36102-29FM
Description:	Human Breast Cancer Stem Cell Freezing Media. This media is available in ready to use format. Media is optimized for freezing Human Breast Cancer Stem Cell Cells. Each lot is qualified to support Human Breast Cancer Stem Cell freezing and thawing for Cell storage.
Storage Conditions:	-20 ⁰ C
Tissue culture:	The Human Breast Cancer Stem Cell Freezing Media was tested and optimized with Human Breast Cancer Stem Cell storage and recovery. The Human Breast Cancer Stem Cell was grown on Celprogen's Human Breast Cancer Stem Cell Extra-cellular Matrix pre-coated flasks and Media. Ideally the cells should be frozen down when they are in their growth phase. Prior to freezing the cells the media should be changed 24 hours before. Freezing cells: a slow and reproducible cooling rate is very important to ensure good recovery of cultures. A decrease of -1 to -3 ⁰ C / minute is recommended for Celprogen's cell culture.
Product Orders:	Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement.
Permits/Forms:	In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits for their country of origin.
Notices & Disclaimers:	<i>CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.</i> The product, Human Breast Cancer Stem Cell Freezing Media, which contains 10% DMSO, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale and may not be transferred to a third party prior to written request and approval by CELPROGEN Inc.