



3914 DEL AMO BLVD
SUITE 901
TORRANCE, CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Product name: Human Immunoglobulins Depleted Solution
--

Catalog number: HS2001-25-1000ml

Description: Human Immunoglobulins Depleted Serum solution (1000ml). This product is optimized for use as an active reagent in tissue culture systems. Each lot is qualified to support Human cell cultures for growth and Proliferation. Formulated for use in a 5% CO₂ humidified incubator.

Storage Conditions: 2- 8 °C for short term storage, -20 °C for long term storage.

Tissue culture: The Human Immunoglobulins Depleted Serum was tested and Optimized with Human cells in Tissue Culture system. The Human Immunoglobulins Depleted Serum was utilized as an active component in Celprogen's Human complete growth media. The Human Immunoglobulins Depleted Serum in complete growth media was optimized on Extra-cellular Matrix and Cell Culture system up to 6 passages within 80 population doublings in a 37°C humidified incubator with 5% CO₂.

Expiration date: 5 years at -20°C

Purity > 98% (Electrophoresis)
Form Liquid
Concentration 5 grams/dL (Nanodrop)
pH 6.5 - 7.4
Sterility No Microbial, Yeast and mold Growth,

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are manufactured under current Good Manufacturing Practices.* The Product, Human Immunoglobulin Depleted Serum, is established and manufactured by CELPROGEN Inc., This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.