



3914 Del Amo Blvd
SUITE 901
TORRANCE, CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Human Primary Coronary Endothelial Freezing Media

Catalog number:	M11002-04FM
Description:	Human Primary Coronary Endothelial Freezing Media. Optimized for freezing Human Primary Coronary Endothelial. Each lot is qualified to support Human Primary Coronary Endothelial freezing and thawing for Cell storage.
Storage Conditions:	-20 ⁰ C
Tissue culture:	The Human Primary Coronary Endothelial Freezing Media was tested and optimized with Human Primary Coronary Endothelial storage and recovery. The Human Primary Coronary Endothelial was grown on Celprogen's Human Primary Coronary Endothelial Extra-cellular Matrix pre-coated flasks and Media. Ideally the cells should be frozen down when they are in their growth phase. Prior to freezing the cells the media should be changed 24 hours before. Freezing cells: a slow and reproducible cooling rate is very important to ensure good recovery of cultures. A decrease of -1 to -3 ⁰ C / minute is recommended for Celprogen's cell culture.
Product Orders:	Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement.
Permits/Forms:	In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.
Notices & Disclaimers:	<i>CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.</i> The product, Human Primary Coronary Endothelial Freezing Media, which contains 10% DMSO, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale and may not be transferred to a third party prior to written request and approval by CELPROGEN Inc.