



3914 DEL AMO BLVD  
SUITE 901  
TORRANCE CA 90503  
www.celprogen.com

Phone: 310 542 8822  
Fax: 310 542 8028  
Email: info@celprogen.com  
stemcells@celprogen.com

## DATA SHEET

### Human Choroidal Endothelial Cell Culture Extracellular Matrix – 10 cm diameter Petri Dish

**Catalog number:** E36052-03-10PD

**Description:** Human Choroidal Endothelial Cell Culture Extracellular Matrix pre-coated 10cm diameter Petri Dish for growth, expansion and proliferation of Human Choroidal Endothelial Cell Cultures. The Human Choroidal Endothelial Cell Culture Extracellular Matrix is available in the following cell culture ware:

#### Well Plates:

6 wells	CAT# E36052-03-6Well plate	5 per pack
12 wells	CAT# E36052-03-12Well plate	5 per pack
24 wells	CAT# E36052-03-24Well plate	5 per pack
96 wells	CAT# E36052-03-96Well plate	5 per pack
384 wells	CAT# E36052-03-384Well plate	5 per pack

#### Flasks:

T25	CAT# E36052-03-T25 flasks	10 per pack
T75	CAT# E36052-03-T75 flasks	5 per pack
T150	CAT# E36052-03-T150 flasks	5 per pack
T225	CAT# E36052-03-T225 flasks	5 per pack

#### Microscope Slides:

8 Well	CAT# E36052-03-8Well slides	5 per pack
10 Well	CAT# E36052-03-10Well slides	5 per pack
18 Well	CAT# E36052-03-18Well slides	5 per pack

#### Microscope Chamber Slides:

8 Well	CAT# E36052-03-8Well slides	5 per pack
4 Well	CAT# E36052-03-4Well slides	5 per pack
2 Well	CAT# E36052-03-2Well slides	5 per pack

#### Microscope Cover Slips:

12mm	CAT# E36052-03-CS12	10 per pack
15mm	CAT# E36052-03-CS15	10 per pack
18mm	CAT# E36052-03-CS18	10 per pack
22mm	CAT# E36052-03-CS22	10 per pack
25mm	CAT# E36052-03-CS25	10 per pack

#### Petri Dishes:

6cm	CAT# E36052-03-PD6	5 per pack
10cm	CAT# E36052-03-PD10	5 per pack



3914 DEL AMO BLVD  
SUITE 901  
TORRANCE CA 90503  
www.celprogen.com

Phone: 310 542 8822  
Fax: 310 542 8028  
Email: info@celprogen.com  
stemcells@celprogen.com

The 6wells, 12wells, 24wells, 96wells, and 384wells are available 5 well plates per pack. The T75, T150, and T225 flasks are available as 5 flasks per pack. The T25 flasks are available as 10 per pack. The 8 Well, 10 Well, and 18 Well Microscope Slides are available as 5 in pack. The 8 Well, 4 Well, and 2 Well Microscope Chamber Slides are available in 5 per pack. The 12mm, 15mm, 18mm, 22mm, and 25mm Microscope Cover Slips are available in 10 per pack. The 6cm and 10cm Petri Dishes are available in 5 per pack. For individual catalog items please visit our web site at [www.celprogen.com](http://www.celprogen.com).

**Storage Conditions:** 2 - 8°C long term storage 24 months (**thaw at room temperature prior to use**)

Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits/Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

#### Notices

**& Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* The product, Human Choroidal Endothelial Cell Culture Extracellular Matrix 10cm diameter Petri Dish is established and manufactured by CELPROGEN Inc., and is for Research Use Only. The product is not for re-sale and may not be transferred to third-party prior to written request and approval by CELPROGEN Inc.