



3914 DEL AMO BLVD.
SUITE 901
TORRANCE, CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Human Pancreatic (Trans-Differentiated) Stem Cell Complete Growth Media with Serum

- Catalog number:** M36004-05S
- Description:** Human Pancreatic (Trans-Differentiated) Stem Cell Complete Growth Media with fetal bovine Serum and Antibiotics: total volume of 500ml. This media is available in ready to use format. Media is optimized for use with Human Pancreatic (Trans-Differentiated) Stem Cell. Each lot is qualified to support Human Pancreatic (Trans-Differentiated) Stem Cell growth and Proliferation. Formulated for use in a 5% CO₂ humidified incubator, and use in low O₂ environment: 5% O₂ in a controlled 37⁰C humidified incubator.
- Storage Conditions:** 2- 8 ⁰C
- Tissue culture:** The Human Pancreatic (Trans-Differentiated) Stem Cell Complete Growth Media with Serum was tested and Optimized with Human Pancreatic (Trans-Differentiated) Stem Cell in tissue culture. The Human Pancreatic (Trans-Differentiated) Stem Cell was grown on Celprogen's Human Pancreatic (Trans-Differentiated) Stem Cell Extra-cellular Matrix pre-coated flasks and Human Pancreatic (Trans-Differentiated) Stem Cell Complete Growth Media with Serum. This media was optimized on Extra-cellular Matrix and Media combination Cell Culture system for up to 12 passages within 80-120 population doublings. This media is formulated to be utilized in a 37⁰C humidified incubator with 5% CO₂. (Protocol is provided with the purchase of product)
- Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).
- Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.
- Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* The product, Human Pancreatic (Trans-Differentiated) Stem Cell Complete Growth Media with Serum, is established and manufactured by CELPROGEN Inc., for Research use only. This product is not for re-sale or may not be transferred to a third party prior to written request and approval by CELPROGEN Inc.