



3914 DEL AMO BLVD.
SUITE 901
TORRANCE, CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Human Pancreatic (Trans-Differentiated) Stem Cell Serum Free Undifferentiation Media

Catalog number: M36004-05U

Description: Human Pancreatic (Trans-Differentiated) Stem Cell Serum Free Undifferentiation Media with Antibiotics: Total Volume of 500ml. Media is available in ready to use format. This media is optimized for use with Human Pancreatic (Trans-Differentiated) Stem Cell . Each lot is qualified to support Human Pancreatic (Trans-Differentiated) Stem Cell undifferentiation, growth and proliferation. Formulated for use in a 5% CO₂ humidified incubator, and use in low O₂ environment: 5% O₂ in a controlled 37⁰C humidified incubator.

Storage Conditions: 2- 8 ⁰C

Tissue culture: The Human Pancreatic (Trans-Differentiated) Stem Cell Serum Free Undifferentiation Media was tested and Optimized with Human Pancreatic (Trans-Differentiated) Stem Cell in Tissue Culture. The Human Pancreatic (Trans-Differentiated) Stem Cell were grown on Celprogen's Human Pancreatic (Trans-Differentiated) Stem Cell Extra-cellular Matrix pre-coated flasks and Human Pancreatic (Trans-Differentiated) Stem Cell Complete Undifferentiation Media. This media was optimized on this Extra-cellular Matrix and Media combination Cell Culture system for up to 12 passages within 80-120 population doublings. The media is formulated to be utilized in a 37⁰C humidified incubator with 5% CO₂. Note: Prior to using serum free media, the cells have to be weaned off their serum conditions. (Protocol is provided with the purchase of the product).

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* The product Human Pancreatic (Trans-Differentiated) Stem Cell Serum Free Undifferentiation Media, is established and manufactured by CELPROGEN Inc., for Research Use Only. This product is not for re-sale or may not be transferred to a third party prior to written request and approval by CELPROGEN Inc.