



3914 DEL AMO BLVD.  
SUITE 901  
TORRANCE, CA 90503  
www.celprogen.com

Phone: 310 542 8822  
Fax: 310 542 8028  
Email: info@celprogen.com  
stemcells@celprogen.com

## DATA SHEET

<b>Human Retinal Stem Cell Serum Free Undifferentiation Media</b>
---

- Catalog number:** M36105-32U
- Description:** Human Retinal Stem Cell Serum Free Undifferentiation Media with Antibiotics: Total Volume of 500ml. Media is available in ready to use format. This media is optimized for use with Human Retinal Stem Cell . Each lot is qualified to support Human Retinal Stem Cell undifferentiation, growth and proliferation. Formulated for use in a 5% CO<sub>2</sub> humidified incubator, and use in low O<sub>2</sub> environment: 5% O<sub>2</sub> in a controlled 37<sup>0</sup>C humidified incubator.
- Storage Conditions:** 2- 8<sup>0</sup>C
- Tissue culture:** The Human Retinal Stem Cell Serum Free Undifferentiation Media was tested and Optimized with Human Retinal Stem Cell in Tissue Culture. The Human Retinal Stem Cell were grown on Celprogen's Human Retinal Stem Cell Extra-cellular Matrix pre-coated flasks and Human Retinal Stem Cell Complete Undifferentiation Media. This media was optimized on this Extra-cellular Matrix and Media combination Cell Culture system for up to 12 passages within 80-120 population doublings. The media is formulated to be utilized in a 37<sup>0</sup>C humidified incubator with 5% CO<sub>2</sub>. Note: Prior to using serum free media, the cells have to be weaned off their serum conditions. (Protocol is provided with the purchase of the product).
- Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).
- Permits & Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.
- Notices & Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* The product Human Retinal Stem Cell Serum Free Undifferentiation Media, is established and manufactured by CELPROGEN Inc., for Research Use Only. This product is not for re-sale or may not be transferred to a third party prior to written request and approval by CELPROGEN Inc.