



3914 DEL AMO BLVD
SUITE 901
TORRANCE CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Product name:	Human Pancreatic Type 1 Langerhans Phenol Free Media with Serum
----------------------	---

Catalog number: M35001-03PNS

Description: Human Pancreatic Type 1 Langerhans Phenol Free and Serum Free Media with fetal bovine Serum and Antibiotics: total volume of 500ml. This media is available in ready to use format. Media is optimized for use with Human Pancreatic Type 1 Langerhans . Each lot is qualified to support Human Pancreatic Type 1 Langerhans Growth and Proliferation. Formulated for use in a 5% CO₂ humidified incubator, and use in low O₂ environment: 5% O₂ in a controlled 37°C humidified incubator. Used for cell based assays, so phenol reagent would not interfere with endpoint colormatric readings.

Storage Conditions: 2- 8 °C

Tissue culture: The Human Pancreatic Type 1 Langerhans Phenol Free Media with fetal bovine Serum and Antibiotics was tested and optimized with Human Pancreatic Type 1 Langerhans Cells in Tissue Culture. The Human Pancreatic Type 1 Langerhans was grown on Celprogen's Human Pancreatic Type 1 Langerhans Extra-cellular Matrix pre-coated flasks. This media was optimized on Extra-cellular Matrix and Media combination Cell Culture system up to 12 passages within 80-120 population doublings in a 37°C humidified incubator with 5% CO₂. (Protocol is provided with the purchase of product)

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices & Disclaimers:

CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans. The Product, Human Pancreatic Type 1 Langerhans Primary Cell Culture Phenol and Serum Free Media, is established and manufactured by CELPROGEN Inc., for Research use only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.