



3914 DEL AMO BLVD
SUITE 901
TORRANCE CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

DATA SHEET

Mouse Mesenchymal Bone Marrow Stem Cell pre-coated - T150 Flask Undifferentiation Extracellular Matrix

Catalog number: U66096-23-T150

Description: Mouse Mesenchymal Bone Marrow Stem Cell pre-coated - T150 Flask Undifferentiation Extracellular Matrix pre-coated for growth, undifferentiation and proliferation of Mouse Mesenchymal Bone Marrow Stem Cell pre-coated - T150 Flask. The Mouse Mesenchymal Bone Marrow Undifferentiation Extracellular Matrix is available in the following cell culture ware:

Well Plates:

6 wells	CAT# 66096-23-6Well plate	5 per pack
12 wells	CAT# 66096-23-12Well plate	5 per pack
24 wells	CAT# 66096-23-24Well plate	5 per pack
96 wells	CAT# 66096-23-96Well plate	5 per pack
384 wells	CAT# 66096-23-384Well plate	5 per pack

Flasks:

T25	CAT# 66096-23-T25 flasks	10 per pack
T75	CAT# 66096-23-T75 flasks	5 per pack
T150	CAT# 66096-23-T150 flasks	5 per pack
T225	CAT# 66096-23-T225 flasks	5 per pack

Microscope Slides:

8 Well	CAT# 66096-23-8Well slides	5 per pack
10 Well	CAT# 66096-23-10Well slides	5 per pack
18 Well	CAT# 66096-23-18Well slides	5 per pack

Microscope Chamber Slides:

8 Well	CAT# 66096-23-8Well slides	5 per pack
4 Well	CAT# 66096-23-4Well slides	5 per pack
2 Well	CAT# 66096-23-2Well slides	5 per pack

Microscope Cover Slips:

12mm	CAT# 66096-23-CS12	10 per pack
15mm	CAT# 66096-23-CS15	10 per pack
18mm	CAT# 66096-23-CS18	10 per pack
22mm	CAT# 66096-23-CS22	10 per pack
25mm	CAT# 66096-23-CS25	10 per pack

Petri Dishes:

6cm	CAT# 66096-23-PD6	5 per pack
10cm	CAT# 66096-23-PD10	5 per pack



3914 DEL AMO BLVD
SUITE 901
TORRANCE CA 90503
www.celprogen.com

Phone: 310 542 8822
Fax: 310 542 8028
Email: info@celprogen.com
stemcells@celprogen.com

The 6wells, 12wells, 24wells, 96wells, and 384wells are available 5 well plates per pack. The T75, T150, and T225 flasks are available as 5 flasks per pack. The T25 flasks are available as 10 per pack. The 8 Well, 10 Well, and 18 Well Microscope Slides are available as 5 in pack. The 8 Well, 4 Well, and 2 Well Microscope Chamber Slides are available in 5 per pack. The 12mm, 15mm, 18mm, 22mm, and 25mm Microscope Cover Slips are available in 10 per pack. The 6cm and 10cm Petri Dishes are available in 5 per pack. For individual catalog items please visit our web site at www.celprogen.com.

Storage Conditions: 2 - 8°C long term storage 24 months (**thaw at room temperature prior to use**)
Room temperature for short term storage 12 months

Product Orders: Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits/Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices

& Disclaimers: *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* The product Mouse Mesenchymal Bone Marrow Stem Cell pre-coated - T150 Flask Undifferentiation Extracellular Matrix is established and manufactured by CELPROGEN Inc., and is for Research Use Only. The product is not for re-sale and may not be transferred to third-party prior to written request and approval by CELPROGEN Inc.